VOCO

Admira® Fusion Flow

EN Instructions for use

MD EU Medical device

In accordance with DIN EN ISO 4049

Product description:

Admira Fusion Flow is a flowable, light-curing, radiopaque nano-hybrid ORMOCER® restorative material. **Admira Fusion Flow** contains 74 % w/w inorganic fillers and is used with a dentine enamel bond.

Shades:

A1, A2, A3, A3.5, A4, B1, B2, C2, Incisal, BL, OA2, WO

Indications:

- Filling of small cavities and extended fissure sealing
- Blocking out undercuts
- Lining or coating cavities
- Fillings of class III - V cavities
- Repairing fillings, veneers and temporary restorations
- Luting translucent prosthetic pieces
- Locking, splinting of loose teeth

Contraindications:

Admira Fusion Flow contains methacrylate and BHT. **Admira Fusion Flow** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:

Admira Fusion Flow is suitable for use in all patients without any age or gender restrictions.

Performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:

Admira Fusion Flow should only be used by a professionally trained dental practitioner.

Use:

Preparation:

Clean the teeth being treated. If applicable, mark occlusal contact points. Bring material to room temperature before use.

Shade selection:

To select the shade, compare the cleaned tooth with the enclosed shade guide prior to anaesthesia. To do so, moisten the **Admira Fusion Flow** sample shades and the tooth and determine the shade, preferably in daylight conditions.

Cavity preparation:

In principle, the cavities should be prepared according to the rules for adhesive restorative treatment; preparation should be minimally invasive in order to protect healthy dental hard tissue. For anterior restorations, angle off the enamel margins. Then clean and dry the cavity. Caries-free lesions in the cervical area do not need to be prepared; thorough cleaning is sufficient in this case.

Creating a dry working field:

Ensure that the work area is sufficiently dry. Use of a rubber dam is recommended.

Matrix placement:

Take appropriate measures for areas where the applied **Admira Fusion Flow** can flow away. Translucent matrices that are wedged in the approximal area are advantageous. Minimal separation facilitates the shaping of the approximal contact and placement of the matrix.

Pulp protection:

A suitable pulp-protective liner that can be covered with a stable cement, if necessary, should be applied in close proximity to the pulp.

Bonding materials:

Admira Fusion Flow is applied with a dentine enamel bond, using the adhesive technique. Any light-curing bonding materials can be used. Observe the relevant instructions for use during preparation (etching technique) and processing.

Application of Admira Fusion Flow:

Lock the enclosed application tip on the syringe with a clockwise turn and directly apply **Admira Fusion Flow**. For fillings with a thickness of over 2 mm, apply and cure in layers. Do not pull the plunger back on the syringe during or after use to guarantee the function of the non-dripping syringe.

The attached application tip is intended for single use only. Remove it after use and seal the syringe again tightly with its original cap. The use of a protective barrier is recommended to protect the syringe against contamination with bodily fluids, dirt from hands or oral tissue.

Light-curing:

Conventional polymerisation devices are suitable for light-curing the material. Curing time if an LED/halogen lamp with a light output of at least 500 mW/cm² is used is 20 s for the shades A1, A2, A3, A3.5, B1, B2, C2, Incisal, BL 40 s for the shades A4, OA2, WO

Bring the light aperture as close as possible to the light-curing unit as otherwise the curing will be less effective. Inadequate curing can cause discolouration and complaints.

Finishing:

The restoration can be finished and polished immediately after removing the matrices (e.g., fine or extra-fine diamond burs, polisher). The tooth should be fluoridated as a final step.

Warnings, precautionary measures:

- Phenolic substances, especially products containing eugenol and thymol, lead to curing failures in restorative composites. Consequently, the use of zinc oxide/eugenol cements or other materials containing eugenol should be avoided in combination with restorative composites.

- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Constituents (in descending order according to content):

Barium aluminium borosilicate glass, ORMOCER® resin, silicon dioxide, initiators, stabilisers, pigments

Storage instructions and application method:

Storage at 4°C-23°C. Seal the syringes immediately after removing the material in order to avoid exposure to light and subsequent polymerisation. Do not use after the expiry date.

Disposal:

Dispose of the product in accordance with local regulations.

Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **Admira Fusion Flow** must be reported to VOCO GmbH and the responsible authority.

Note:

The Summary of Safety and Clinical Performance of **Admira Fusion Flow** can be found in the European database on medical devices (EUDAMED - <https://ec.europa.eu/tools/eudamed>). Detailed information can also be found at www.voco.dental.

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VC 60 AA2816 E1 1121 99 © by VOCO
ORMOCER® is a registered trademark of Fraunhofer Gesellschaft zur Förderung der angewandten Forschung e.V., München, Germany.