Breeze™

Self-Adhesive Resin Cement

Instructions for Use

English

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Hrvatski

Český

Dansk

Nederlands

Eesti

Suomi

Français

Deutsch

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Magyar

Italiano

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Latviešu

Lietuviu

Norsk

Polski

Português

Românește

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Srpski

Slovenský

Slovensko

Español

Svenska

Türkçe

Manufactured by

PENTRON

800-551-0283 / 714-516-7557 / www.pentron.com

Breeze cement is a self-adhesive resin cement indicated for the cementation of most types of restorations. Its high bond strength and low expansion makes Breeze cement the perfect choice for most of your cementation needs.

Available in 3 shades, this new self adhesive cement is an ideal choice for esthetic composite and cementable all-ceramic restorations. By eliminating the need for etching and bonding, Breeze cement reduces both sensitivity and chair time. Breeze cement is packaged in our popular double barrel syringe with automix tips, providing versatility, speed and accuracy during placement; no hand-mixing or trituration devices are required.

INDICATIONS:

- Cementation of composite, cementable all-ceramic, PFM and alloy crowns, bridges, inlays, and onlays.
- Cementation of fiber reinforced composite posts.
- Cementation of titanium, titanium alloy, stainless steel and zirconia posts.

CONTRAINdications:

- Cementation of porcelain veneers.
- Cementation of feldspathic porcelain restorations.
- All-ceramic restorations that require bonding.

Storage and Handling

Storage: Refrigeration required.

Store at 2° - 8°C (36°F - 46°F)

Use: Use at room temperature.

NOTE: Refer to outer package for expiration date. Do not use after expiration date.

Preparation

1. Trial seat the restoration prior to cementation to ensure proper fit and clean thoroughly.
2. Prepare the restoration according to the manufacturer's instructions.
Note: When placing an all-ceramic restoration, apply silane to the internal surface of the restoration and air dry.
3. Prepare and clean the tooth to be restored. Rinse the prepared tooth thoroughly and remove excess water by briefly drying it with a 2 second stream of air. Do not over dry the tooth.

Important Note: Although there is no separate etching or bonding procedure when using Breeze cement, the tooth surface should be free from any contamination such as saliva, blood, eugenol, or any other substance that may interfere with the setting and/or bonding of resin adhesives.

Cementation

1. Prepare the syringe of Breeze cement for initial use.
 - a) Remove the cap from the syringe. Observe the level of cement base and catalyst in the two orifices of the syringe. If the level of material

is not even, bleed excess material from the syringe until the base and catalyst both flow evenly.

- b) Attach the mixing tip to the syringe.

Note: Mixing tips with intra-oral tips are enclosed, if preferred for cement application. Order Item Codes N72AA (50 pk) or N72AB (100 pk) if additional intra-oral tips are desired. Mixing Tips are for single patient use only, to prevent cross-contamination between patients.

2. Gently dispense the cement directly into the restoration.

Note: If cementing a post, use a lentulo spiral to completely coat the canal walls with cement. Use of the intra-oral tips noted above is recommended.

3. Place the restoration firmly and carefully on the prepared tooth and ensure that it is fully seated within a minimum of 1 minute from initial mixing.

Important Note: When placing multiple units, it may be necessary to replace the mixing tip during use. The working time of Breeze cement is a minimum of 1 minute which begins as soon as material is first dispensed through the mixing tip. As a result, working time decreases with each unit that is placed.

4. Light cure the margins of the restoration for one second to gel the cement and facilitate clean up.

5. Remove excess cement with a scaler or explorer.

Note: Breeze Self-Adhesive Resin Cement will completely set within 4 minutes. If the restoration is translucent, it may be light cured for 40 seconds to achieve an immediate set.

Working Characteristics

Mode of Cure:	Dual cure
Method of Mixing:	Auto-Mix
Mixing Ratio (base:catalyst)	1:1
Mixing Time	Auto-mix
Working Time (minimum)*	1 minute
Intra-oral Setting Time (from initial mix)*	Within 4 minutes

* These times may vary based on storage conditions, temperature, humidity, and age of product.

CAUTION: Uncured methacrylate resin may cause contact dermatitis and damage the pulp. Avoid contact with skin, eyes and soft tissue. Wash thoroughly with water after contact.

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Self-Adhesive Resin Cement

LIMITED WARRANTY —

Pentron's technical advice, whether verbal or in writing, is designed to assist dentists in using Pentron's product. The dentist assumes all risk and liability for damages arising out of the improper use of Pentron's product.

In the event of a defect in material or workmanship, Pentron will replace the defective product or part thereof, or reimburse the actual cost of the defective product. In order to take advantage of this limited warranty, the defective product must be returned to Pentron.

PENTRON

Manufactured by Pentron Clinical

1717 West Collins Ave. Orange, CA 92867 USA

800-551-0283 / 714-516-7557 / www.pentron.com

EC REP

European Union Representative
SpofaDental a.s.

Markova 238

CZ-506 01 Jičín

Phone: +420 493 583 111

CE
0086

DMC1637 Rev0
08/2014