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REUSABLE SURGICAL INSTRUMENTS

This procedure applies to all reusable medical devices (hereafter referred to as devices) from the company MEDIN, a.s. for which sterilization before use is required. If any surgical instrument has its own instructions for use manual, then the instructions given in the relevant IFU take precedence over the instructions given in this manual.

1. Caution

- 1.1. Devices covered by this instruction for use manual are supplied non-sterile and must be prepared (disinfected, cleaned and sterilized) according to the procedure below.
- 1.2. The devices are supplied in protective packaging, not intended for sterilization. Remove the device from packaging before first use and remove the coverings from the sharp parts of the device.
- 1.3. If the product's packaging is damaged, it must be returned to the manufacturer for inspection and repackaging.
- 1.4. When cleaning and disinfecting stainless steel equipment, avoid cleaning agents containing chloride ions, which may cause corrosion.
- 1.5. Cleaning solutions containing hydrogen peroxide or stronger alkaline solutions may lead to discolouration and damage to the identification of the product.

- 1.6. For surgical instruments with carbide parts (e.g. needle holders), avoid alkaline solutions, hydrogen peroxide solutions and oxygenated water for cleaning.
- 1.7. Do not use abrasive cleaning materials like wire wool or abrasive detergents for cleaning. If possible, always clean instruments in a surgical instrument washer-disinfector according to ISO 15883-1.
- 1.8. Special attention is required when cleaning long, narrow cannulas, blind holes, and complex parts.
- 1.9. Sterilization parameters only apply to sufficiently cleaned devices.
- 1.10. Sterilization parameters only apply to well-assembled, maintained, calibrated and compliant reusable surgical instruments by ISO 15883-1 and 17665-1.
- 1.11. Only healthcare professionals are allowed to process the implants before use.

2. Limitations on processing

- 2.1. Repeated cleaning, disinfection and sterilization cycles have minimal effect on the instruments.
- 2.2. The life of all reusable equipment is limited by wear and damage. Signs of wear may include but are not limited to contamination, corrosion, pitting corrosion, discolouration, scratches, deformation, bending, kinking, breakage, cracks, and missing or illegible markings. The instrument should always be inspected for completeness, damage or excessive wear after each preparation before sterilization to ensure its functionality.

3. Treatment at the point of use

- 3.1. After using the reusable instruments, the coarse residues, such as tissue and blood debris, should be prevented from drying on the instrument's surface. The instruments should be decontaminated and cleaned as soon as possible after use.
- 3.2. To avoid contamination of unused surgical instruments, personnel and the surrounding environment, separate contaminated surgical instruments from non-contaminated instruments. Contaminated surgical instruments must be separated from medical waste.
- 3.3. Used surgical instruments should be placed in a disinfectant solution in a decontamination bath or decontamination container or covered with a drape moistened with demineralized or potable water to prevent blood or tissue debris from drying.

4. Transportation

- 4.1. Avoid contact between unused and contaminated instruments during transport to prevent contamination.

5. Preparation before cleaning

- 5.1. If possible, open the instruments. Separate pointed and sharp surgical instruments.
- 5.2. Wipe the oil from the locks of new surgical instruments with a cloth. A large amount of oil in the locks could lead to staining on the surface of the surgical instruments after sterilization.
- 5.3. Remove sharp, removable parts of surgical instruments before manual cleaning.
- 5.4. If possible, disassemble the instruments before cleaning. For further information on the possible disassembly and assembly of individual instruments, please refer to the separate instruction for use manual or contact your MEDIN sales representative.

6. Automated cleaning

- 6.1. Always prefer validated automated cleaning and disinfection methods over manual cleaning for reasons of process safety.
- 6.2. We recommend washing instruments before sterilization using automatic washers with thermal disinfection capability.
- 6.3. Place the instruments in the washer so they cannot move freely or overlap during the washing process and prevent possible damage to their surface.
- 6.4. The use of a specific detergent and disinfectant is not prescribed. When selecting a cleaning solution, consider the recommendations of the manufacturer of the washer.
- 6.5. When using detergent, follow the manufacturer's instructions for the relevant product, observing the prescribed concentration, solution temperature, water quality and exposure time.

7. Manual cleaning

- 7.1. Wash the tools by hand under running water and remove coarse dirt mechanically with a brush. Then, wash the instruments in a fresh cleaning and disinfecting solution in an ultrasonic bath and rinse again under running water.
- 7.2. The use of a specific cleaning and disinfecting agent is not prescribed.
- 7.3. When using a cleaning and disinfecting agent, follow the manufacturer's instructions for the relevant product, observing the prescribed concentration, solution temperature, water quality and exposure time. High concentrations or prolonged exposure may cause damage to the instrument.

8. Thermal disinfection

- 8.1. For cleaning in a washer, disinfect at a minimum temperature of 90°C for at least 5 minutes (according to ISO 15883-1).

9. Drying

- 9.1. If the washer does not have a drying cycle, dry each tool thoroughly.

10. Inspection and Maintenance

- 10.1. Inspect all instruments visually for condition and function before sterilization. Instruments should show no signs of wear, as described above. Non-functional instruments showing wear or damage must be discarded and disposed of; reuse is strictly forbidden.
- 10.2. Make sure that the surgical instruments have an undamaged surface and are well-adjusted and functional. Check the proper function, including sharpness of cutting surgical instruments, bending of flexible parts of surgical instruments, movement of joints, hinges and other moving parts.
- 10.3. Ensure there is no visible contamination on the instruments before the sterilization process. If you notice any visible contamination during the visual inspection, please repeat the cleaning steps mentioned above.
- 10.4. Lubricate all moving parts of the instrument. Make sure to use the following substances for lubrication: paraffin/white oil-based, compliant with the applicable European Pharmacopoeia, biocompatible, suitable for steam sterilization and steam permeable.
- 10.5. Unless specified otherwise, reassemble disassembled instruments before sterilization.
- 10.6. Contact your MEDIN sales representative for further detailed instructions on functional checks.