PARASORB® Cone Genta

DE	Gebrauchsanweisung: Bitte aufmerksam lesen!
EN	Instructions for use: Please read carefully!
FR	Instructions d'utilisation: À lire attentivement!
ES	Instrucciones de utilización: ¡Rogamos leer atentamente!
IT	Istruzioni per l'uso: Leggere attentamente!
EL	Οδηγίες χρήσεως: να διαβάζεται με προσοχή!
LT	Naudojimo instrukcija: prašome įdėmiai perskaityti!
RU	Инструкция по применению: пожалуйста, прочитайте внимательно!
RO	Instrucțiuni de utilizare: Vă rugăm cititi cu atenție!
NO	Bruksanvisning: Les nøye!
DA	Brugsanvisningen: Skal læses grundigt!
FI	Käyttöohjeisiin: Luettava huolellisesti.
SV	Bruksanvisningen: Läs noggrant.
NL	Gebruiksaanwijzing: Gelieve aandachtig te lezen!
PT	Instruções de utilização: Por favor, leia atentamente!
CS	Návodu k použití: Čtěte pozorně!
PL	Instrukcjach używania: Należy uważnie przeczytać!

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 HIBC-Code

EN Instructions for use – Please read carefully!

The information given in this package leaflet is updated regularly.
Please read instructions carefully prior to use.

Collagen cone with antibiotic protection for surgical use
- haemostyptic, sterile -

QUALITATIVE AND QUANTITATIVE COMPOSITION

1 cone contains:
Collagen from equine tendons, 22.4 mg
Gentamicin sulfate 16 mg, corresponding to 8.8-11.44 mg gentamicin

PHARMACEUTICAL FORM

Resorbable local haemostyptic agent with antibiotic protection

THERAPEUTIC INDICATIONS

- PARASORB® Cone Genta is inserted for local haemostasis and as a tamponade into clean, clean contaminated and contaminated wound cavities.
- PARASORB® Cone Genta is used for haemostasis following extractions or maxillary surgical intervention with a high risk of infection. In patients at risk from bleeding a combination with fibrin sealant is indicated.
The information in the relevant Instructions for Use must be observed.
- PARASORB® Cone Genta is also used for reducing alveolar ridge atrophy following tooth extractions.

POSODOLOGY AND METHOD OF ADMINISTRATION

After removing the focus of inflammation PARASORB® Cone Genta (according to requirement) is inserted into the defect cavity. PARASORB® Cone Genta is loosely placed or if there is a bone defect mixed with tiny bone fragments and put in place. To reduce alveolar ridge atrophy, following tooth extraction PARASORB® Cone Genta is introduced into the alveolus. Because collagen is resorbable no removal is needed.

Haemostasis takes place by a physical method. The contact of collagen with blood leads to platelet aggregation. Large numbers of platelets precipitate on the collagen matrix, disintegrate and release coagulation factors, which together with plasma factors facilitate the formation of fibrin. The collagen matrix ensures additional strengthening of the coagulum and thus inhibits connective tissue healing of the extraction alveolus. Furthermore the gentamicin contained in PARASORB® Cone Genta produces enhanced vascularisation of the coagulum and granulation tissue. In this way atrophy of the alveolar ridge is reduced and the quality of newly-formed bone improved (ArabOghli et al., 2006; Europeo). PARASORB® Cone Genta is able, because of its structure, to absorb large quantities of fluids. Through this purely mechanical process, by means of secretion absorption, shed materials such as bacteria and fibrin layers are absorbed.

The dosage of PARASORB® Cone Genta is governed by the size of the defect and the weight of the patient. Dependent upon this, the maximum permissible quantity of gentamicin sulfate for patients up to 50 kg body weight is 600 mg gentamicin sulfate or for patients over 50 kg body weight 1000 mg gentamicin sulfate. PARASORB® Cone Genta is applied in a dry state by exerting slight pressure in order to achieve better adhesion. Unless otherwise directed, PARASORB® Cone Genta is inserted directly into the extraction alveolus. PARASORB® Cone Genta adapts to the extraction alveolus, brings about haemostasis and is resorbed. Because of collagen's affinity for bloody surfaces, when PARASORB® Cone Genta is applied, dry gloves and instruments should be used.

The release of the incorporated gentamicin happens concurrently with resorption of the collagen cone, so that the cone is protected from external contamination and also from contamination by the introduction and spread of germs during surgery. In the event of diffuse bleeding, pressure when applying the cone at the affected sites in the wound area can be advantageous. Overloading the wound with collagen can lead to the formation of a seroma. Care must be taken that only a single layer of collagen is applied to the wound area. With a combination of PARASORB® Cone Genta and a fibrin sealant the latter is applied on the side which will be placed on the area to be treated.

CONTRAINDICATIONS

PARASORB® Cone Genta must not be used if there is a known sensitivity to collagen and/or gentamicin or other aminoglycoside antibiotics. There are no data about the safe use of PARASORB® Cone Genta during pregnancy and lactation available. Aminoglycoside antibiotics can pass through the placenta and reach the foetus. There is therefore a risk of intrauterine damage. Since aminoglycoside antibiotics also pass into breast milk they should not be used during lactation, or breastfeeding should be discontinued.

SPECIAL WARNINGS AND PRECAUTIONS FOR USE

PARASORB® Cone Genta should only be used under strict indication guidelines in patients with limited renal function also patients with autoimmune diseases or neuromuscular diseases such as Parkinson's or myasthenia gravis.
Simultaneous systemic administration of aminoglycoside antibiotics should be avoided or only given if there is strict monitoring of serum-gentamicin levels and renal function. Cross-allergies between aminoglycosides exist. There are no records about the use of PARASORB® Cone Genta in patients with immunological diseases or connective tissue diseases available. Although there has been no evidence found that the use of animal collagen leads to a deterioration in disease PARASORB® Cone Genta should only be used under strict indication guidelines in these patients. There are no records about its use in paediatrics available.

To guarantee full functionality of the product, the usual application techniques for a haemostypticum should be observed. For safe application the following instructions for its application must be observed:

- Moistening or wetting the haemostypticum prior to implantation can lead to impairment in its action both in terms of haemostyptic properties and also with respect to a loss in the inherent protective characteristics of PARASORB® Cone Genta because of a premature release of the highly water-soluble gentamicin sulphate. Moistening causing a reduction in blood absorption of the cone leads to a reduction in the physical properties of platelet aggregation and therefore considerably affects haemostasis.
- The desired platelet aggregation is also achieved by lining the wound with PARASORB® Cone Genta.

UNDESIRABLE EFFECTS

Under some circumstances at the start of tissue resorption by the PARASORB® Cone Genta, increased secretion may occur. With simultaneous systemic administration of aminoglycoside antibiotics or in patients with limited renal function an increase in nephrotoxicity and ototoxicity may occur.

INTERACTION

Although after implantation of PARASORB® Cone Genta only very low plasma levels are reached, attention should be paid to the interactions of gentamicin with other drugs. A simultaneous administration of aminoglycosides with loop diuretics such as furosemide or ethacrynic acid must be avoided since loop diuretics themselves have an ototoxic effect. If i.v. diuretics are given at the same time these can increase the concentration of gentamicin in serum or tissue and boost toxicity. Simultaneous or successive systemic administration or topical application of potentially neurotoxic and/or nephrotoxic substances such as Cisplatin, other aminoglycosides, Streptomycin, Cephaloridine, Viomycin, Polymyxin B or Polymyxin E can increase the toxicity of gentamicin.

Simultaneous local use of beta-lactam antibiotics can lead to changed action (inactivation). A tendency to neuromuscular blockade can be increased by the simultaneous administration of muscle relaxants, e.g. d-Tubocurarine, Suxamethonium or Pancuronium, as well as ether. A neuromuscular blockade can best be remedied by calcium salt. If a neuromuscular blockade occurs with simultaneous administration of Suxamethonium in the case of acquired or genetic cholinesterase deficit, artificial respiration must be given and cholinesterase administered.

EFFECTS ON ABILITY TO DRIVE AND USE MACHINES

To date none known.

OVERDOSE

The recommended application quantity of PARASORB® Cone Genta causes a negligible increase in serum-gentamicin levels. Regarding dosage over and above this there are no known findings available. With a deviation from this dosage in each case the serum gentamicin levels and kidney function must be monitored. In cases of severe intoxication, a peritoneal dialysis or haemodialysis may be considered.

INCOMPATIBILITIES

To date none known.

PLEASE NOTE

Because of the affinity of collagen for bloody surfaces PARASORB® Cone Genta should only be applied using dry gloves and instruments.

SPECIAL PRECAUTIONS FOR STORAGE

Do not store at temperatures exceeding 25° C.
Protect against moisture. Keep out of the reach and sight of children.

SHELF LIFE

The use-by date and batch number are printed on the sales packaging and the blister cover. After the use-by date has expired PARASORB® Cone Genta must no longer be used. PARASORB® Cone Genta once removed from the protective packaging must not be re-sterilised and must therefore be discarded. The contents of unused, opened or damaged packaging must not be re-sterilised and must therefore be discarded. PARASORB® Cone Genta is intended for single use only. PARASORB® Cone Genta, pack size MK11, is packaged in double sterilised packaging. The interior package is sterile inside and outside and therefore can be passed on in a sterile state.

PACK SIZES

PARASORB® Cone Genta	PARASORB® Cone Genta
10 pcs, size Ø: 1.2 cm; H: 1.6 cm	5 pcs, size Ø: 1.2 cm; H: 1.6 cm
REF MK10	REF MK11

DESCRIPTION OF THE SYMBOLS USED ON THE PACKAGE

-  Reference number
-  Batch number
-  Use by year - month
-  Consult instructions for use
-  Do not reuse
-  Do not resterilize
-  Do not use if package is damaged
-  Sterilised using ethylene oxide
-  Upper limit of temperature
-  Protect against humidity

CE 0197

CE marking and identification number of the notified body. Product conforms to the essential requirements of the Council Directive 93/42/EEC concerning medical devices.

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COMPOSITION QUALITATIVE ET QUANTITATIVE

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RESORBA Medical GmbH