VOCO

Grandio®

Table 1

Light intensity (LED-/halogen lamp)	Polymerisation time (per layer)	
≥ 500 mW/cm ²	20 s A1, A2, A3, A3.5, A4, B1, B2, B3, C2, C3, D2, D3, Incisal, BL	40 s OA2, OA3.5
≥ 800 mW/cm ²	10 s A1, A2, A3, A3.5, A4, B1, B2, B3, C2, C3, D2, D3, Incisal, BL	20 s OA2, OA3.5

EN Instructions for use

MD EU Medical device
In accordance with DIN EN ISO 4049

Product description:

Grandio is a light-curing, radiopaque nano-hybrid restorative material. **Grandio** contains 87% by weight inorganic fillers and is applied with a dentine enamel bond. It is available in screw syringes and caps for direct application.

Shades:

A1, A2, A3, A3.5, A4, B1, B2, B3, C2, C3, D2, D3, Incisal, OA2, OA3.5, BL

Indications:

- Class I to V fillings
- Reconstruction of traumatically affected anteriors
- Facetting of discoloured anteriors
- Correction of shape and shade for better aesthetic appearance
- Locking, splinting of loose teeth
- Repairing veneers
- Filling deciduous teeth
- Core-build-up under crowns
- Composite inlays

Contraindications:

Grandio contains methacrylate and BHT. **Grandio** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:

Grandio is suitable for use in all patients without any age or gender restrictions.

Performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:

Grandio should only be used by a professionally trained dental practitioner.

Use:

Preparation:

Clean the teeth being treated. If applicable, mark occlusal contact points. Bring material to room temperature before use.

Shade selection:

To select the shade, compare the cleaned tooth with the enclosed shade guide prior to anaesthesia. To do so, moisten the **Grandio** sample shades and the tooth and determine the shade, preferably in daylight conditions.

Cavity preparation:

In principle, the cavities should be prepared according to the rules for adhesive restorative treatment; preparation should be minimally invasive in order to protect healthy dental hard tissue. For anterior restorations, angle off the enamel margins. Then clean and dry the cavity. Caries-free lesions in the cervical area do not need to be prepared; thorough cleaning is sufficient in this case.

Ensure that the working field is sufficiently dry. The use of a rubber dam is recommended. It is recommended that shaping aids be used in cases where a portion of the cavity is in the proximal region. It helps to use translucent female parts, wedged in the proximal region. Keeping separation to the minimum makes it easier to shape the proximal contacts and attach a female part. Suitable protection for the pulp should be applied close to the pulp and covered with a layer of stable cement if necessary.

Bonding materials:

Grandio is applied with a dentine enamel bond, using the adhesive technique. Any light-curing bonding materials can be used. Observe the relevant instructions for use during preparation (etching technique) and processing.

Application:

Apply the selected shade of **Grandio** (see shade selection) in layers of no more than 2 mm thickness; adapt using a suitable instrument and then cure.

Light-curing:

Conventional polymerisation devices are suitable for light-curing the material. The polymerisation time per layer (Polymerisation time) as a function of the light intensity (Light intensity) can be found in Table 1 (Table 1).

Bring the light aperture as close as possible to the light-curing unit as otherwise the curing will be less effective. Inadequate curing can cause discolouration and complaints.

Finishing:

Restorations can be finished and polished immediately after removing the female parts (e.g. fine and extra-fine diamond burs, abrasive discs). The tooth should then be fluoridized.

Grandio can be used for direct or indirect inlays using the known procedures. The physical stabilisation values are improved by the usual external (additional) curing.

Warnings, precautionary measures:

– Phenolic substances, especially products containing eugenol and thymol, lead to curing failures in restorative composites. Consequently, the use of zinc oxide/ eugenol cements or other materials containing eugenol should be avoided in combination with restorative composites.

– Each **Grandio** Caps unit should be used for one patient only and should be applied exerting uniform but not excessive pressure.

– Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Constituents (in descending order according to content):

Strontium aluminium borosilicate glass, barium aluminium borosilicate glass, silicon dioxide, BisGMA, TEGDMA, initiators, stabilisers, pigments

Storage instructions and application method:

Storage at 4°C-28°C. Seal the syringes immediately after removing the material in order to avoid exposure to light and subsequent polymerisation. Do not use after the expiry date.

Disposal:

Dispose of the product in accordance with local regulations.

Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **Grandio** must be reported to VOCO GmbH and the responsible authority.

Note:

The Summary of Safety and Clinical Performance of **Grandio** can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>). Detailed information can also be found at www.voco.dental.

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© by VOCO