VOCO

# Rebilda Post

0482

## Instructions for use

The radiopaque, translucent **Rebilda Post** is made of glassfibre-reinforced composite and exhibits elasticity behaviour similar to that of dentine. A durable, highly aesthetic metal-free restoration can be achieved when used in the adhesive technique. **Rebilda Post** is available in four sizes with corresponding drills.

Size	10 (red)	12 (green)	15 (black)	20 (yellow)
Length	19 mm	19 mm	19 mm	19 mm
Coronal diameter	1,00 mm	1,20 mm	1,50 mm	2,00 mm
Apical diameter	0,50 mm	0,65 mm	0,80 mm	1,00 mm

### Indication:

Endodontic post build-up for support and anchorage of coronal restorations in situations with insufficient tooth substance.

### Contraindications:

- Allergies to methacrylates
- Patients with bruxism or closed bite
- Less than 1.5 mm circular coronal tooth substance

### Application instructions:

#### A) Pre-treatment of root canal and Rebilda Post

- Successful completion of endodontic treatment
- Utilisation of a rubber dam is recommended
- Selection of the appropriate post size based on x-ray
- Correct determination of the required working length (apically, a minimum of 4 mm endodontic filling must be retained, ½ of the post should be located in the root canal and ½ in the core build-up)
- Removal of the root canal filling material with a Gates-Glidden bur, Peeso reamer and/or other hot instrument down to the proposed depth
- Precision drilling with an appropriate instrument (1000 - 2000 min⁻¹ instrument rpm; low speed contra angle hand piece). Avoid preparation beyond the apex and/or lateral perforation of the root. [Drills and posts must always be used as a system (cp. colour coding) and may not be combined with other brands]
- To ensure optimal drilling performance, intermittently rinse the root canal and drill to remove dentine residues.
- Clean the post with alcohol before the Try-in
- Check the position in the mouth (the endodontic post should exactly fill the root canal; it should not be wedged and should be easily removable)
- Shorten the post to the appropriate length with diamond tools extra-orally (Do not use forceps or scissors to avoid the risk of delamination)
- Clean the post with alcohol and dry
- Carry out silanisation of post (e.g. Ceramic Bond), allow to act for 60 s and then dry with oil-free air
- Flush the root canal with a 2.5 - 5 % sodium hypochlorite solution (NaOCl). [Do not use hydrogen peroxide (H₂O₂), other disinfecting agents or EDTA solutions, since residues can impair the adhesive strength and/or curing reaction of the luting composite]
- Flush immediately with water and dry with paper points

#### B) Adhesive luting with Rebilda DC

(in combination with a dual-curing adhesive)

1. When luting **Rebilda Post**, apply the dual-curing adhesive according to its instructions for use.

#### 2. Application of Rebilda DC cartridge:

When using mixing tip type 19, apply **Rebilda DC** directly into the root canal. If necessary, the application needle can be pre-bent by hand for this purpose. When not using mixing tip type 19, proceed as follows:  
Wet **Rebilda Post** with **Rebilda DC**. Do not apply **Rebilda DC** directly into the root canal, do not use a Lentulo! Insert **Rebilda Post** in a twisting motion while allowing a slight amount of excess cement. Remove the excess cement if **Rebilda DC** is not being used as core build-up material.

#### Application of Rebilda DC QuickMix syringe:

Applying **Rebilda DC** directly into the root canal with the assistance of a thin intraoral tip type 1 is recommended (application beginning from apical direction, keep the cannula immersed in the luting composite during the application). Do not use a Lentulo! **Rebilda Post** can additionally be wetted with **Rebilda DC**. Insert **Rebilda Post** in a twisting motion while allowing a slight amount of excess cement. Remove the excess cement if **Rebilda DC** is not being used as core build-up material.

3. To fix the post and cure the dual-curing adhesive applied to the occlusal surface, light-cure for 40 s.

If using a luting composite from another manufacturer, proceed according to the respective information provided by the manufacturer.

Avoid the use of all materials that contain eugenol, including some normally appropriate endodontic restorative, as they can impede curing of the luting composite.

#### C) Core build-up

- Follow instructions for use for the selected core build-up material (e.g. **Rebilda DC**).

#### D) Removal of the Rebilda Post in the event of a revision:

- Take an x-ray to determine the location, length and size of the post
- Remove the restoration and core build-up in order to facilitate direct, straight access to the post
- Drill a central pilot hole in the middle of the post with a round bur
- Take an x-ray, if necessary, to check the course of the drill
- Drill **Rebilda Post** out with the appropriate drill (Attention: Due to the similar behaviour of the rotating instrument in the glassfibre-reinforced endodontic post and in the radicular dentine, it is necessary to use extreme caution to avoid endodontic perforation. If necessary, check intermediate steps with an x-ray.)

#### E) Cleaning, disinfection and sterilisation (Endodontic posts and drills are not delivered sterile):

- Clean and autoclave (134 °C for 18 min.) drill before each use. Use a separate autoclave bag for each instrument
- Use a suitable agent for cleaning and disinfecting rotating instruments (in accordance with EN 1040 and EN 1275). Follow the respective manufacturer's instructions. Some aldehyde-free cleaning or disinfecting solutions can cause discolouration of the colour-coded rings. Use a brush or ultrasonic bath (without contact to other instruments) for cleaning. Rinse well with water and dry afterwards. Verify that cleaning has been thorough and repeat if necessary
- Always store drills dry and in a dust-free container. Visually check the drill before each use. Incompletely cleaned drills have to be cleaned again. Separate out damaged and/or dull drills and discontinue using them
- The posts can be disinfected with medicinal alcohol

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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