1 Intended use

Intended purpose

Sealing of pits, fissures and foramina caeca as a protection against caries

Patient target group

- Patients with permanent teeth
- Patients with deciduous teeth

Intended users / Special training

- Dentists
- No special training required

Use

For dental use only.

Description

Helioseal® F Plus is a light-curing, white-shaded fissure sealant (acc. to EN ISO 6874:2015) featuring fluoride release. Helioseal F Plus is used to seal pits, fissures and foramina caeca as a protection against caries.

Technical data

The product contains a photoinitiator and cures with light as follows:

Light intensity / wave length	Exposure time
> 500 mW/cm ² / 400-500 nm	20 s
> 1000 mW/cm ² / 400-500 nm	10 s
2000 mW/cm ² / 400-500 nm	5 s

The curing depth achieved with an exposure time of 10 s (at a light intensity of > 1000 mW/cm² and a wave length of 400-500 nm) is > 3 mm.

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Areas of application

Helioseal F Plus is used to seal pits, fissures and foramina caeca.

Contraindications

The use of the product is contraindicated if the patient is known to be allergic to any of its ingredients.

Limitations of use

- If a dry working field cannot be established.
- If the stipulated working procedures cannot be applied.
-  Cavifils and cannulas are intended for single use only.

Do not re-use

Helioseal® F Plus

Helioseal® F Plus

[en]	Instructions for Use	Light-cured fissure sealant	[fi]	Käyttöohjeet	Valokovettainen fissuurapainotte
[de]	Gebrauchsinformation	Lichthärtender Fissurenversiegler	[no]	Bruksanvisning	lyshærdende fissurforsæglings
[fr]	Mode d'emploi	Matériau de scellement de sillons photopolymérisable	[nl]	Gebruiksaanwijzing	Lichthardtende fissuursealant
[it]	Istruzioni d'uso	Sigillante per fessure fototindente	[el]	Οδηγίες Χρήσεως	Φωτοπολυμεριζόμενο εμφορστικό υλικό
[es]	Instrucciones de uso	Sellador de fisuras fotopolimerizable	[tr]	Kullanma Talimatı	Işıklı polimerize olan fissür dolgusu
[pt]	Instruções de Uso	Selante de fissuras fotopolimerizável	[ru]	Инструкция по применению	Светоотверждаемый материал для запечатывания fissур
[sv]	Bruksanvisning	Ljushärdande fissurförsegling	[pl]	Instrukcja stosowania	Światłoutwardzalny materiał szlącący do uszczelniania bruz
[da]	Brugsanvisning	lyshærdende resin til fissurforsæglings	[mk]	Упатство за употреба	Светлосно-полимеризирачки запевач за fissури
[sl]	Navodila za uporabo	Sredstvo za tesnjenje fissur, ki se strjuje s svetlobo	[bg]	Инструкции за употреба	Фотополимеризиращ силант за fissури
[hr]	Uputa za uporabu	Svijetlosno polimerizirajući materijal za pečaćenje fissura	[sq]	Udhëzime përdorimi	Restaurues plasarëjtje shi i foto polimerizueshem
[cs]	Návod k použití	Světlem tuhnoucí materiál na pečetení fissur	[ro]	Instrucţiuni de utilizare	Material de sigilare pentru şanţuri şi foseţe, fotopolimerizabil
[sk]	Návod na používanie	Svetlom vytvrdzované pečatidlo na pečatenie rýh	[uk]	Інструкція щодо використання	Фотополімерний матеріал для герметизації fissур
[hu]	Használati utasítás	Fényre kötő barázdázó anyag	[et]	Kasutamisi juhend	Valguskõvastuv lõhede täitematerjal
[sr]	Упутство за употребу	Светлосно-полимеризујући запивач fissура	[lv]	Lietošanas instrukcija	Ģaismā cietināms plaisu silants
[mk]	Упатство за употреба	Светлосно-полимеризирачки запевач за fissури	[lt]	Naudojimo instrukcijos	Šviesą kietinamas vagelių silantas

Side effects

In rare cases, components of the product may lead to sensitization. The product must not be used in such cases.

Interactions

Substances containing eugenol/clove oil may inhibit the polymerization of methacrylate-based materials. Consequently, the application of such materials together with the product must be avoided.

Clinical benefit

Protection against caries

Composition

UDMA, HEMA phosphate, aromatic aliphatic UDMA, Al fluorosilicate glass, silicon dioxide, polyacrylate

2 Application

I. Cleaning the enamel surface

Thoroughly clean the enamel surface to be sealed.

II. Isolation

Adequate relative or absolute isolation with rubber dam is required.

III. Conditioning

Condition the enamel surface with an etching gel containing phosphoric acid according to the instructions for use of the product used.

IV. Thorough rinsing and drying

- Thoroughly rinse off the conditioner with water spray.
- Dry the surface with oil-/moisture-free air.
- Etched enamel has a matte white appearance; avoid contamination with saliva of the surface to be sealed. In the case of contamination, the conditioning step must be repeated.

V. Applying and light-curing the sealant

- Apply the product directly from the single-use cannula or the Cavifil and distribute it evenly on the tooth surface.
- Ensure sufficient exposure to the curing light to prevent incomplete polymerization of the sealant.
- For recommendations regarding exposure time and light intensity see section "Technical data".
-  The instructions for use of the curing light must be observed.
- Check the sealant and the occlusion.

Application notes

- Use a suitable hygienic protective sleeve for the syringe as indicated by the manufacturer.
- Avoid any contact of the reusable syringe not covered with a protective sleeve with the mouth of the patient.
- Dispose of contaminated syringes.
- Do not disinfect syringes with oxidizing disinfection agents.
- To avoid the formation of air bubbles in the syringe, the plunger of the syringe must not be pulled back during and after the application.

4 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclar.com, and your responsible competent authority.
- The current Instructions for Use are available in the download section of the Ivoclar Vivadent AG website (www.ivoclar.com).

Warnings

- Observe the Safety Data Sheet (SDS) (available in the download section of the Ivoclar Vivadent AG website www.ivoclar.com).
- Unpolymerized material should not come in contact with skin/mucous membrane or eyes. In an uncured state, the product may have a mild irritating effect and, in rare cases, may lead to a sensitization against methacrylates.
- Commercial medical gloves do not provide protection against the sensitizing effects of methacrylate.

Disposal information

Remaining stocks must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks. Some of these risks are listed below:

- Loss of retention may lead to accidental swallowing or inhalation of material and dental re-treatment.

4 Shelf life and storage

- Storage temperature 2 - 28 °C (36 - 82 °F) 
- Close syringe immediately after use. Exposure to light may lead to premature polymerization.
- Do not use the product after the indicated date of expiration.
- Date of expiration: see note on Cavifils, syringes and packages.
- Before use, visually inspect the packaging and the product for damage.

5 Additional information

Keep material out of the reach of children!

Not all products are available in all countries.

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the product for its suitability and use for any purpose not explicitly stated in the Instructions.

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