VOCO
AlignerFlow® LC

EN Instructions for use
MD EU Medical device

Product description:
AlignerFlow LC is a flowable, light-curing and fluorescent nanohybrid composite for aligner attachments within the scope of aligner treatments. AlignerFlow LC contains 83% by weight inorganic fillers and is used with a dentine enamel bond.

Shades:
A1, A2, A3

Indications:
- Fabrication of aligner attachments
- Retention of aligners
- Bonding of lingual retainers
- Occlusal build-ups

Contraindications:
AlignerFlow LC contains methacrylate and BHT. AlignerFlow LC should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:
AlignerFlow LC is suitable for use in all patients without any age or gender restrictions.

Performance features:
The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:
AlignerFlow LC should only be used by a professionally trained dental practitioner.

Use:
Preparation:
Clean the teeth being treated. Bring material to room temperature before use. Ensure that the work area is sufficiently dry.

Bonding materials:
AlignerFlow LC is used in the adhesive technique with a dentine enamel bond. Any light-curing bonding material can be used. Follow the respective instructions for use with regard to preparation (etching technique) and application. Enamel which is unprepared must always be etched at the beginning of the bonding process.
For the application of AlignerFlow LC on glass ceramic, zirconium dioxide, aluminium oxide ceramic or metal and any associated use of a primer or similar, please follow the instructions for use of the applied adhesive.

Application:
Tips with two different diameters are supplied with the material. The type 44 tip is particularly suitable for the application of AlignerFlow LC into the aligner template for the fabrication of attachments. The type 46 tip with the larger diameter is recommended for the bonding of retainer wires.
Lock the enclosed application tip on the syringe with a clockwise turn. Do not pull the plunger back on the syringe during or after use to guarantee the function of the non-dripping syringe.
The attached application tip is intended for single use only. Remove it after use and seal the syringe again tightly with its original cap. The use of a protective barrier is recommended to protect the syringe against contamination with bodily fluids, dirt from hands or oral tissue.

a) Fabrication of aligner attachments / Retention of aligners
Apply AlignerFlow LC directly into corresponding recesses of the attachment template. Insert the attachment template into the mouth and cure each attachment through the transparent template.
Tip: Properties of the attachment template, e.g. material composition, surface quality, 3D print quality, can affect the processing of AlignerFlow LC. To test the expected compatibility with a template material, e.g. the easy removability of the attachment from the template, we recommend testing extraorally.
Following removal of the attachment template, finish the attachments in accordance with the generally applicable standards if and as necessary.
Once the aligner treatment is completed, the attachments can be removed using a grinding/finishing diamond burr. Subsequent polishing of the tooth surface is then recommended.

b) Attachment of lingual retainers
Place the prefabricated retainer wire in the desired position on the tooth surfaces, fixing in place temporarily if necessary, and apply AlignerFlow LC to the wire and surrounding tooth surfaces. Cover the ends of the wire with AlignerFlow LC too. Cure each applied portion of the material separately.

c) Occlusal build-ups
In principle, the preparation of dental hard tissue should follow the rules for adhesive restorative treatment; preparation should be minimally invasive in order to protect healthy dental hard tissue. Apply AlignerFlow LC in layers of no more than 2 mm thickness and then light-cure.

Light-curing:
Conventional polymerisation devices are suitable for light-curing the material. Depending on the light output, polymerisation times are as follows:

Table with 2 columns: Light intensity, Time. Rows: LED/halogen lamp, ≥ 1000 mW/cm², 500 - 1000 mW/cm².

Bring the light aperture of the lamp as close as possible to the surface of the attachment template or, in the case of retainers and occlusal build-ups, to the surface of the material, as otherwise the curing will be less effective.
Inadequate curing can cause discolouration and complaints.

Warnings, precautionary measures:
- UV light can cause damage. Patients and dental personnel should always take note of the corresponding safety measures. This includes wearing appropriate protective UV goggles, among other things. Do not look into the light source and do not aim the light directly into the eyes. Only use the light at the treatment site.
- Phenolic substances, especially products containing eugenol and thymol, lead to curing failures in restorative composites. Consequently, the use of zinc oxide/eugenol cements or other materials containing eugenol should be avoided in combination with restorative composites.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Constituents (in descending order according to content):
Barium aluminium borosilicate glass, silicon dioxide, DDDMA, BisGMA, TEGDMA, BisEMA, fumed silica, initiators, stabilisers, pigments

Storage instructions and application method:
Storage at 4°C-28°C. Seal the syringes immediately after removing the material in order to avoid exposure to light and subsequent polymerisation. Do not use after the expiry date.

Disposal:
Dispose of the product in accordance with local regulations.

Reporting obligation:
Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of AlignerFlow LC must be reported to VOCO GmbH and the responsible authority.

Note:
The Summary of Safety and Clinical Performance of AlignerFlow LC can be found in the European database on medical devices (EUDAMED - https://ec.europa.eu/tools/eudamed). Detailed information can also be found at www.voco.dental.

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VC 60 AA1526 E1 1024 99 © by VOCO