IPS Empress® Direct Opaque

[en] Instructions for Use

Opaquer

Rx ONLY

C € 0123

Date information prepared:
2024-07-04 / Rev. 0
web

Manufacturer:
Ivoclar Vivadent AG
Benderstrasse 2
9494 Schaan/Liechtenstein
www.ivoclar.com

ivoclar

English

1 Intended use

Intended purpose

Esthetic characterization of direct and indirect anterior and posterior restorations and masking of discolorations.

Patient target group

Patients with permanent teeth.

Intended users

Dentists

Special training

No special training required.

Use

For dental use only.

Description

IPS Empress Direct Opaque is a light-curing opaquer for use in dental restorative procedures.
IPS Empress Direct Opaque cures with light in the wavelength range of 400–500 nm.

Indications

- None

Types of restorations:

Masking of

- exposed metal surfaces when repairing defective ceramic and composite veneers intraorally;
- discoloration on teeth.

Contraindications

The use of the products is contraindicated if the patient is known to be allergic to any of their ingredients.

Limitations of use

- If a dry working field cannot be established.
- If the stipulated working technique cannot be applied.
- When used as a restorative without a capping layer.
- Apply the product at ambient temperature. Cool temperatures render the material difficult to extrude.
- The use of different application cannulas can render the material difficult to extrude.
- The syringe may be used for multiple applications, but should be protected by a single-use hygienic sleeve. Use a suitable hygienic protective sleeve for the syringe as indicated by the manufacturer.
-  The application cannula of the syringe must not be used for more than one patient due to hygienic reasons (prevention of cross-contamination between patients).

Side effects

In rare cases, components of IPS Empress Direct Opaque may lead to sensitization. The products must not be used in such cases.

In order to avoid irritation of the pulp, provide areas close to the pulp with suitable pulp/dentin protection. Selectively apply a calcium hydroxide-containing material in areas close to the pulp and cover with a suitable cavity liner.

Interactions

- Phenolic substances such as eugenol/clove oil inhibit the polymerization of methacrylate-based materials. Consequently, the application of such materials together with IPS Empress Direct Opaque must be avoided.
- Discoloration may occur in combination with cationic mouthwashes, plaque disclosing agents and chlorhexidine.
- Disinfectants with an oxidative effect (e.g. hydrogen peroxide) may interact with the initiator system, which in turn may impair the curing process.

Clinical benefit

- Restoration of esthetics.

Composition

Barium glass, Bis-GMA, UDMA, D3MA, ytterbium trifluoride, silicone oxide, Si-Zr mixed oxide, Ba-Al fluorosilicate glass, titanium oxide
Total content of inorganic fillers: 23 vol%.

Particle size of inorganic fillers: between 0.03 µm and 16.3 µm.

2 Application

Shade determination

The product is available in the universal shade Opaque.

Isolation

Adequate relative or absolute isolation is required.

A) Application for the repair of indirect restorations

Preparatory steps

Proceed as follows:

1. Roughen the metal surface (using an intraoral air abrasion unit or diamond or corundum grinding instruments).
2. Bevel the margins of the defective ceramic or composite surface (1–2 mm) under water cooling using a diamond grinding instrument. Dry with oil-free air.
3. Condition the surface to be repaired. Apply a thin layer of bonding agent/primer to the pre-treated surfaces with a brush or microbrush, according to the instructions for use of the manufacturer of the product used.
4. Apply the adhesive to the entire surface to be repaired (metal ceramic/composite resin) according to the instructions for use of the product used.

Application of the product:

Proceed as follows:

1. Apply the product in a layer of max. 0.5 mm thickness exclusively to the prepared metal surface.
2. Light-cure the product according to the indicated exposure time and light intensity in Table 1.
 - Hold the light emission window of the curing light as close to the material as possible. If the surface is still not sufficiently masked, apply a second layer of the product and light-cure again.
 - If the diameter of the light guide tip is smaller than the restoration diameter, perform an overlapping polymerization in order to ensure that all areas of the restoration are covered.

NOTICE! The instructions for use of the curing light must be observed.

3. If the restoration includes proximal areas, use a matrix and secure it with wedges.
4. Repair the restoration by applying increments of the composite material and light-cure the material according to the instructions for use of the manufacturer of the respective product.
5. Finish and polish the restoration.

B) Application for the masking of tooth discolorations

Proceed as follows:

1. Condition the surface and apply the adhesive according to the instructions for use of the product used.
2. Apply the product in very thin layers of max. 0.5 mm with a brush or with another suitable instrument.
3. Light-cure the product according to the indicated exposure time and light intensity in Table 1.

- Hold the light emission window of the curing light as close to the material as possible. If the surface is still not sufficiently masked, apply a second layer of the product and light-cure again.
- If the diameter of the light guide tip is smaller than the restoration diameter, perform an overlapping polymerization in order to ensure that all areas of the restoration are covered.

NOTICE! The instructions for use of the curing light must be observed.

4. Coat with a layer of suitable composite restorative or build-up material according to the instructions of the manufacturer of the product used.
5. Finish and polish the restoration or prepare for indirect restoration.

Application notes

- Only cannulas intended for the product must be used. In case other cannulas are used, there is a risk that the cannula may suddenly become loose during application.
- Do not expose the product to intensive light while applying it, as this might lead to premature polymerization.
- Do not disinfect syringes with oxidizing disinfection agents.
- Use a suitable hygienic protective sleeve for the syringe as indicated by the manufacturer.
- Avoid any contact of the reusable syringe not covered with a protective sleeve with the mouth of the patient.
- Dispose of contaminated syringes.

3 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclar.com, and your responsible competent authority.
- The current Instructions for Use are available on the website (www.ivoclar.com/eIFU).
- Explanation of the symbols: www.ivoclar.com/eIFU
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>.
- Basic UDI-DI: 76152082AFILL015JM

Warnings

- Observe the Safety Data Sheet (SDS) (available at www.ivoclar.com).
- Avoid contact of the uncured product with the skin/mucous membrane or eyes. The unpolymerized product may have a slight irritating effect and may lead to a sensitization against methacrylates. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Disposal information

Remaining stocks or removed restorations must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following clinical residual risks are known:

- Failure of the adhesive bond (loss of the filling)
- Postoperative sensitivity
- Inclusion of air bubbles during the placement of the filling
- Chipping, fractures
- Ingestion of material

4 Shelf life and storage

- Storage temperature 2–28 °C (36–82 °F)
- Close syringe immediately after usage. Exposure to light leads to premature polymerization.
- Do not use the product after the indicated date of expiration.
- Expiry date: See information on syringes and packages.
- Before use, visually inspect the packaging and the product for damage. In case of any doubts, please contact Ivoclar Vivadent AG or your local dealer.

5 Additional information

Keep material out of the reach of children!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the materials for their suitability and use for any purpose not explicitly stated in the Instructions.

Table 1

Light intensity	Exposure time
500–900 mW/cm ²	40 s
1,000–1,300 mW/cm ²	20 s

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