EN

CORECEM® FOR DENTAL USE ONLY

DESCRIPTION

CORECEM® "shade A1" and **CORECEM ILLUSION®** are dual-cure flowable composite resins used for post cementation and core build-up during an endodontic treatment or during a restoration procedure. The term **CORECEM®** will be used in the rest of the document. Both formulations have identical mechanical properties with the only difference being the pigment used for coloration and for the Illusion Technology that enables the composite to change color.

CORECEM® consists of a light curable resin (called a base composite) and a catalyst, supplied in two barrels of an auto-mixing syringe (see figure 1). When the base and catalyst are mixed through the static mixing tip, the mixture begins the process of self-curing (also known as self-curable or auto-polymerization). It can also be cured at any time using a curing light.

CORECEM® is supplied in auto-mixing 1:1 MINIMIX (5 ml) or AUTOMIX (25 ml) syringes, provided with root canal tip for the mixing of the base and the catalyst.

COMPOSITION

CORECEM® is composed of a light curable resin (called a base composite) and a catalyst.

CORECEM® ILLUSION® is a special formulation that is tooth colored at body temperature but turns blue when sprayed with cold water.

INTENDED USE

CORECEM® is used as a core build-up material to replace missing tooth structure above the gingiva. It is also used as a cement to bond endodontic fiber posts inside the tooth where root canals treatment has been accomplished. The fiber post acts as an anchor (retention) to retain the core build-up where large amounts of tooth are missing

These products should only be used in hospitals, clinics or dental office with qualified dental personnel (dentists).

INDICATIONS FOR USE

CORECEM® is indicated in case of missing substantial amounts of original tooth structure. It can be used on root canal treated teeth (non-vital), in conjunction with a fiber post, or it can be used on vital teeth that do not need root canal treatment.

CONTRAINDICATIONS

- Do not use if coronal tooth structure is less than 1.5 mm.
- Do not use on patient with known allergy to methacrylate.

WARNINGS

None known to date.

PRECAUTIONS

- CORECEM®** is not recommended to be used for people allergic to methacrylic resins. Additionally, there is no data for the use of **CORECEM®** for children, pregnant or breast-feeding women. Thus, it is not recommended for such use.
- The product is solid non-sterile. It cannot be sterilized.
- Do not use any resin to adjust viscosity of composite restorative material.
- Do not use in combination with substances containing eugenol because eugenol can impair the hardening of the composite and causes discoloration. Do neither store the composite material in proximity of eugenol containing products, nor let the composite allow to come into contact with materials containing eugenol.
- Store unopened material at a dark and cold place. Temperature should not exceed 20 °C (68 °F). Opened cartridges should be used up within 3 months.
- The shelf-life for **CORECEM** Minimix 9 g is 24 months (two years) and for **CORECEM** Automix 42 g is 18 months. Do not use after expiration date.
- Store used **CORECEM®** syringe in the dark with used mixing cannula in place.
- Do not pull back on the plunger after dispensing. This may cause constriction of the syringe.
- Contact of resin pastes with skin and gingival tissue should be avoided, especially by anyone having known resin allergies to methacrylates.
- Use of a rubber dam is recommended.

ADVERSE REACTIONS

Possible allergic reaction to methacrylate.

STEP-BY-STEP INSTRUCTIONS

Preparing the mixing syringes

- Remove the cap of the **CORECEM®** syringe and throw it away (do not use it again!). It is replaced by a special

1:1 mixing cannula. Turn the cannula 90° until it locks in position.

- Discard the first 2-3 mm (about the size of a pea) of the extruded material. This must be done for each new mix.
- Do not pull back on the plunger after dispensing. This may cause constriction of the syringe.
- Working time in the self-cured mode (23 °C (74 °F)): 1:30 minutes from start of mixing.
- Setting time in the self-cured mode (37 °C (99 °F)): 5:30 minutes from start of mixing.

Post cementation
- Prepare the post space according to manufacturer's instructions. Etching and Primer are required.
- Dispense a small amount 2-3 mm to eliminate the previous mixing and voids.
- Apply **CORECEM®** inside the root canal using the root canal tip (recommended). Filling should occur from bottom upward to prevent voids.
- Seat the post, light-cure for at least 20 s. (400 mW/cm²).
- Core build-up**
- Fill in **CORECEM®** into a plastic form or fill in directly into the preparation. Filling should occur from bottom upward to prevent voids.
- Light-cure for at least 40 s 400 mw/cm² or allow the system to self-cure for at least 3 minutes.
- Remove the matrix when the material is completely set.
- Start the final crown preparation.
- For **CORECEM® ILLUSION** use cold water spray during trimming. After finishing, blue color will become dentin shade like A1.

STORAGE

- Store unopened material at a dark and cold place (< 20 °C – 68 °F)
- The shelf-life for **CORECEM®** Minimix 9 g is 24 months (two years) and for **CORECEM®** Automix 42 g is 18 months. Do not use after expiration date
- Opened **CORECEM®** can be stored along with its syringe and mixing cannula in the dark. Temperature should not exceed 20 °C (68 °F). Opened cartridges should be used up within 3 months
- Do not store the composite material in proximity of eugenol containing products.

DISPOSAL

Dispose in a safe manner in accordance with local/national regulations.

SAFETY AND LIABILITY

WARRANTY

RTD warrants this product will be free from defects in material and manufacture. RTD makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and RTD's sole obligation shall be repair or replacement of the RTD product.

LIMITATION OF LIABILITY

Except where prohibited by law, RTD will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

For dental use only.

**** Refer to the manufacturer's instructions for use.***

***** Please refer to your autoclave manufacturer for specific sterilization instructions if needed.***

TERMS AND SYMBOLS IDENTIFICATION

	Manufacturer
	Batch code
	Catalogue number
	Number of unit / Content
	Upper limit of temperature
	Use-by date
	Keep away from sunlight

	MANUFACTURER:
	RTD
	3 rue Louis Néel, Technoparc Espace Gavanière
	38120 St Egrève - France
	RTD-info@dentsplysirona.com; www.rtddental.com

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