VOCO

Aqua Ionofil Plus

0482

EN

Instructions for use

MD EU Medical device

In accordance with DIN EN ISO 9917-1

Product description:
Aqua Ionofil Plus is a radiopaque, water-mixable glass ionomer filling cement with chemical adhesion to enamel and dentine. **Aqua Ionofil Plus** is available in shade A3 and has enamel-like translucency.

Indications:
– Restoration of deciduous teeth (especially class I)
– Repair of V-shaped defects and enamel as well as root erosions in the dental neck area including class V fillings
– Fillings for cavities caused by carious lesions
– Fillings of class III and small fillings of class I (incl. extended fissure sealing) as well as core build-up under crowns and as a cavity lining

Contraindications:
Aqua Ionofil Plus contains polyacrylics. **Aqua Ionofil Plus** should therefore not be used in patients with a known hypersensitivity (allergy) to this ingredient.

Patient target group:
Aqua Ionofil Plus is suitable for use in all patients without any age or gender restrictions.

Performance features:
The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:
Aqua Ionofil Plus should only be used by a professionally trained dental practitioner.

Use:
Preparation:
Prepare the cavity using a minimally invasive approach (removing only carious tooth substance). The layer thickness of the restorations must not be less than 1 mm. Avoid thinly tapering margins. Clean the cavity, rinse with water and dab or blow dry. Do not overdry the dentine: it should have a moist, shiny appearance. If the residual dentine thickness is insufficient, the area close to the pulp should be covered with a suitable calcium hydroxide product.

Attention! Mix Aqua Ionofil Plus only with water.

Shake powder well before use.

Allow **Aqua Ionofil Plus** to reach room temperature (max. 23 °C) before application. In case of higher temperatures it is recommended to use a cooled glass plate.
Dispense **Aqua Ionofil Plus** powder onto a glass plate or suitable mixing pad and mix the powder portion by portion with the water. Use 1 drop of water for 1 level measuring scoop of powder (mixing ratio powder : water = 5,4 - 6,3 g : 1).
When dispensing the water, hold the dropper bottle in an absolutely vertical position, making sure that the metal cannula protrudes approx. 2 - 3 mm from the dropper. Do not use drops for mixing which escape when tilting the bottle.
Use a strong plastic spatula or another abrasion-resistant spatula for mixing.

Mixing time	30 - 40 s
Working time (15 - 23°C)	2 - 3 min.
Setting time (37°C)	2 - 3 min.
Total setting time	5 - 6 min.

The mixing time is approx. 30 - 40 s, but should not exceed 1 min. Working time from end of mixing procedure is approx. 2 - 3 min (at room temperature). After the filling has been inserted into the cavity, the setting time of **Aqua Ionofil Plus** is approx. 2 - 3 min.
Ensure that the cavity remains dry during the filling procedure.
During the setting phase, especially immediately after application, glass ionomer cements are sensitive to humidity (breath, saliva) and also to drying-out.
Therefore, if possible, cover immediately with a strip or matrix or apply a protective varnish (see instructions for use **Final Varnish LC**).
The rough excess material can be removed with a sharp instrument after **Aqua Ionofil Plus** has set (no earlier than 4 - 6 min after application depending on mixing temperature and mixing ratio). Finish with an ultra-fine-grained diamond or a flexible disk.
Afterwards, apply a coating of varnish.

Warnings, precautionary measures:
– Avoid contact with mucous membranes and eyes.
– Please also do not use any insulating material for the dentine in order to preserve full adhesion (with the exception of pulp protection).
– Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Constituents (in descending order according to content):
Fluorosilicate glass, polyacrylic acid, water, tartaric acid, colour pigments

Storage instructions and application method:
Storage at 4-23°C. Reclose powder and liquid containers well immediately after use. Do not use after the expiry date.

Disposal:
Dispose of the product in accordance with local regulations.

Reporting obligation:
Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **Aqua Ionofil Plus** must be reported to VOCO GmbH and the responsible authority.

Note:
The Summary of Safety and Clinical Performance of **Aqua Ionofil Plus** can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>). Detailed information can also be found at www.voco.dental.

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VC 60 AA1509 E1 0724 99 © by VOCO