VOCO

Ionoseal®

EN Instructions for use

MD EU Medical device

Product description:
Ionoseal is a light-curing, glass ionomer composite cement for linings, extended fissure sealing and care of smaller lesions.

- Indications:**
- Linings under amalgams and composites
 - Extended fissure sealing
 - Care of smaller lesions

Contraindications:
Ionoseal contains methacrylate and BHT. **Ionoseal** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:
Ionoseal is suitable for use in all patients without any age or gender restrictions.

Performance features:
The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:
Ionoseal should only be used by a professionally trained dental practitioner.

Use:

1. Apply the preparation in accordance with the rules of applied restorative treatment. Make sure that the tooth surfaces are perfectly dry and avoid contamination!
2. Cover areas in the vicinity of the pulp with a calcium hydroxide preparation. It is recommended that an adhesive bond be achieved (especially in the case of cavities without enamel limitation on all sides).
3. **Ionoseal tube:** The required amount of **Ionoseal** can be applied directly from the tube onto the application instrument. Avoid contamination, i.e. disinfect after use. In the case of linings with a thickness greater than 1 mm, apply in layers and cure.
Syringe: In order to guarantee that the syringe does not drip or run, it is important to ensure that the plunger is not retracted during or after use. The attached application tip is intended for single use only. Remove it after use and seal the syringe again tightly with its original cap. The use of a protective barrier is recommended to protect the syringe against contamination with bodily fluids, dirt from hands or oral tissue.
4. Conventional polymerisation devices are suitable for light-curing the material. The light output should be a minimum of 500 mW/cm² for halogen polymerisation devices and 300 mW/cm² for LED devices. Bring the light source as close as possible to the surface (approx. 2 mm) and cure for at least 20 s. At greater distances, increase the curing times.
5. Finish the **Ionoseal** surface with rotary instruments and place the fillings as per the manufacturer's specifications.

Warnings, precautionary measures:

- **Ionoseal** adheres to materials, e.g. crown and bridge materials, temporary methacrylate-based sealing materials, and certain impression materials (e.g. polyether), which can result in damage to or dislodgement of the base fill. In these cases, it is recommended that a bonding agent be applied.
- Due to the high stability of light-cured **Ionoseal**, it is not possible to remove excess material with a probe as the entire base fill may become dislodged (use rotary instrument).
- Always allow pulp protection agents or dentine isolating agents to set or dry prior to the use of **Ionoseal** (see the product's instructions for use), in order not to reduce the bonding adhesion to dentine.
- Once dispensed, **Ionoseal** may not be exposed to direct light (surgical lamps, sunlight) and must be applied as quickly as possible.
- Phenolic substances, especially products containing eugenol and thymol, interfere with the curing of filling composites. The use of eugenol-containing materials in combination with **Ionoseal** should thus be avoided.
- Curing times that are too short or an inadequate light intensity will also result in incomplete curing. Thus, the exposure lamp and light guide must be inspected regularly.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for the intended purpose.

Constituents (in descending order according to content):

Fluorosilicate glass, barium sulphate, BisGMA, UDMA, HEDMA, pyrogenic silica, TEGDMA, initiators, stabilisers, dyes

Storage:

Store **Ionoseal** at 4°C-23°C away from light and tightly sealed. The material must be warmed up to room temperature prior to use. Do not use after the expiry date.

The vertical storage of opened tubes in the stand tray (removable) is especially practical. It allows for **Ionoseal** to always be immediately available.

Disposal:

Dispose of the product in accordance with local regulations.

Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition, and a serious risk to public health that arise or could have arisen in association with the use of **Ionoseal**, must be reported to VOCO GmbH and to the responsible authority.

Note:

The Summary of Safety and Clinical Performance of **Ionoseal** can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>). Detailed information can also be found at www.voco.dental.

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VC 60 AA1326 E1 0623 99 © by VOCO