VOCO

Ionolux®

Handmix

EN **Instructions for use**
MD **EU Medical device**
In accordance with DIN EN ISO 9917-2

Product description:

Ionolux is a resin-reinforced light-curing radiopaque glass ionomer restorative material. The material has outstanding aesthetics, good handling properties and high chemical adhesion.

Shades:
A1, A2, A3

Indications:

- Class III and V restorations, primarily cervical fillings and root caries
- Restoration of deciduous teeth
- Small class I fillings
- Temporary fillings
- Core build-ups
- Linings

Contraindications:

Ionolux contains polyacrylic acid, BHT and methacrylates. **Ionolux** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:

Ionolux is suitable for use in all patients without any age or gender restrictions.

Performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:

Ionolux should only be used by a professionally trained dental practitioner.

Use:

Preparation:

In principle, cavity preparation should be minimally invasive, to conserve healthy tooth substance. There is no need for undercuts. Class V cavities should have a minimum depth of 1 mm. Remove preparation residues with water spray, then dry with air so that the surface remains slightly glossy (do not overdry). Place a matrix, if necessary. If the residual dentine thickness is insufficient, the area close to the pulp should be covered with a calcium hydroxide product.

Shade selection:

With the help of a VITA® shade guide on a moist tooth.

Mixing:

Shake the powder and the liquid bottle thoroughly for at least 3 s prior to dispensing. Dispense a level spoonful of powder onto a suitable mixing pad. Close the bottle immediately after dispensing to protect the powder against moisture. Dose 2 drops of liquid per level scoop of powder (powder/liquid ratio 3.2 : 1 g/g). Hold the bottle vertically when dispensing the drops. The dropper should be free from dried residues. Close the liquid bottle immediately after use. Work the first half of the powder into the liquid, then the second; for this, use spatulas of solid plastic or other abrasion-proof materials. Spread out the paste several times to achieve a homogeneous consistency.

Process the cement at room temperature (15 °C - 23 °C).

Application:

Apply the mixed **Ionolux** into the cavity with a suitable instrument and contour. Then light-cure with a suitable polymerisation lamp. Fillings of more than 2 mm thickness should be applied and light-cured in layers.

Light-curing:

The material can be light-cured with market-available polymerisation lamps. Their output should not be lower than 500 mW/cm². The curing time is 20 s per layer. Hold the tip of the light-curing lamp as close as possible to the surface of the filling. If the distance is more than 5 mm, the curing depth will be smaller.

Powder/liquid ratio [g/g]	3.2 : 1.0
Dosage [measuring spoon/drop]	1 : 2
Mixing time [s]	30
Working time [min]	>3:00
Light-curing [s]	20

Finishing:

After light-curing, remove excess material and finish and polish the filling as usual.

Warnings, precautionary measures:

- In case of contact with oral mucosa or skin, remove immediately with a cotton swab soaked in alcohol. Rinse with water.
- In case of eye contact, rinse well with plenty of water and consult an ophthalmologist.
- Phenolic substances, especially preparations containing eugenol or thymol, inhibit the curing of **Ionolux**. Therefore, do not use zinc oxide eugenol cements or other eugenol-containing substances in combination with **Ionolux**.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Constituents (in descending order according to content):

Fluorosilicate glass, polyacrylic acid, HEMA, water, UDMA, GlyDMA, initiators, stabilisers, pigments

Storage instructions and application method:

Storage at 4°C - 23°C. Tightly close powder and liquid container immediately after use. Store dry and in an area protected from light. Do not use after the expiry date.

Disposal:

Dispose of the product in accordance with local regulations.

Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **Ionolux** must be reported to VOCO GmbH and the responsible authority.

Note:

The Summary of Safety and Clinical Performance of **Ionolux** can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>). Detailed information can also be found at www.voco.dental.

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VC 60 BB1990 E1 0822 99 © by VOCO
VITA® is a registered trademark of VITA Zahnfabrik H. Rauter GmbH & Co.