# Evetric® Bond

## [en] Instructions for Use

Light-curing, single-component dental adhesive

Rx ONLY

C € 0123

Date information prepared:  
2024-10-08 / Rev.1

Manufacturer:  
Ivoclar Vivadent AG  
Bendererstrasse 2  
9494 Schaan/Liechtenstein  
[www.ivoclar.com](http://www.ivoclar.com)

ivoclar

English

## 1 Intended use

### Intended purpose

Bonding of dental restorations to dental hard tissue, sealing of dentin lesions.

### Patient target group

- Patients with permanent teeth
- Patients with deciduous teeth

### Intended users / Special training

- Dentists
- No special training required

### Application

For dental use only.

### Description

Evetric® Bond is a light-curing, self-etching single-component dental adhesive for enamel and dentin. Evetric Bond can be light-cured for 10 s at a light intensity of 500 to 1,300 mW/cm<sup>2</sup>.

### Indications

Missing tooth structure in anterior and posterior teeth

### Areas of application:

- Direct-placed light-curing composite restorations in permanent teeth
- Direct-placed core build-ups with light-, self- and dual-curing composites in permanent teeth.
- Repairs of fractured composite restorations in permanent teeth.
- Desensitization of hypersensitive cervical areas in permanent teeth.
- Restorations of deciduous teeth

The limitations of use must be observed.

### Contraindications

- The use of the product is contraindicated if the patient is known to be allergic to any of its ingredients
- Direct pulp capping

### Limitations of use

- Do not use Evetric Bond as a primer for ceramic restorative materials. A suitable ceramic primer must be used.
- Do not use Evetric Bond for applications in which sufficient illumination cannot be ensured (e.g. cementation of root canal posts).
- Evetric Bond must not be used if the stipulated working technique cannot be employed.
- The disposable applicators must not be reused.
- Light intensities > 1,300 mW/cm<sup>2</sup> must not be used for deciduous teeth.

## Side effects

In rare cases, components of Evetric Bond may lead to sensitization. In such cases, Evetric Bond should no longer be used.

## Interactions

- Peroxide-based bleaching agents and disinfectants may inhibit the polymerization of adhesives. Such materials must not be used ≤ 14 days prior to the application of Evetric Bond.
- Phenolic substances (e.g. eugenol) may inhibit the polymerization of Evetric Bond.
- Before any treatment with Evetric Bond, completely rinse off all iron- and aluminium-based styptic agents with water.
- Fluoridation agents must not be used prior to the application of Evetric Bond.

## Clinical benefit

- Reconstruction of chewing function
- Restoration of esthetics

## Composition

HEMA, Bis-GMA, ethanol, water, D3MA, MDP, methacrylate-modified polyacrylic acid, silicon dioxide, camphor quinone, ethyl p-dimethyl aminobenzoate, 2-dimethyl aminoethyl methacrylate.

## 2 Application

### 2.1 Isolation

Adequate relative or absolute isolation using suitable auxiliaries is required.

### 2.2 Pulp protection / cavity liner

In very deep cavities, areas close to the pulp must be selectively coated with a calcium hydroxide liner and subsequently covered with a pressure-resistant cement (e.g. a glass ionomer cement).

### 2.3 Application of the adhesive

- Dispense the desired amount of Evetric Bond into a mixing well and apply it using a disposable applicator.
- Protect the adhesive from light.
- Carefully close the bottle immediately after each use.
- Starting with the enamel, completely coat the tooth surfaces to be treated with Evetric Bond.
- The adhesive must be scrubbed into the tooth surface for at least 20 s. This time must not be shortened. Applying the adhesive on the tooth surface without scrubbing is inadequate.
- Disperse Evetric Bond with oil- and moisture-free compressed air until a glossy, immobile film layer results.
-  A new disposable applicator must be used for each application.

### Application notes

Do not use the bottle intraorally. For hygienic reasons, suitable application aids should be used.

### Important information

- Exposure to intensive light should be avoided during application.
- Prevent any contamination of the applicator with blood, saliva, sulcus fluid or water during application. In the event of contamination, the cavity must be rinsed again and the procedure must be repeated with a new applicator.

### 2.4 Light-curing the adhesive

Evetric Bond can be light-cured for 10 s at a light intensity of 500 to 1,300 mW/cm<sup>2</sup>.

If the diameter of the light guide tip is smaller than the restoration diameter, perform an overlapping polymerization in order to ensure that all areas of the restoration are covered.

**The instructions for use of the curing light must be observed.**

### Safety note

Avoid direct exposure of the gingiva, mucous membrane or skin to the light emitted by the curing light.

### 2.5 Application of the restorative composite

Proceed according to the instructions for use of the restorative composite used.

## **2.6 Procedure for intraoral repairs of fractured composite fillings/ composite restorations with light-curing materials**

- Roughen the surface to be repaired using diamond finishers and thoroughly clean with water spray and dry with oil- and water-free compressed air.
- Proceed according to the direct application for Evetric Bond.

## **3 Safety information**

- Evetric Bond is irritant. Avoid contact with the skin, mucous membrane and eyes.
- After skin contact, immediately rinse with copious amounts of water.
- If the material should come in contact with the eyes, immediately rinse with copious amounts of water and consult a physician/ ophthalmologist.
- In rare cases, contact with the skin may lead to sensitization towards the ingredients.
- Commercial medical gloves do not provide protection against the sensitizing effects of methacrylate.
- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/ Liechtenstein, website: [www.ivoclar.com](http://www.ivoclar.com), and your responsible competent authority.
- The current Instructions for Use are available in the download section of the Ivoclar Vivadent AG website ([www.ivoclar.com](http://www.ivoclar.com)).
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>.
- Basic UDI-DI: 76152082ABOND005HS

### **Warnings**

Observe the Safety Data Sheet (SDS) (available in the download section of the Ivoclar Vivadent AG website [www.ivoclar.com](http://www.ivoclar.com)).

### **Disposal information**

Remaining stocks or removed restorations must be disposed of according to the corresponding national legal requirements.

### **Residual risks**

Users should be aware that any dental intervention in the oral cavity involves certain risks. Some of these risks are listed below:

- Failure of the adhesive bond (loss of the filling)
- Postoperative sensitivity

## **4 Shelf life and storage**

- Storage temperature 2–28 °C (36 – 82 °F)
- If the indicated storage temperature is exceeded, proper functioning of the product cannot be ensured.
- Do not use the product after the indicated date of expiration.
- Expiration date: See note on packages.

Before use, visually inspect the packaging and the product for damage. In case of any doubts, please contact Ivoclar Vivadent AG or your local dealer.

## **5 Additional information**

Keep material out of the reach of children!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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