**RelyX™ U200****Self-adhesive Resin Cement****Самоадгезивни Композиционни Цемент****Самодејзивен композитен цимент****Kompozitni cement za adhezivno cementiranje protetiskih nadomjestaka****Önkondicionáló rezincment****Samopregajajuci cement do osadzania prac protetycznych****Ciment râsnic autoadeziv****Instructions for Use****Инструкция по применению****Указания за употреба****Upotrebe za upotrebu****3M Deutschland GmbH****Dental Products****Caris-Strur-1****41453 Neuss - Germany****3M, ESPE, Cojet, Lava, RelyX и Sol-lex  
are trademarks of 3M or 3M Deutschland GmbH.  
© 2014, 3M. All rights reserved.****CE 0123****4407587872****Használati utasítások****Instrukcja užívania****Instrucțiuni de utilizare****Notes****Pre-treatment of Zirconia and Aluminum Oxide Ceramic Restorations****Please follow the instructions for use of the restoration material. If the manufacturer has not provided deviating instructions, we recommend the following procedure:****► Use hydrofluoric acid to etch the glass ceramic restoration surface to be luted.****► Rinse thoroughly with water for 15 seconds and dry with water-free and oil-free air.****► Apply an appropriate silane in accordance with the instructions for use, e.g.: RelyX™ Ceramic Primer: allow to react for 5 sec, then air dry so that the solvent evaporates completely.****Pre-treatment of Maryland and Inlay/Onlay Bridges****Please follow the instructions for use of the restoration material. If the manufacturer has not provided deviating instructions, we recommend the following procedure:****► Blast the restoration surface to be luted with aluminum oxide ≤40 µm.****► Clean the blasted surface with alcohol and dry it with water-free and oil-free air.****► Use hydrofluoric acid to etch the glass ceramic restoration surface to be luted.****► Rinse thoroughly with water for 15 seconds and dry with water-free and oil-free air.****► Apply an appropriate silane in accordance with the instructions for use.****► For Maryland and Inlay/Onlay bridges, please refer to "Pre-treatment of Maryland and Inlay/Onlay bridges".****Pre-treatment of Composite Restorations****Please follow the instructions for use of the restoration material. If the manufacturer has not provided deviating instructions, we recommend the following procedure:****► Use hydrofluoric acid to etch the restoration surface to be luted with aluminum oxide ≤40 µm.****► Rinse thoroughly with water for 15 seconds and dry with water-free and oil-free air.****► Apply an appropriate silane in accordance with the instructions for use.****► For Maryland and Inlay/Onlay bridges, please refer to "Pre-treatment of Maryland and Inlay/Onlay bridges".****Pre-treatment of Glass Fiber-reinforced Posts****Please follow the instructions for use of the post to be used. If the manufacturer has not provided deviating instructions, we recommend the following procedure:****► Clean the restorative surface of the blasting agent with water-free and oil-free air.****► When using RelyX™ Fiber Post application of a silane is not required. Other fiber-reinforced posts should be treated in accordance with the relevant instructions for use.****Pre-treatment of Impant Abutments****Please follow the recommendations for the relevant restorative material (e.g., metal, zirconia).****Pre-treatment of Maryland and Inlay/Onlay Bridges****► Coat (silicate) the restoration surface to be luted with the micro-blasting device Cojet Prep and the blast-coating agent Cojet Sand from a distance of 2-10 mm and vertically to the surface for 15 sec, see instructions for use for Cojet Prep and Cojet Sand.****► Clean the restorative surface of the blasting agent with water-free and oil-free air.****► Use hydrofluoric acid to etch the glass ceramic restoration surface to be luted.****► Rinse thoroughly with water for 15 seconds and dry with water-free and oil-free air.****► Apply an appropriate silane in accordance with the instructions for use.****► For etched glass ceramics, please refer to "Pre-treatment of Etchable Glass Ceramics".****Dosing****When dosing the first time:****► Remove the clicker dispenser from the foil pouch and note the removal date.****► Remove the clicker and the clicker cap.****► Remove a sufficient amount of paste on a mixing pad by pressing down the clicker lever. This leaves the filling quantity of base paste and catalyst pastes in the clicker. Discard the paste which has been dispensed.****Precautionary Measures****For Patients and Dental Personnel****► Base paste: contains substances that may cause an allergic reaction by skin contact in certain individuals. Avoid use of this product in patients with known acrylic and/or peroxide (especially peroxide) allergies.****► If prolonged contact with oral tissue occurs, flush with large amounts of water. If a allergic reaction occurs, seek medical attention as needed, remove the product if necessary and discontinue future use.****For Dental Personnel****► This product contains substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, avoid use in certain individuals. If a skin contact occurs, avoid use until the uncurled product. If skin contact occurs, wash with soap and water.****► The use of protective gloves and a no-touch technique is recommended.****Acrylates may penetrate commonly used gloves. If the product contacts the glove, remove and discard the glove, wash the hands immediately with soap and water and then re-choose.****► If allergic reaction occurs, seek medical attention as needed.****3M SDSs can be obtained from www.mmm.com, or contact your local subsidiary.****Precautions during Processing****► Seal temporary restorations using a euugenol-free product (e.g., with RelyX™ Temp NE). The use of temporary cements containing euugenol can inhibit the polymerization process of RelyX U200 during the luting cementing.****► Do not use hydrogen peroxide ( $H_2O_2$ ) as the restorations are adversely affected by the adhesive strength and setting reaction of RelyX U200.****► To avoid pulpal irritation, cover areas in close proximity to the pulp by applying small amounts of sealing calcium hydroxide material. To avoid pulp increases during the cementing of the restorations, the pulp protection should be applied prior to taking an impression for the final restoration.****Tooth Preparation****Preparation of the Cavity/Tooth Stump****► Prior to final preparation, clean the prepared stump or the cavity thoroughly with pumice slurry, rins with a spray, and lightly air dry in only 2-3 bursts of water-free and oil-free air, or use cotton pellets to dry it off.****► Application:****► Protect the working area from contamination with water, blood, saliva and sulfuric fluid during the application and setting phase.****► Apply the mixed cement evenly to the entire cavity and as appropriate to the adhesive strength and setting reaction of RelyX U200.****► Do not use substances such as desensitizers, disinfectants, astringents, dentin sealants, rinsing solutions containing EDTA, etc., after the final cleaning with pumice slurry and water. Their residues may have a detrimental effect on the bonding strength and setting reaction of the cement.****Application in the Root Canal****► Do not use Lentulo-Spirals to insert the cement in the root canal as this can excessively accelerate setting.****► Squeeze the tip to the point.****► As an alternative, the cement can be inserted into the root canal with the aid of an applicator tip, e.g., Centrix® Needles Tubes (Centrix is a registered trademark of the U.S.-based company Centrix); this will minimize the entrainment of air bubbles.****► Rinse immediately after the use of a rubber dam and keep the cavity free of any contamination during cementation.****► Etch the enamel surface of the cavity with 37% phosphoric acid for 15–20 sec. Then rinse thoroughly with water and dry with water-free and oil-free air. In the case of exposed dentin, make sure to selectively etch the enamel to avoid post-operative sensitivity.****Trial Fitting of the Restoration****► Remove the temporary restoration and thoroughly clean the tooth surface of any residue from the temporary cement.****► Make a trial insertion of the final restoration and check its fit and contact points.****► If a low viscosity silicone is used for the trial fit, the silicone residue must then be removed completely.****► For glass ceramic restorations, do not check the occlusion until the restoration has been cemented. Breakage could occur if done prior to being cemented in place.****► Avoid any contamination whatsoever of the surfaces to be treated during pre-treatment and until the final cementation.****Pre-treatment of Restorations****Pre-treatment of Metal Restorations****Please follow the instructions for use of the restoration material. In the absence of ceramic and composite work, we recommend the following procedure:****► Blast the restoration surface to be luted with aluminum oxide ≤40 µm.****► Clean the blasted surface with alcohol and dry it with water-free and oil-free air.****► For Maryland and Inlay/Onlay bridges, please refer to "Pre-treatment of Maryland and Inlay/Onlay bridges".****Pre-treatment of Etchable Glass Ceramic Restorations****Please follow the instructions for use of the restoration material. If the manufacturer has not provided deviating instructions, we recommend the following procedure:****► Use hydrofluoric acid to etch the glass ceramic restoration surface to be luted.****► Rinse thoroughly with water for 15 seconds and dry with water-free and oil-free air.****► Apply an appropriate silane in accordance with the instructions for use, e.g.: RelyX™ Ceramic Primer: allow to react for 5 sec, then air dry so that the solvent evaporates completely.****Pre-treatment of Ceramic Restorations****Please follow the instructions for use of the restoration material. If the manufacturer has not provided deviating instructions, we recommend the following procedure:****► Use hydrofluoric acid to etch the glass ceramic restoration surface to be luted.****► Rinse thoroughly with water for 15 seconds and dry with water-free and oil-free air.****► Apply an appropriate silane in accordance with the instructions for use.****Pre-treatment of Movable Prototypes and Maryland Bridges****Please follow the instructions for use of the restoration material. If the manufacturer has not provided deviating instructions, we recommend the following procedure:****► Use hydrofluoric acid to etch the glass ceramic restoration surface to be luted.****► Rinse thoroughly with water for 15 seconds and dry with water-free and oil-free air.****► Apply an appropriate silane in accordance with the instructions for use.****Pre-treatment of Composite Restorations****Please follow the instructions for use of the restoration material. If the manufacturer has not provided deviating instructions, we recommend the following procedure:****► Use hydrofluoric acid to etch the glass ceramic restoration surface to be luted.****► Rinse thoroughly with water for 15 seconds and dry with water-free and oil-free air.****► Apply an appropriate silane in accordance with the instructions for use.****Pre-treatment of Metal Restorations****Please follow the instructions for use of the restoration material. In the absence of ceramic and composite work, we recommend the following procedure:****► Blast the restoration surface to be luted with aluminum oxide ≤40 µm.****► Clean the blasted surface with alcohol and dry it with water-free and oil-free air.****► For Maryland and Inlay/Onlay bridges, please refer to "Pre-treatment of Maryland and Inlay/Onlay bridges".****Pre-treatment of Etchable Glass Ceramic Restorations****Please follow the instructions for use of the restoration material. If the manufacturer has not provided deviating instructions, we recommend the following procedure:****► Use hydrofluoric acid to etch the glass ceramic restoration surface to be luted.****► Rinse thoroughly with water for 15 seconds and dry with water-free and oil-free air.****► Apply an appropriate silane in accordance with the instructions for use.****Pre-treatment of Ceramic Restorations****Please follow the instructions for use of the restoration material. If the manufacturer has not provided deviating instructions, we recommend the following procedure:****► Use hydrofluoric acid to etch the glass ceramic restoration surface to be luted.****► Rinse thoroughly with water for 15 seconds and dry with water-free and oil-free air.****► Apply an appropriate silane in accordance with the instructions for use.****Pre-treatment of Movable Prototypes and Maryland Bridges****Please follow the instructions for use of the restoration material. If the manufacturer has not provided deviating instructions, we recommend the following procedure:****► Use hydrofluoric acid to etch the glass ceramic restoration surface to be luted.****► Rinse thoroughly with water for 15 seconds and dry with water-free and oil-free air.****► Apply an appropriate silane in accordance with the instructions for use.****Pre-treatment of Composite Restorations****Please follow the instructions for use of the restoration material. If the manufacturer has not provided deviating instructions, we recommend the following procedure:****► Use hydrofluoric acid to etch the glass ceramic restoration surface to be luted.****► Rinse thoroughly with water for 15 seconds and dry with water-free and oil-free air.****► Apply an appropriate silane in accordance with the instructions for use.****Pre-treatment of Metal Restorations****Please follow the instructions for use of the restoration material. In the absence of ceramic and composite work, we recommend the following procedure:****► Blast the restoration surface to be luted with aluminum oxide ≤40 µm.****► Clean the blasted surface with alcohol and dry it with water-free and oil-free air.****► For Maryland and Inlay/Onlay bridges, please refer to "Pre-treatment of Maryland and Inlay/Onlay bridges".****Pre-treatment of Etchable Glass Ceramic Restorations</**

- (sk) Samoadhézivný živícový cement
- (sl) Samo-adhezívni kompozitní cement
- (cs) Samoadhézivní kompozitní pryskyřičný cement
- (tr) Self-Adhesive Rezin Simon
- (et) Isekinnitruu kompositisidmik
- (lv) Pašīlēmōjs svecu cements
- (lt) Savaiminio sukimimo dervinis cementas fiksavimui
- (uk) Самоадгезивний композитний цемент

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© 2014, 3M All rights reserved.