

Dispenser 1:1/2:1

Dispenser 4:1/10:1

3M Deutschland GmbH
Health Care Business
Carl-Schurz-Str. 1
41453 Neuss – Germany

3M, ESPE, Garant and Protemp are trademarks of 3M or 3M Deutschland GmbH.
Used under license in Canada. © 2021, 3M. All rights reserved.

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310019238/02

ENGLISH

General
These Instructions for Use should be kept for the duration of product use. The product may only be used when the product labeling is clearly readable. For details on all additionally mentioned products please refer to the corresponding instructions for use.

Intended Purpose
Intended purpose: application device for dispensing dental materials from corresponding double chamber syringes.
Intended users: educated dental professionals who have theoretical and practical knowledge on usage of dental products.
Patient target group: all patients requiring dental treatment, unless the patient's condition limits the use.

Contraindications
None

Potential Undesirable Side Effects and Complications
None known.

Use of the Mixing Tip
Please refer to the Instructions for Use of the Impression Material to choose the correct Mixing Tip.

Cleaning and Disinfection of the Garant Dispenser
Cautions
Make sure that the cleaning and disinfection agents you have chosen do not contain any of the following materials:
• Oxidizing agents (e.g., hydrogen peroxides)
• Oils
• Glutaraldehydes

Use only validated methods for cleaning and disinfection.
Use solely cleaning and disinfection agents tested for efficacy and compatibility with the medical device (e.g., FDA certification). CavWipes™ have been tested for use as a cleaning and disinfection agent. Always observe all applicable legal and hygiene regulations for dental offices and/or hospitals.

General
The following procedure should be applied after use or before any subsequent use. The Garant Dispenser must first be cleaned manually and, after cleaning, be manually disinfected in accordance with the following instructions.

Manual Cleaning by Wiping
► Clean the Garant Dispenser after every application with ready-to-use disinfectant wipes (e.g., CavWipes™) based on the active substances alcohols and quaternary compounds or equivalent cleaning agents. Follow the cleaning agent manufacturer's directions for use.
► Continue wiping the dispenser with the wipes until no contamination is visible.
► Special care must be exercised in cleaning the grooves.

Manual Disinfection by Wiping (Intermediate Level Disinfection)
► Disinfect the dispenser with ready-to-use disinfectant wipes (e.g., CavWipes™) based on the active substances alcohols and quaternary compounds or equivalent disinfection agents. Follow the cleaning agent manufacturer's directions for use, observing in particular the contact times.
► Disinfect the Garant Dispenser for three minutes using at least one new disinfectant wipe. Keep all surfaces of the entire Garant Dispenser moist during the three-minute period.
► Special care must be exercised in disinfecting the grooves.

Inspection, Maintenance, Testing
► Examine the Garant Dispenser for damage, discoloration, and contamination before every use.
► Do not use damaged devices under any circumstances. If the Garant Dispenser is not visibly clean, repeat the reprocessing procedure.

Notes
• Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.
• Garant Mixing Tips Purple, Yellow, Green, Blue, White, Garant Intra-oral Tips White, Yellow and Bite are single-use devices and would be damaged when removing set impression material.
• When reattaching Garant Intra-oral Tips, the proper fit of the tip on the mixing tip cannot be ensured.
• Please keep the single-use accessory in the original bag to guarantee the traceability of the batch.

Disposal
Dispose of the contents or container in accordance with the applicable regulations. Please pay special attention to the disposal of contaminated waste to avoid health risks because of improper handling.

Reference Number and Symbol Title	Symbol	Description of Symbol
ISO 15223-1 5.1.1 Manufacturer		Indicates the medical device manufacturer as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.
ISO 15223-1 5.1.3 Date of Manufacture		Indicates the date when the medical device was manufactured.
ISO 15223-1 5.1.4 Use-by date		Indicates the date after which the medical device is not to be used.
ISO 15223-1 5.1.5 Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified.
ISO 15223-1 5.1.6 Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified.
ISO 15223-1 5.4.2 No re-use		Indicates a medical device that is intended for one use or for use on a single patient during a single procedure.

Reference Number and Symbol Title	Symbol	Description of Symbol
ISO 15223-1 5.4.4 Caution		Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
CE Mark		Indicates conformity to European Union Medical Device Regulation or Directive.
Medical Device		Indicates the item is a medical device.
Rx Only		Indicates that U.S. Federal Law restricts this device to sale by or on the order of a dental professional. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1).
Recycle		Indicates packaging is recyclable with mixed paper.
Low-density polyethylene		Indicates plastic component made of low-density polyethylene is recyclable.
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.

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