D.T. LIGHT-POST® ILLUSION® X-RO® FOR DENTAL USE ONLY

DESCRIPTION

D.T. LIGHT-POST ILLUSION® X-RO® are dental posts used for reconstruction of devitalized teeth. They are made of X-RO yellow translucent fiber found in a colored resin Epoxy-vinyl ester matrix, which becomes translucent in the mouth. The post color will re-appear on command under cold air, water spray if re-treatment is required.

The D.T. LIGHT-POST ILLUSION® X-RO® is a smooth-surfaced post, apically tapered with a length of 20 mm. The optimized cylindrical-conical design limits the removal of dentin, thereby increasing the strength and clinical quality of the restoration.

The radiopacity level of the D.T. LIGHT-POST ILLUSION® X-RO® is 340% Al Equivalent (ISO 4049). The devices are available in 4 sizes and are used in association with D.T. LIGHT-POST® drills. An universal pre-shaping drill to perform a first hole and a finishing drill associated to each post size according to an O-ring color part (See Figure 1).

D.T. LIGHT-POST® drills are associated with:

- D.T. LIGHT-POST posts
- D.T. LIGHT-POST ILLUSION X-RO posts

	■ 0,5	■ 1	■ 2	■ 3
Ø apical tip	0,80	0,90	1,00	1,20
ØPost head	1,25	1,52	1,80	2,18
Double taper	.02/04	.02/06	.02/08	.02/10

COMPOSITION

D.T. LIGHT-POST ILLUSION® X-RO® posts are fiber-reinforced polymer, made of XRO fiber (80%) and epoxy resin matrix (20%) which contains colored pigment.

D.T. LIGHT-POST® drills are made of stainless steel

INTENDED USE

D.T. LIGHT-POST ILLUSION® X-RO® posts have been designed to support and secure the coronal restoration in case of insufficient residual tooth structure (<4 mm).

D.T. LIGHT-POST® drills are intended to prepare the tooth for insertion of a post that will support and secure the coronal restoration, where insufficient residual tooth structure (<4 mm) exists.

These products should only be used in hospitals, clinics or dental office with qualified dental personnel (dentists).

INDICATIONS FOR USE

D.T. LIGHT-POST ILLUSION® X-RO® posts and associated drills are indicated in case of insufficient residual tooth structure (<4 mm) to support and secure the coronal restoration.

CONTRAINDICATIONS

D.T. LIGHT-POST ILLUSION® X-RO® posts and associated drills are contraindicated if coronal tooth structure is less than 1.5 mm.

WARNINGS

None known to date.

PRECAUTIONS

There is no data for the use of D.T. LIGHT-POST® ILLUSION® X-RO® posts and associated drills for pregnant women, breast-feeding women and children under 18 years old (for posts only). It is not recommended for such use.

If the device is used in a way different than the one recommended in the IFU, loosening, breakage of the post or extraction of the tooth might occur.

D.T. LIGHT-POST® ILLUSION® X-RO® posts:

- The post must be immersed or soaked in isopropanol for 1 minute and dried gently with air syringe for 15 s at least prior to bonding before insertion in the canal. Avoid touching the posts with your fingers after that.

The post can neither be sterilized nor disinfected in thermosinfector and the use of autoclave is not recommended.

- Posts are single use devices. They must not be re-used to avoid the risk of contamination.
- Shortening the post with a diamond disc should be done outside of the mouth. Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post.

RTD warrants this product will be free from defects in material and manufacture. RTD makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and RTD's sole obligation shall be repair or replacement of the RTD product.

The post can neither be sterilized nor disinfected in thermosinfector and the use of autoclave is not recommended.

- Posts are single use devices. They must not be re-used to avoid the risk of contamination.
- Shortening the post with a diamond disc should be done outside of the mouth. Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post.

RTD warrants this product will be free from defects in material and manufacture. RTD makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and RTD's sole obligation shall be repair or replacement of the RTD product.

The post can neither be sterilized nor disinfected in thermosinfector and the use of autoclave is not recommended.

- Posts are single use devices. They must not be re-used to avoid the risk of contamination.
- Shortening the post with a diamond disc should be done outside of the mouth. Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post.

RTD warrants this product will be free from defects in material and manufacture. RTD makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and RTD's sole obligation shall be repair or replacement of the RTD product.

The post can neither be sterilized nor disinfected in thermosinfector and the use of autoclave is not recommended.

- Posts are single use devices. They must not be re-used to avoid the risk of contamination.
- Shortening the post with a diamond disc should be done outside of the mouth. Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post.

RTD warrants this product will be free from defects in material and manufacture. RTD makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and RTD's sole obligation shall be repair or replacement of the RTD product.

The post can neither be sterilized nor disinfected in thermosinfector and the use of autoclave is not recommended.

- Posts are single use devices. They must not be re-used to avoid the risk of contamination.
- Shortening the post with a diamond disc should be done outside of the mouth. Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post.

RTD warrants this product will be free from defects in material and manufacture. RTD makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and RTD's sole obligation shall be repair or replacement of the RTD product.

The post can neither be sterilized nor disinfected in thermosinfector and the use of autoclave is not recommended.

- Posts are single use devices. They must not be re-used to avoid the risk of contamination.
- Shortening the post with a diamond disc should be done outside of the mouth. Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post.

ADVERSE REACTIONS

- Debonding of the post
- Root fracture
- Fracture of the post

STEP-BY-STEP INSTRUCTIONS

- Root canal preparation: remove gutta-percha to preplanned depth with a Gates-Glidden, Peeso reamer or Largo. Radiographic verification is recommended.
- Select the correct D.T. LIGHT-POST® ILLUSION® X-RO® size. Determine the post size according to the anatomical situation of the tooth using the radiograph and the plastic calibrating card.
- Determine the length of the coronal build-up, which will account for 1/3 of the final post length. Remove the root filling with the D.T. Preshaping Drill (rotation speed 1000-2000 rpm) as deep as necessary for the post to be inserted to 2/3 of the root length. At least 4 mm of root canal filling should remain in the apical region. In case of curved canal, this calculated length is not suitable.
- Shape the canal with the D.T. Finishing Drill corresponding to the selected D.T. LIGHT-POST® ILLUSION® X-RO® post size, (rotation speed 1000-2000 rpm).
- Try in the post, to check that the post is properly seated in the canal
- Shorten the post to its final length with a diamond disc outside of the mouth. Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post.
- Clean the post with alcohol.

Apply etchant in the post space and to the exposed dentin for 15 seconds. Rinse for 10 seconds. Remove the excess water with paper points but leave surface moist. Do not put the etchant in contact with the gingiva.

- Apply 2 coats* of the primer with a thin brush (e.g. Compositush RTD) in the post space. Remove excess using paper points and gently air dry all surfaces. Etchamel and dentin surfaces should have a moist, glossy appearance. If not, repeat application. Light-cure the primer! In the post space for 10 to 20 seconds.
- Apply a single coat of primer* to the post. Gently air dry 5 seconds and light-cure the post for 10 to 20 seconds (if needed) outside the mouth.

11.a. If using a dual cure resin cement for both cementation and core build up such as Corecem™ (recommended technique), mix the cement and apply it in the post space from bottom to the top using the root canal tips. Seat the post immediately. Light-cure for 40 to 60 seconds. Continue applying the core build up using a plastic form, or directly on the tooth.

11.b. If using separate cement and a core build up resins: Mix the cement* and apply it to the post and then in the post space using a Lentulo spiral or a syringe tip. Seat the post immediately. Remove excess cement with appropriate instruments. If using dual-cure cement, light-cure for 40 to 60 seconds, applying gentle pressure on the post with the tip of the curing light probe. Apply 2 coats of bonding agent* to the exposed post(s), cement and involved tooth structure. Remove excess and air-dry with jet of air. Light-cure for 20 seconds (if needed).

11.c.i. If using a self-etching cement*, there is no need to etch and prime the post space. Mix the cement and apply it in the post space from bottom to the top using a root canal tip. Seat the post immediately. Remove excess cement. Light cure if needed.

11.c.ii. After etching, rinsing, drying of involved tooth structure, apply 2 coats of bonding agent* to the exposed post(s) and coronal part. Remove excess and air dry with jet of air. Light cure for 20 seconds (if needed).

12. Directly model the build-up using a composite core build-up material*. RTD recommends that the coronal end of the post be covered by core composite. In cases where a prosthetic crown is not indicated, a minimum of 1 mm of composite restorative should cover the post and that should be an appropriately-shaded, polishable, aesthetic composite. If the ILLUSION post cannot be adequately covered with composite, then a non-ILLUSION version can be utilized.

STORAGE

Do not store the post in direct sun light or artificial light.

DISPOSAL

Dispose in a safe manner in accordance with local/national regulations.

SAFETY AND LIABILITY

WARRANTY

RTD warrants this product will be free from defects in material and manufacture. RTD makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and RTD's sole obligation shall be repair or replacement of the RTD product.

LIMITATION OF LIABILITY

Except where prohibited by law, RTD will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

For dental use only.

* Refer to the manufacturer's instructions for use.

** Please refer to your autoclave manufacturer for specific sterilization instructions if needed.

TERMS AND SYMBOLS IDENTIFICATION

Manufacturer
Batch code
Catalogue number
Number of unit / Content
Keep away from sunlight

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