Plaque Test

EN Instructions for Use

- Disclosing liquid to make bacterial biofilm visible

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English

Description

Plaque Test is a disclosing liquid used to reveal plaque on the surfaces of natural and restored teeth as well as implants. Under a polymerization light, e.g. Bluephase®, plaque appears fluorescent yellow after Plaque Test has been applied.

Composition

Glycerin, Aqua, CI 45350, Potassium Phosphate, Ethylparaben, Sodium Hydroxide

Indication

Exposure of plaque on tooth surfaces

Contraindication

The application of Plaque Test is contraindicated if a person is known to be allergic to any of its ingredients.

Side effects

None known to date

Application

Use Plaque Test to selectively stain biofilm on teeth. This measure is particularly recommended to reveal retention niches for plaque during orthodontic treatment or if restorations or implants are present. The fluorescent pigment does not stain teeth, soft tissue and restorations. Only the plaque-coated areas glow under the light of a polymerization unit. This method can also be used in the determination of plaque indices.

Step by Step

1. Apply Plaque Test generously on all tooth surfaces using a cotton swab.
2. The mouth should be rinsed so that Plaque Test adheres only to areas covered with plaque.
3. Shine a polymerization light, e.g. Bluephase® on the teeth. Plaque covered with the disclosing liquid has a fluorescent yellow appearance.
4. Remove plaque by brushing the teeth or by professional tooth cleaning using Proxyl® pastes.

Note

The glycerine gives Plaque Test a sweet taste. The disclosing liquid is sugarless.

Storage

- Close bottle immediately after use.
- Shelf life: see label / packaging.
- Do not use Plaque Test after the indicated date of expiration.
- Storage temperature: 2–28 °C/36–82 °F.

Store out of the reach of children!

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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