Politip

[en] Instructions for Use

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1 Intended use

Intended purpose

Intraoral finishing and polishing of dental composite and amalgam restorations

Patient target group

- Patients with direct and indirect composite restorations
- Patients with amalgam restorations

Intended users

- Dentists (clinical procedure)

Special training

None

Use

For dental use only.

Description

Dental instruments for intraoral finishing (grey instrument) and polishing (green instrument) of all common dental composites and amalgam.

Indications

None

Areas of application:

- Intraoral finishing and polishing of dental composite restorations
- Intraoral finishing and polishing of dental amalgam restorations

Contraindications

The use of this product is contraindicated if the patient is known to be allergic to any of its ingredients.

Limitations of use

None

Reprocessing restrictions

Maximum 10 reprocessing cycles. Before use, visually inspect the product for damage.

Side effects

None known to date

Interactions

None known to date

Clinical benefit

- Restoration of esthetics

Composition

Polishing body: synthetic rubber, silicon carbide particles and titanium oxide

Handle: stainless steel

2 Application

Finishing / Polishing the restoration

- ✓ The product has been properly reprocessed. See chapter "Cleaning, disinfection and sterilization".
- ✓ The handpiece is technically and hygienically impeccable, properly maintained, and cleaned.
 1. Secure the product in the handpiece.
 2. Finish / polish the surface using gentle rotary motions without jamming or leverging, applying light pressure, ample water spray (> 50ml/min.) and efficient suction (approx. 2N; 7,500-10,000 RPM with accurate rotation).
 3. Remove the product from the handpiece.
 4. Reprocess the product. See chapter Cleaning, disinfection and sterilization.

3 Cleaning, disinfection and sterilization

NOTICE! Observe the country-specific regulations and guidelines regarding the hygiene standards and reprocessing of medical devices in dental practices. Wear protective goggles and gloves that comply with the requirements of the PPE regulation (EU) 2016/425 when handling any used and contaminated instruments.

NOTICE! Do not use strong acids or bases for cleaning.

Pre-cleaning the product

1. Pre-clean the product under running water (of drinking water quality) for at least 60 seconds using a brush with plastic bristles.
2. In case of coarse contamination: Spray the product with a water spray gun under high pressure for approx. 60 seconds.

Cleaning and disinfecting the product

Clean and disinfect contaminated products using automated or manual reprocessing procedures. Automated reprocessing in a washer-disinfector is preferred over manual reprocessing.

A) Clean and disinfect the product using machine cleaning

NOTICE! Make sure that the effectiveness of the washer-disinfector has been verified according to ISO 15883. The responsibility for the validation, performance qualification, and periodic routine testing of the reprocessing procedures lies with the operator.

NOTICE! Use a thermal program (temperature 90–95 °C / 194–203 °F) according to the validated directions stated in these Instructions. Observe the instructions of the manufacturer.

NOTICE! Only use cleaning / disinfection solutions that are classified by the manufacturer as "suitable for rubber polishers or synthetics/silicones". Observe the instructions of the manufacturer. The information on machine cleaning in these Instructions for Use refers to the cleaning agent neodisher MediClean Dental/Dr. Weigert.

1. Place the product in the small parts sieve/loading tray of the washer-disinfector.
2. Pre-rinse the product for 5 min with cold water.
3. Clean the product for 10 min at 55 °C / 131 °F using the cleaning agent.
4. Perform the final rinse of the product for 2 min using deionized, low-germ (max 10 cfu/ml) and low-endotoxin (max. 0.25 EU/ml) water.
5. Disinfect the product for 5 min at 90 °C / 194 °F using distilled water (A0 value > 3000).
6. Dry the product for 15 min at >120 °C / >248 °F.
7. Dry the product with a clean, lint-free cellulose cloth.
8. Inspect the product for integrity and cleanliness.
9. If any macroscopically visible residual contamination is present: Repeat the cleaning and disinfection procedure.

B) Clean and disinfect the product manually

NOTICE! Only use cleaning / disinfection solutions specified by the manufacturer as "suitable for rubber polishers or synthetics/silicones".

Observe the instructions of the manufacturer. The information on manual cleaning in these Instructions for Use refers to the cleaning agent ID 215/Dürr Dental (concentration 2 %) and the disinfection agent ID 212/Dürr Dental (concentration 2 %).

NOTICE! Observe the following parameters:

Step	Sonication time in the ultrasonic bath	Reaction time in the ultrasonic bath	Rinsing
Cleaning	1 min	–	30 s, drinking water
Disinfection	2 min	5 min	30 s, deionized, low-germ (max. 10 cfu/ml) and low-endotoxin (max. 0.25 EU/ml) water

Table 1

1. Prepare the cleaning / disinfection solution.
2. Pour the solution into the ultrasonic bath.
3. Immerse the product with the cleaning basket in the solution.
4. Sonicate the product in the ultrasonic bath, see Table 1.
5. Allow the product to remain in the ultrasonic bath, see Table 1.
6. Remove the product with the cleaning basket from the ultrasonic bath.
7. Rinse the product, see Table 1.
8. Dry the product with a clean, lint-free cellulose cloth.
9. Inspect the product for integrity and cleanliness.
10. If any macroscopically visible residual contamination is present: Repeat the cleaning and disinfection procedure.

Packaging, sterilizing and storing the product

NOTICE! Observe the following parameters:

	Method	Conditions	Drying time
1	Steam sterilization (autoclave) Fractionated vacuum	134 °C / 273.2 °F for 5 min	Local practice
2	Steam sterilization (autoclave)* Fractionated vacuum	132 °C / 269.6 °F for 3 min	10 min
3	Steam sterilization (autoclave)** Fractionated vacuum	134 °C / 273.2 °F for 3 min	Local practice
4	Steam sterilization (autoclave)*** Fractionated vacuum	134 °C / 273.2 °F for 18 min	Local practice

Table 2

* Recommended for the USA; ** Recommended for the UK; *** Recommended for Switzerland and France

For countries outside the United States:

NOTICE! Use a sterile barrier system that is compliant with ISO 11607-1 and designated by the manufacturer for steam sterilization and a sealing seam device. Observe the instructions of the manufacturer. The information on packaging in these Instructions for Use refers to the foil paper packaging steriCLIN and the sealing seam system HAWO, type 880 DC-V.

NOTICE! Make sure that the packaging is sufficiently large so that the seam is not stretched when the product is inserted.

For the United States:

NOTICE! Make sure that the sterilizer and any sterilization accessories and the sterile barrier system being used are approved by the FDA for the intended sterilization.

1. Insert the product in the sterile barrier system (ISO 11607-1).
2. Seal the sealing seam using the heat sealing procedure.
3. Visually inspect the sealing seam for defects.
4. In the event of sealing defects: Repackage the product.
5. Repeat the steps 1 to 4.
6. Place the double-bagged product in the sterilizer.
7. Steam-sterilize the product, see Table 2.
8. Store the product for a maximum of 48 h.

4 Safety information

- The product has been developed solely for use in dentistry. Process according to the Instructions for Use.
- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, www.ivoclar.com and your responsible competent authority.
- Wear a face mask, protective gloves and goggles.
- Check to make sure that the packaging and the product are intact and undamaged before use. If in doubt, contact Ivoclar Vivadent AG or your local sales partner.

Warnings

CAUTION

Improper use

Damage of the restoration and the surrounding soft tissue
Damage of the instruments

- a) Use only instruments that are technically impeccable and maintained.
- b) Securely lock the product in the handpiece.
- c) Comply with the maximum rotary speed.
- d) Comply with the maximum contact pressure.
- e) Avoid local overheating.
- f) User water spray.

CAUTION

Improper reprocessing

Contamination, damage of the product

- a) Follow the reprocessing instructions.

Supporting documents

Document	To be found at:
Current version of Instructions for Use	www.ivoclar.com/eIFU
Structure of Instructions for use and Warnings	www.ivoclar.com/eIFU
Explanation of symbols	www.ivoclar.com/eIFU

Disposal information

Remaining stocks must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following known residual clinical risks exist:

- Ingestion of material
- Aspiration of material

5 Additional information

Keep material out of the reach of children!

The manufacturer accepts no liability for damages resulting from a failure to observe the Instructions for Use. Furthermore, the user is responsible for testing the products for their suitability for any purpose not explicitly stated in the Instructions.

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