**Intended use****Intended purpose**

Absolute isolation

**Use**

For dental use only.

**Description**

OptraDam® Plus is an anatomically shaped rubber dam for the absolute isolation of the working field. It provides retraction of lips and cheeks during dental treatment and offers the clinician increased visibility and accessibility. In addition, the device facilitates the evacuation of water spray. OptraDam Plus features high flexibility and elasticity in all directions of movement, provides patient comfort and assists in keeping the patient's mouth open. OptraDam Plus is available in regular and small sizes.

It is suitable for use in therapeutic dental procedures, for which an absolutely dry and isolated working field is required (e.g. for the adhesive cementation of restorations, direct restorative procedures and root canal treatments).


**Indication**

Isolation of single and multiple teeth during dental treatment

**Contraindications**

The use of OptraDam Plus is contraindicated in patients with latex allergy.

**Limitations of use**

 OptraDam Plus is for single patient use only. The product is not intended to be reprocessed or reused. The material is not suitable for sterilization. Cleaning and disinfecting agents may have a detrimental effect on the material and lead to undesirable reactions in the patient.

**Side effects**

 OptraDam Plus contains natural rubber latex which may cause allergic reactions.

**Composition**

OptraDam Plus is made of latex and POM and it is powdered with corn starch.

**Application**

When OptraDam Plus is used, the small intraoral ring of the device is positioned in the area of the gingivobuccal fold, while the outer ring remains outside the mouth.

The elastic component between the two rings thus embraces the lips of the patient and provides retraction due to the restoring force of the rings (see Figure 1).

# OptraDam® Plus

<b>EN Instructions for Use</b>	<b>NO Bruksanvisning</b>	<b>SK Návod na používanie</b>	<b>LV Lietošanas instrukcija</b>
- Anatomically shaped rubber dam	- Anatomičk formet kofferdam	- Anatomický tvárnany koferdám	- Anatomišķas formas gumijas koferdams
<b>DE Gebrauchsinformation</b>	<b>NL Gebruiksaanwijzing</b>	<b>HU Használati utasítás</b>	<b>LT Naudojimo informacija</b>
- Anatomisch ausgeförmter Kofferdam	- Anatomisch gewordene cofferdam	- Anatómiaiakkal formált kofferdám	- Anatomiškai suformuota guminė užtvara
<b>FR Mode d'emploi</b>	<b>EL Οδηγίες Χρήσεως</b>	<b>SR Упутство за употребу</b>	<b>MK Упатство за употреба</b>
- Digue anatomique en latex	- Ανατομικά διαμέρισμένος ελαστικός απορρυτήρας	- Анатомски обликован гумени кофердам	- Анатомична обликуван кофердам
<b>IT Istruzioni d'uso</b>	<b>TR Kullanma Talimatı</b>	<b>BG Инструкции за употреба</b>	<b>SQ Udhëzime përdorimi</b>
- Digue anatomica	- Anatomičk olarak sekilendirilmiş rubber dam	- Анатомично обликуван кофердам	- Дигеј гоме me формë anatomike
<b>ES Instrucciones de uso</b>	<b>RU Инструкция по применению</b>	<b>PL Instrukcja stosowania</b>	<b>RO Instructiuni de utilizare</b>
- Dique de goma anatómico	- Коффердам анатомической формы	- Anatomicznie ukształtowany koferdam	- Diga din cauciuc cu forma anatomica
<b>PT Instruções de Uso</b>	<b>PT Instruções de Uso</b>	<b>SL Navodila za uporabo</b>	
- Dique de borracha com conformação anatómica	- Dique de borracha com conformação anatómica	- Anatomisk oblikovana gumijsasta opna	
<b>SV Bruksanvisning</b>			
- Anatomiškt utformad koferdam			

Fig. 1: OptraDam Plus in place

1. Select the appropriate sized OptraDam Plus. The size (R for regular or S for small) is indicated on the packaging and on the back of the latex bag.
2. The rubber dam punch is used to make the appropriate holes, depending on the teeth to be exposed. Check the quality of your rubber dam punch regularly. Blunt punches perforate the latex only incompletely, which results in the latex material tearing and thus in insufficient isolation of the teeth.
3. OptraDam Plus may be inserted from both the 12 o'clock position with the patient lying down, or from the 8 o'clock position with the patient sitting upright (or, the 4 o'clock position if the clinician is left-handed).
4. To ease placement, the operator should grasp the intraoral ring by reaching through the outer ring. The intraoral ring should be held between the thumb and middle finger and pressed together slightly (see Figure 2). When grasping the intraoral ring, make sure that it is not completely pulled out of the latex bag.


Fig. 2: Correct hold on the intraoral ring with thumb and middle finger.  
The ring is slightly compressed in the process.

5. The slightly compressed intraoral ring is then inserted into the buccal corridor on one side, so that the elastic component embraces the corner of the mouth and the outer ring is still located extraorally (see Figure 2).
6. As soon as OptraDam Plus is securely in place on one side, the intraoral ring is inserted in the same manner on the other side by slightly bending it (see Figure 2).
7. Subsequently, the intraoral ring is placed behind the lower and upper lips and thus the OptraDam Plus assumes its final, stable position (see Figure 3). The intraoral ring can be inserted with more ease if the patient's mouth is relaxed. If, in individual cases, the intraoral ring should show a tendency to slip out of the gingivobuccal fold if the mouth is closed completely, positioning of the intraoral ring deeper inside the vestibule is usually sufficient to solve the problem. Sometimes, the use of another size may also be a solution.
8. The dam is secured around the teeth by fitting the septum of the dam interproximally and in the sulcus using dental floss.

Fig. 3: After having positioned the intraoral ring behind the corners of the mouth, it is pushed behind the lower and upper lip.

9. The securely positioned OptraDam Plus improves the overall view and accessibility of the working field. The saliva ejector is introduced below the integrated frame and positioned in the corner of the mouth. Lateral movements of the mandible and opening of the mouth to various degrees may provide additional space and access to the treatment area. If the interproximal spaces are large and retention of the rubber dam is insufficient, a ligature of dental floss may be used additionally to hold the latex bag in place.
10. Removal of the rubber dam: Carefully cut the interdental latex septa. Using a paper towel, grasp the outer ring in the lower area and lift it slightly so that the intraoral ring is loosened from the lower gingivobuccal fold. The intraoral ring can then be easily removed from the upper gingivobuccal fold and hygienically disposed of with the paper towel.

#### Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: [www.ivoclarvivadent.com](http://www.ivoclarvivadent.com), and your responsible competent authority.
- The current Instructions for Use are available in the download section of the Ivoclar Vivadent AG website ([www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)).

#### Disposal information

- Remaining stocks must be disposed of according to the corresponding national legal requirements.

#### Shelf life and storage

- Storage temperature 2–28 °C (36–82 °F)

- Keep out of direct sunlight.

- Do not use the product after the indicated date of expiration.

- Expiration date: See note on packaging

#### Additional information

Keep out of the reach of children!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

**Ivoclar Vivadent Pty. Ltd.**

1 – 5 Overseas Drive, P.O. Box 367, Noble Park, Vic. 3174, Australia

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[www.ivoclarvivadent.us](http://www.ivoclarvivadent.us)

**Ivoclar Vivadent Shanghai Trading Co., Ltd.**

2/F Building 1, 881 Wuding Road, Jing An District, 200040 Shanghai, China

Tel. +86 21 6032 1657, Fax +86 21 6176 0968

[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

**Ivoclar Vivadent Marketing Ltd.**

Calle 134 No. 7-B-83, Of. 520, Bogotá, Colombia

Tel. +57 1 627 3399, Fax +57 1 633 1663

[www.ivoclarvivadent.co](http://www.ivoclarvivadent.co)

**Ivoclar Vivadent SAS**

B.P. 118, 74410 Saint-Jorioz, France

Tel. +33 4 50 88 64 00, Fax +33 4 50 68 91 52

[www.ivoclarvivadent.fr](http://www.ivoclarvivadent.fr)

**Ivoclar Vivadent GmbH**

Dr. Adolf-Schneider-Str. 2, 73479 Ellwangen, Jagst, Germany

Tel. +49 7961 889 0, Fax +49 7961 6326

[www.ivoclarvivadent.de](http://www.ivoclarvivadent.de)

**Ivoclar Vivadent Marketing (India) Pvt. Ltd.**

503/504 Raheja Plaza, 15 B Shah Industrial Estate

Veera Desai Road, Andheri (West), Mumbai, 400 053, India

Tel. +91 22 2673 0302, Fax +91 22 2673 0301

[www.ivoclarvivadent.in](http://www.ivoclarvivadent.in)

**Ivoclar Vivadent Marketing Ltd.**

The Icon, Horizon Broadway BSD, Block M5 No. 1

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[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

**Ivoclar Vivadent Ltd.**

12 Omega St, Rosedale, PO Box 303011 North Harbour

Auckland 0751, New Zealand

Tel. +64 9 914 9999, Fax +64 9 914 9990

[www.ivoclarvivadent.co.nz](http://www.ivoclarvivadent.co.nz)

**Ivoclar Vivadent Polska Sp. z o.o.**

Al. Jana Pawła II 78, 00-175 Warszawa, Poland

Tel. +48 22 635 5496, Fax +48 22 635 5469

[www.ivoclarvivadent.pl](http://www.ivoclarvivadent.pl)

**Ivoclar Vivadent LLC**

Prospekt Andropova 18 korp. 6/, office 10-06, 115432 Moscow, Russia

Tel. +7 499 418 0300, Fax +7 499 418 0310

[www.ivoclarvivadent.ru](http://www.ivoclarvivadent.ru)

**Ivoclar Vivadent Marketing Ltd.**

Qlaya Main St., Siricon Building No.14, 2<sup>nd</sup> Floor, Office No. 204

P.O. Box 300146, Riyadh 11372, Saudi Arabia

Tel. +966 11 293 8345, Fax +966 11 293 8344

[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

**Ivoclar Vivadent S.L.U.**

Carretera de Fuencarral nº24, Portal 1 – Planta Baja  
28108-Alcobendas (Madrid), Spain

Tel. +34 91 375 78 20, Fax +34 91 375 78 38

[www.ivoclarvivadent.es](http://www.ivoclarvivadent.es)

**Ivoclar Vivadent AB**

Dalvägen 14, 169 56 Solna, Sweden

Tel. +46 8 514 939 30, Fax +46 8 514 939 40

[www.ivoclarvivadent.se](http://www.ivoclarvivadent.se)

**Ivoclar Vivadent Liaison Office**

: Tesvikiye Mahallesi, Sakayik Sokak, Nisantas' Plaza No:38/2

Kat:5 Daire:24, 34021 Sisli – Istanbul, Turkey

Tel. +90 212 343 0802, Fax +90 212 343 0842

[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

**Ivoclar Vivadent Limited**

Compass Building, Feldspar Close, Warrens Business Park

Enderby, Leicester LE19 4SD, United Kingdom

Tel. +44 116 284 7880, Fax +44 116 284 7881

[www.ivoclarvivadent.co.uk](http://www.ivoclarvivadent.co.uk)

**Ivoclar Vivadent, Inc.**

175 Pineview Drive, Amherst, N.Y. 14228, USA

Tel. +1 800 533 6825, Fax +1 716 691 2285

[www.ivoclarvivadent.us](http://www.ivoclarvivadent.us)