Kerr™

Instructions For Use

Demi^{Plus}™

LED Light Curing System
Powered by Demetron™ Technology

Table of Contents

en - ENGLISH.....	3
fi - SUOMI	9
es - ESPAÑOL	14
fr - FRANÇAIS.....	19
tr - TÜRKÇE	24
de - DEUTSCH	29
hu - MAGYAR	34
ro - ROMÂNĂ	39
nl - NEDERLANDS	44
no - NORSK.....	49
pl - POLSKI	54
pt - PORTUGUÊS	60
it - ITALIANO	67
cs - ČESKY	70
sv - SVENSKA	75
et - EESTI.....	80
da - DANSK.....	85
bg - БЪЛГАРСКИ	90
el - ΕΛΛΗΝΙΚΑ	96
sk - SLOVENSKY.....	102
sl - SLOVENSKO	107
ru - РУССКИЙ.....	112
lt - LIETUVIŲ.....	118
sr - SRPSKI	123
is - ÍSLENSKA.....	128
uk - УКРАЇНСЬКА	133
ko - 한국어	138
ja - 日本語	143
zh - 简体中文	148
zh-t - 繁體中文	153
id - INDONESIA	158
vi - TIẾNG VIỆT	163
he - עברית.....	169
ar - العربية.....	174

Illustrations

Figure 1



Figure 2

Figure 3

Figure 4

Demi Plus LED Light Curing System

OPERATOR MANUAL

GENERAL INFORMATION

The Demi Plus is a wand-style curing light used for the polymerization of light-cured dental material and is powered by a rechargeable lithium-ion battery pack. The unit employs proprietary Periodic Level Shifting (PLS) technology.

INDICATIONS FOR USE

The Demi Plus is an LED visible light curing unit intended for polymerization of light cure materials by dental professionals.

CONTRAINDICATIONS

None Known.

WARNINGS

- I. Demi Plus emits high intensity light and must be only used as indicated in this manual.
- II. Do not use the Demi Plus curing light without a Light Guide.
- III. Do not look directly at the light emitted from this curing light. Suitable protective eyewear is required for medical staff and a patient when using Demi Plus.
- IV. Persons having a history of photosensitive reactions or who are using photosensitizing drugs should not be exposed to light emitted from the Demi Plus.
- V. The Demi Plus requires special precautions regarding Electromagnetic Compatibility (EMC) and must be installed and put into service according to the Guidance and Manufacturer's Declaration for EMC included in this manual.
- VI. Portable and mobile Radio Frequency (RF) communications equipment can affect the Demi Plus. Reference the Guidance and Manufacturer's Declaration for EMC included in this manual.
- VII. Use only the Kerr-authorized AC adapter and power cords.
- VIII. Use of other accessories that are not authorized for use with this device may cause malfunction and compromise patient safety.
- IX. No modification of this equipment is allowed.
- X. For the batteries air shipping and transportation, consult and follow the International Air Transport Association guidelines.

PRECAUTIONS

- I. This product is intended to be used by qualified dental professionals on a general patient population.
- II. The Demi Plus device components have not been cleaned and disinfected prior to shipping. The fiber-optic Light Guide is provided in a non-sterile condition. Perform the necessary cleaning, disinfection and sterilization steps prior to patient contact according to "Cleaning, Disinfection and Sterilization" Section of this manual.
- III. Do not place the Light Guide tip of the curing light directly on or towards unprotected gingiva, oral mucous membrane or skin. Prolonged exposure to the light emitted from the Demi Plus might result in a soft tissue injury.
- IV. Observe the manufacturer's specified curing times, particularly when using Demi Plus on areas near the pulp tissue. Prolonged exposure to the light emitted from the Demi Plus might result in irreversible pulpitis.
- V. Regularly inspect your curing light to ensure it is providing the expected light irradiance. It is recommended that curing capability be tested using a Hardness Disk to ensure a complete cure of the material.
- VI. The Demi Plus Charger should not be stored within 1.5 meters from the patient to prevent liquids from entering Charger vents.
- VII. Demi Plus disposable barrier is designed for single use only, to prevent cross contamination between patients.

PACKAGED CONTENTS

- 1 Demi Plus curing light body
- 1 Protective light shield
- 1 Plug-in power supply
- 1 Battery charger
- 1 Battery pack
- 1 Operator manual
- 1 Demi/Demi Plus disposable barrier

PRODUCT SPECIFICATIONS

Type of Information	Specifications
Light Emitting Diodes	An array of four LEDs die-mounted on a substrate.
Output peak wavelength range	450-470 nanometers
Baseline output	1,100 mW/cm ²
Curing Device Dimensions	Length: 9.25 in (23.5 cm)
	Width: 1.2 in (2.9 cm)
Charging Dock Dimensions	Length: 6.0 in (16.0 cm)
	Width: 4.63 in (11.1 cm)
Weight	Curing Device with Battery: 5.5 oz (155 g)
	Charging Dock: 12.4 oz (352 g)
Battery Performance	Type: Rechargeable Lithium-Ion Battery Output: 3.6 Volts nominal @ 2.6A-H capacity Performance: Curing device will generate 500, 5-sec cures with a fully charged battery. Battery pack life: 300 full charge/recharge cycles.
Battery Overcurrent Protection	Electronic and resettable polyfuse
Standard Light Guide(s)	P/N 921746, 13 mm-8 mm extended light guide
AC supply connection	100-240V AC, 0.8-0.4A, 47-63Hz
Power input	12V DC, 2.5A, nominal
Equipment class	Class II
Safety	IEC 60601-1
EMC (Electro-Magnetic Compliance)	IEC 60601-1-2
Protection from electric shock	Type BF
Protection from ingress of liquids	Curing Device with battery: IPX0
	Charging Dock: IPX0
Operation	Duty Cycle 20 sec ON / 1 min OFF
Operating Environment	Ambient temperature: 60°F to 95°F (16°C to 35°C)
	Relative humidity: 10% - 80%, non-condensing
	Atmospheric pressure: 0.5-atm to 1.0-atm (500 hPa - 1060 hPa)
Transport and Storage Environment	Ambient temperature: -4°F to 95°F (-20°C to 35°C)
	Relative humidity: 10% - 85%
	Atmospheric pressure: 0.5-atm to 1.0-atm (500 hPa - 1060 hPa)

ADVERSE EVENTS

If a serious incident occurs with this medical device, report it to the manufacturer and to the competent authority for that country in which the user and/or patient is established.

STEP-BY-STEP INSTRUCTIONS

ATTACH THE ADAPTER TO THE POWER SUPPLY

The power supply is a universal-input unit that accepts 100V - 240V AC nominal.

Note: Only Kerr-supplied power supply, shown in (Figure 3, ① part number 952269) shall be used to power the Demi Plus battery charger.

- I. Select the proper adapter (as shown in Figure 4, ① ② ③ ④). Their applications are as follows:
 - ① 120V US — Type A
 - ② Euro — Type C
 - ③ UK — Type G
 - ④ Australian — Type I
- II. Each adapter also has a slot into which the retaining latch on the power supply fits.
- III. With the slot oriented toward the power supply, place the adapter onto the power supply so that the tabs on the adapter are aligned with their respective slots on the power supply.
- IV. Press down on both ends of the adapter (Figure 3, ③), then slide the adapter into the power supply until a click is heard and the adapter is locked in place (Figure 3, ③). If the adapter feels loose, remove the adapter and re-install.
- V. (Figure 3, ④) shows the adapter locked in place on the power supply.
- VI. To remove the adapter from the power supply, press on the latch retainer button. (Figure 3, ②)
- VII. Slide the adapter away from the power supply and remove the adapter.
- VIII. Retain unused adapters for future use.

Demi Plus PARTS

The following parts are identified by Figure 1:

- Lithium-Ion Battery Pack ⑤ is used to power the curing light.
 - Light Guide ⑥. The output of the Curing Light LED is emitted from the end of this piece.*
 - Light Shield ⑦. This shield protects the eyes of the operator from the light emitted from the Light Guide when in use.
- *Other light guides that may yield different outputs are available for the Demi Plus.

INSTALLATION OF THE BATTERY PACK, LIGHT SHIELD AND LIGHT GUIDE ONTO THE CURING LIGHT BODY

- I. The battery pack will only fit onto the curing light body in one orientation. If you have difficulty fitting the battery pack to the curing light body, do not force it.
- II. Orient the battery pack (Figure 1) so that the broad flat area ⑨ is toward the side of the curing light body that contains the switches.
- III. Slide the battery pack and curing light body together until a click is heard and the battery pack is locked in place.
- IV. To remove the battery pack (Figure 1), press the textured indents ⑧ on both sides of the battery pack and slide the battery free.
- V. Install the light shield (Figure 1, ⑦) to the curing light body until the light shield is locked in place.
- VI. Install the light guide (Figure 1, ⑥) to the curing light body until a click is heard and the light guide is locked in place.

CHARGE THE CURING LIGHT / BATTERY PACK

Note: The battery pack contains lithium-ion cells. A new battery pack should be fully charged prior to its first use. If the battery pack is below room temperature, allow the battery pack to come up to room temperature before charging. Plug the power supply into the charger. The socket is located at the back of the charger (Figure 2, ②).

- I. Plug the power supply into an electrical outlet.
- II. Insert the curing light into the opening(s) of the charger.
 - a. The charger can charge any combination of two battery packs or curing lights (in differing states of charge) at the same time.
 - b. The lighted rings (Figure 2, ①) that surround the openings in the charger indicate the state of charge of each battery pack as follows:
 - Green light — battery pack is charged.

- Yellow light — battery pack is charging.
- Blinking Yellow Light (with battery pack or curing light in the charger) — indicates a fault condition. Please call Customer Service.

DEMI PLUS CONTROLS

The following controls are used by the operator during the curing process (Figure 1):

- I. Curing Mode LEDs ①. Three LEDs represent the three curing modes; from left to right these are: 5 second, 10 second, and 20 second cycles.
- II. Mode Selection Switch ②. Pressing this switch sequentially toggles the unit through the three curing modes. If the unit is in Sleep mode, pressing this switch awakens the unit in the curing mode last used.
- III. Trigger Switch ③. Pressing (and releasing) this switch initiates the chosen curing cycle. If the unit is in Sleep mode, pressing this switch awakens the unit in the curing mode last used.
- IV. Battery Power LED ④. This LED indicates the state of charge of the battery in the following way:
 - a. If the light is unlit, the battery pack is either CHARGED or the unit is in Sleep mode.
 - b. If the light is yellow, the battery pack has a LOW charge.
Note: The curing light will operate normally in this mode. Recharge the battery pack at the next available opportunity.
 - c. If the light is red, the battery pack is DISCHARGED (dead-battery) and will need to be recharged. So that no curing cycles are interrupted, when the battery voltage drops below the dead-battery threshold during a curing cycle, the curing cycle will complete and the battery status LED will not change to red until the curing cycle is completed.

DEMI PLUS CURING LIGHT OPERATING STATES

There are three distinct operational states

- I. **Idle:** The unit is Idle when it is not Curing or in Sleep mode and the battery is CHARGED or LOW. The user may cycle between 5 second, 10 second, and 20 second curing modes by pressing the mode switch. Curing cycles are also initiated when the unit is in the Idle state.
- II. **Curing:** A Curing cycle is initiated by pressing, and releasing, the trigger switch while the unit is Idle. A 5 second, 10 second, or 20 second cycle will be started, depending on the status indicated by the Curing Mode LEDs. Once a Curing cycle is started, the Curing Light LED will be turned on and a tone sounds. Depending upon which mode was selected, tone will sound every five seconds until the cycle is complete, at which time a final, longer duration, tone sounds.
- III. **Note:** Pressing the Trigger Switch for a second time during a Curing cycle immediately interrupts the Curing cycle and causes the final tone to sound, regardless of how much time has elapsed since beginning the cycle (pressing the Mode Selection Switch during a Curing cycle has no effect).
- IV. **Sleep:** The Demi Plus goes to Sleep mode to conserve battery life after approximately 8 minutes of no activity. It can be awakened by pressing either the mode or trigger switch once, at which time the unit will be in the Idle state. The unit will reflect the last cure mode chosen. The mode or trigger switches can then be pressed to switch between modes or to initiate a cure. While in Sleep mode, all LEDs are turned off and the curing light goes into a low power state.

OPERATE THE DEMI PLUS CURING LIGHT

- I. Place Demi Plus disposable barrier on the Demi Plus curing light.
- II. Provide suitable protective eyewear for the patient and medical personnel performing the procedure.
- III. Select the desired curing cycle by pressing the Mode Selection Switch. User may choose between 5, 10 and 20 seconds.
Note: Always consult the Instructions for Use of the applied restorative material when selecting the curing cycle.
- IV. Position the tip of the Light Guide directly over the material to be polymerized. Avoid direct contact of the Light Guide tip with the uncured material.
- V. Press and release the Trigger Switch to activate the curing cycle. A tone signal sounds every 5 seconds during the curing cycle. Once the selected curing cycle has elapsed, the curing light will switch off automatically. If needed, the curing light might be switched off during the curing cycle by pressing the Trigger Switch again.

- VI. When the clinical procedure with Demi Plus is finished, perform device reprocessing according to the instructions in the "Cleaning, Disinfection and Sterilization" Section of this manual. Remove the Demi Plus disposable barrier from Demi Plus prior to reprocessing.

CLEANING, DISINFECTION AND STERILIZATION

Preparation for the device reprocessing:

- I. Reprocess all Demi Plus components immediately after each use to prevent drying of soil and contaminants.
- II. Carefully discard examination gloves, rinse and sanitize hands using an appropriate hand sanitizer solution and use a new pair of examination gloves. Follow standard precautions for personal protection using cleaning agents and disinfectant solution as recommended by the manufacturer.
- III. Remove the Light Guide and Light Shield from the Curing Light Body prior to cleaning.
- IV. Do not disassemble the Battery Pack from the Curing Light Body prior to cleaning.
- V. Inspect the Light Guide and Light Shield prior to cleaning for the appearance of defects such as cracks, chippings, deformations which are indicators that the instrument is not in condition to be re-processed with the required level of confidence.
- VI. The use of automated cleaning devices or automated disinfecting devices is NOT recommended for the cleaning of any Demi Plus components.

Precautions for electronic components:

- I. Demi Plus Curing Light Body and Charger are not autoclavable as it will damage the circuitry.
- II. Do not spray the Demi Plus components with any liquids as it may damage the circuitry. Wipe the plastic surfaces with a cloth lightly moistened with cleaning/disinfecting solution.
- III. Do not allow liquids to collect in the Charger or in the light guide opening on the Curing Light Body as it may come in contact with the connectors and damage the circuitry.
- IV. Do not immerse the Curing Light Body and the Charger in any solution as it will void the warranty.

CLEANING AND DISINFECTION OF THE CURING LIGHT BODY, LIGHT SHIELD AND CHARGER

The Curing Light Body, Light Shield and Charger must be reprocessed after each use following the cleaning and disinfection steps below to avoid the risk of cross-contamination between patients.

All plastic parts may be safely cleaned and disinfected with an intermediate-level disinfectant which formulation contains a quaternary ammonium compound and 20% or less concentration of alcohol. Do not use denatured alcohol, Lysol®, phenol, ammonia complex or iodine complex solutions.

Cleaning:

- I. Clean all the surfaces with a cloth lightly moistened with a medical-grade intermediate level-disinfectant cleaner. Please refer to manufacturer's instructions regarding concentrations, temperature and recommended contact time.
- II. Use a cotton swab and a small, soft brush moistened with a fresh disinfectant cleaner to remove any soil that may have accumulated in crevices and hard to reach areas.
- III. Wipe the device components with additional cloths lightly moistened with disinfectant cleaner until no visible soil is detected on the used cloth.
- IV. Visually inspect all components to ensure cleanliness. If any contamination is visible, repeat the cleaning steps. Use one more cloth moistened with disinfectant cleaner to ensure that no residual contaminants are left on the device components.
- V. Finally, to remove detergent residue, use a clean, lint-free cloth lightly moistened with distilled water and thoroughly wipe all component surfaces. Allow to dry.

Disinfection:

- I. Use a new cloth lightly moistened with a medical-grade intermediate-level disinfectant and wipe all the components surfaces. Please refer to manufacturer's instructions regarding concentrations, temperature and recommended contact time.
- II. Use a cotton swab and a small, soft brush moistened with disinfectant to apply the substance to crevices and hard to access areas. Ensure that all the surfaces remain visibly wet for the period of time specified by disinfectant manufacturer.
- III. To remove disinfectant residue, use a clean lint-free cloth lightly moistened with distilled water and thoroughly wipe all component surfaces. Device components are ready for reuse when all surfaces are visibly dry.

CLEANING AND STERILIZATION OF THE FIBER-OPTIC LIGHT GUIDE

The Light Guide must be reprocessed after each use following the cleaning and sterilization steps below to avoid the risk of cross-contamination between patients. Prior to cleaning, disassemble the Light Guide from the curing light body. Inspect the Light Guide for any defects. Discard it if the following conditions are present: broken fibers (internal black spots or darkened areas), severely chipped or etched fiber ends.

Cleaning:

- I. Inspect the Light Guide for any residues of cured resin material at the tip of the Light Guide. If any resin material is found, remove it completely.
- II. Clean the Light Guide with a medical-grade non-enzymatic cleaning solution using a soft bristled brush until all visible soil is removed. Ensure to remove any soil that may have accumulated in crevices. Please refer to manufacturer's instructions regarding cleaning solution concentrations, temperature and recommended contact time.
- III. Rinse with distilled water for a minimum of 30 seconds. Visually inspect to ensure cleanliness. If any contamination is visible, repeat the cleaning steps until there is no visible contamination.
- IV. Dry with a clean, lint free cloth until there is no visible moisture.

Sterilization:

- I. Inspect the Light Guide for any defects. If damage is detected, do not proceed with sterilization process and discard the Light Guide.
- II. Place the Light Guide in a sterilization pouch suitable for steam sterilization. Ensure that used sterilization pouches comply with national requirements such as ISO 11607-1 and for US – use FDA-cleared accessories.
- III. Sterilize the Light Guide using one of the sterilization cycles listed below:

Type of Information Specifications		
Cycle	Gravity Displacement*	Pre-Vacuum
Temperature (°C)	121°C (250°F)	132°C (270°F)
Exposure Time (minutes)	30 minutes	4 minutes
Drying Time (minutes)	30 minutes	20 minutes

- IV. When the product is taken out of the autoclave, inspect the sterilization pouch for any damage. In case of a perforated pouch, do not use the product, as the sterilization could be compromised. Re-pack the product and repeat the sterilization procedure.
- V. After sterilization, place the pouches containing the Light Guides in a dry and dark place. Follow the instructions provided by the pouch manufacturer regarding storage conditions and maximum allowed time in storage.

NOTE: The instructions provided in the section "Cleaning, disinfection and sterilization Instructions" have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. Health Care facilities are responsible for making sure that the sterilization equipment is calibrated according to the manufacturer's manuals and specifications. In addition, health care facilities are responsible for training their staff on infection control, proper sterilization, cleaning and disinfection procedures.

STORAGE AND DISPOSAL

See "Transport and Storage Environment" under "Product Specification" Section.

BATTERY DISPOSAL: The battery pack contains a lithium-ion (Li-ion) battery. Recycle or dispose of batteries according to local and national regulations.

ELECTRONICS DISPOSAL: Recycle or dispose of electronics according to local and national regulations.

ACCESSORIES AND REPLACEMENT PARTS

Part No.	Fiber Optic Light Guides
20941	8 mm Curved Light Guide
21020	8 mm Curved Turbo Light Guide
921551	11 mm Curved Turbo Light Guide
921676	4 mm Curved Turbo Light Guide
921677	2 mm Light Guide
921746	8 mm Curved Turbo+ Extended
952213	8 mm Curved Turbo+ Light Guide
Part No.	Eye Protection
20816	Protective Light Shield
Part No.	Accessories
20399	Disposable Hardness Disk (pack of 6)
910726	Handheld LED Radiometer
Part No.	Replacement Parts
921918-1	Replacement Kit, Battery, Demi Plus
921919-1	Replacement Kit, Charger, Demi Plus
921920-1	Replacement Kit, Curing Light Body, Demi Plus
952269	Plug-in Universal Power Supply

WARRANTY

Kerr hereby warrants that for a period of three years from the date of purchase of Demi Plus (only one year for batteries), this instrument shall be free from defects in material and workmanship. For this warranty to apply, the customer must notify Kerr of the defect before the expiration of the warranty period. The product must be returned in order for this warranty to apply. If the product upon inspection by Kerr proves defective during the warranty period, Kerr, at its option, will either: (1) repair the product at no charge for parts or labor or (2) replace the product with a comparable product which may be new or refurbished. This warranty does not apply to any defect, failure or damage caused by customer's improper use or modification or inadequate maintenance of the product.

THE WARRANTY STATED HEREIN IS THE SOLE WARRANTY APPLICABLE TO KERR PRODUCTS. KERR EXPRESSLY DISCLAIMS ANY AND ALL OTHER WARRANTIES EXPRESSED OR IMPLIED, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. KERR'S LIABILITY WITH RESPECT TO ITS PRODUCTS IS EXPRESSLY LIMITED TO THE REMEDIES SET FORTH ABOVE. THE REMEDIES ARE THE CUSTOMER'S EXCLUSIVE REMEDIES. KERR SHALL UNDER NO CIRCUMSTANCES BE LIABLE FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES.

SYMBOLS DESCRIPTION

Full explanation of symbols used on Kerr packaging located at:

<http://www.kerrdental.com/symbols-glossary>

Symbol	Description
	Catalogue number
	Manufacturer
	European conformity
	Refer to instruction manual/booklet. Blue background and white person/booklet.
	Caution
	Authorized representative in the European Community
	Caution: Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner.
	Do not re-use
	Type BF applied part
	Class II equipment
	Radiated RF
	Waste Electronic and Electronic Equipment (WEEE)
	Serial number
	Temperature storage limitations
	Humidity limitation
	Atmospheric pressure limitation
	Keep dry
	1. Sterilizable in a steam sterilizer (autoclave) at temperature specified
	2. Medical Device
	China RoHS
	Importer

SYMBOLS DESCRIPTION

1. Sterilizable in a steam sterilizer (autoclave) at temperature specified
2. Medical Device

The Demi Plus is intended for use in the electromagnetic environment specified below. The customer or user of the Demi Plus shall insure it is used in such an environment.

GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC EMISSIONS

EMISSIONS TEST	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT -GUIDANCE
RF emissions CISPR 11	Group I	The Demi Plus uses RF energy only for its internal function. Therefore, its RF emissions are very low and not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Demi Plus is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power distribution network that supplies buildings or for domestic purposes.
Harmonic emissions	Class A	
Voltage fluctuations/flicker emissions IEC 60601-3-3	Complies	

GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC IMMUNITY

IMMUNITY TEST	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT -GUIDANCE
Electrostatic discharge (ESD) IEC 61000-4-2	±8kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV, 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±0.5 kV, ±1 kV line(s) to line(s) ±0.5 kV, ±1 kV, ±2 kV line(s) to ground	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and Voltage variations on power supply input lines IEC 61000-4-11	Voltage dips: 0% UT for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT for 1 cycle and 70% UT for 25/30 cycles at single phase (0°) Voltage interruptions: 0% UT, 250/300 cycle	Mains power quality should be that of a typical commercial or hospital environment. Note: UT is the AC mains voltage prior to the application of the test level.
Power frequency magnetic field IEC 61000-4-8	30 A/m; 50 Hz or 60 Hz	Power frequency magnetic field should be at levels characteristic of a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	3* Vrms 0.15 MHz – 80 MHz 80% AM at 1 kHz	Portable and mobile RF communications equipment should be used no closer to any part of the Demi Plus including the power supply and cable than the recommended separation distances calculated from the equation applicable to the frequency of the transmitter. Field strength from fixed RF transmitters, as determined by an electromagnetic site survey, ¹ should be less than the compliance level in each frequency range ² . Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz – 2.7 GHz, 80% AM at 1 kHz	

Note 1: At 80 Mhz and 800 Mhz the higher range applies

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

Note 3: *6 Vrms in the ISM bands. The ISM band is 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz

Note 4: UT is the a.c. mains voltage prior to application of the test level.

¹ Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radio, amateur radio, AM and FM broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. In the measured field strength in the location in which the Demi Plus is used exceeds the applicable RF compliance levels above, the Demi Plus should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Demi Plus.

² Over the frequency range of 150 KHz to 80 Mhz, field strengths should be less than 10V/m.

The Demi Plus is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the Demi Plus can help prevent electromagnetic distance between portable and mobile RF communications equipment (transmitters) and the Demi Plus as recommended below according to the maximum power output of the communications equipment.

RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT AND THE DEMI PLUS

Radiated Maximum output power of transmitter (W)	Separation distance according to frequency (m)		
	150 KHz to 80 Mhz, $d = 0.4\sqrt{P}$	80 MHz to 800 MHz, $d = 0.4\sqrt{P}$	800 MHz to 2.5 GHz, $d = 0.7\sqrt{P}$
0.01	0.1	0.1	0.1
0.1	0.2	0.2	0.3
1	0.4	0.4	0.7
10	1.3	1.5	2.2
100	4.0	4.0	7.0

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 Mhz and 800 Mhz the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

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**Kerr Corporation**
Kerr Australia Pty. Ltd.
**Kerr Italia S.r.l.**