# OptraGloss®

## [en] Instructions for Use

Universal polishing system for ceramic materials and composite resin

**Rx ONLY**

**C € 0123**

Date information prepared:  
2024-05-14 / Rev. 3

Manufacturer:  
Ivoclar Vivadent AG  
Bendererstrasse 2  
9494 Schaan/Liechtenstein  
[www.ivoclar.com](http://www.ivoclar.com)

**ivoclar**

English

### 1 Intended use

#### Intended purpose

- Extraoral and intraoral polishing of dental ceramic restorations
- Extraoral and intraoral polishing of dental composite resins

#### Patient target group

- Patients with indirect ceramic restorations
- Patients with direct composite restorations
- Patients with indirect composite restorations

#### Intended users / Special training

Dentists

#### Use

For dental use only.

#### Description

OptraGloss® is a universal polisher for extraoral and intraoral use and is suitable for the high-gloss polishing of restorations made of dental ceramics in two steps and for the high-gloss polishing of all popular dental composite resins in one step.

OptraGloss consists of the following components:

- Diamond pre-polisher (PP, dark blue) for the pre-polishing of ceramics in the shapes Flame, Cup and Lens
- Diamond high-gloss polisher (HP, light blue) for the high-gloss polishing of composite resins and ceramic materials in the shapes Flame, Cup, Lens and Spiral Wheel

OptraGloss is suitable for:

- the clinical polishing of glass-ceramic, feldspathic ceramic and oxide ceramic restorations
- the clinical polishing of composite resin restorations

#### Indications

None

#### Areas of application

- Clinical polishing of glass-ceramic, feldspathic ceramic and oxide ceramic restorations
- Clinical polishing of composite resin restorations

#### Contraindications

The use of the product is contraindicated if the patient is known to be allergic to any of its ingredients.

#### Limitations of use

None

#### Processing restrictions

The polishers are designed for multiple use and can be reprocessed up to 10 times if no apparent defect is found after the procedure.

#### Side effects

There are no known side effects to date.

#### Interactions

There are no known interactions to date.

#### Clinical benefit

Restoration of esthetics

#### Composition

- All shapes except the Spiral Wheel and Flame:  
Diamond material, dicyclohexylmethane-4,4'-diisocyanate, diethylmethylbenzenediamine, titanium dioxide  
The shanks are made of stainless steel.
- Spiral Wheel shape:  
Diamond material, dicyclohexylmethane-4,4'-diisocyanate, diisodecyl phthalate, titanium dioxide.  
The shanks are made of stainless steel.
- Flame shape:  
Diamond material  
The shanks are made of stainless steel.

### 2 Application

#### Ceramics

- 1<sup>st</sup> step: For the pre-polishing of ceramic restorations, use OptraGloss PP (dark blue) first.
- 2<sup>nd</sup> step: For high-gloss polishing, use OptraGloss HP (light blue).

#### Composite resin

- 1<sup>st</sup> step: For the high-gloss polishing of composite resin restorations, use OptraGloss HP (light blue).

#### Application notes

- Max. rotation speed: 10,000 rpm.
- Use with copious water spray (>50 ml/min).
- Apply medium contact pressure.
- Always bring the polisher to rotation before you apply it on the surface to be treated.
- Polishing should be performed using light circular movements to avoid the creation of grooves.
- Ensure adequate evacuation during polishing.
- Polishing can be performed intraorally (with water cooling) or extraorally (without water cooling).

### 3 Cleaning, disinfection and sterilization

General notes:

- Please observe the country-specific regulations and guidelines regarding the hygiene standards and reprocessing of medical devices in dental practices.
- Protective gloves and goggles complying with the requirements of the regulation on personal protective equipment PPE (EU) 2016/425 must be worn when handling any used and contaminated instruments.
- New polishers are not sterile and must be sterilized before they are used for the first time.
- The polishers must be cleaned and disinfected immediately after each clinical use. The polishers must be sterilized immediately before each use.
- Automated reprocessing in a washer-disinfector unit should generally be preferred over manual reprocessing.
- The label on the cleaning and disinfecting solutions must specifically state "suitable for rubber polishers or synthetics/silicones" (polishers are NOT included in rotary instruments for many manufacturers).
- A washer-disinfector (according to EN ISO 15883) whose effectiveness has been verified must be used to ensure reliable automated cleaning and disinfection. It is the responsibility of the operator to define, document and implement the validation, continued compliance with the performance requirements and regular routine control of the reprocessing performance of the washer-disinfector.

#### Cleaning preparation

- After use, pre-clean the polishers under running water (of drinking water quality) for 60 seconds using a brush (plastic bristles).
- Make sure to thoroughly brush the areas of the polishing head that are difficult to access.
- In case of coarse contamination additionally spray the cup-shaped polishers with a water spray gun under high pressure until the contamination is removed (approx. 60 seconds).

#### Automated cleaning and disinfection:

- Thermal disinfectant: Manufacturer's specification according to DIN EN ISO 15883. Cleaning program as indicated by the manufacturer's operating instructions.
- Place the polishers in a suitable small parts sieve or position them on a loading rack in such a way that all surfaces of the polishers are cleaned and disinfected.
- The polishers must be cleaned and disinfected in the cleaning and disinfection unit using a thermal program (temperature 90–95 °C / 194–203 °F) according to the directions of the device manufacturer.

#### Program:

##### *Holding time – disinfection*

- Pre-rinse for 5 min using cold water.
- Cleaning agent: neodisher MediClean Dental/Dr. Weigert (10 min, 55 °C / 131 °F).
- Clean for 10 min at 55 °C / 131 °F using distilled water.
- Rinse for 2 min using distilled water.
- Disinfect for e.g. 5 min at 90 °C / 194 °F (A0 value > 3000) using distilled water.
- Subsequently, dry for 15 min at up to 120 °C / 248 °F.
- Only use deionized, low-germ (max. 10 cfu/ml) and low-endotoxin (max. 0.25 EU/ml) water for the final rinse.

#### Manual cleaning and disinfection:

##### – Cleaning

- Clean using a suitable cleaning agent, e.g. ID 215/Dürr Dental:
  - Concentration 2%.
- Prepare the cleaning solution according to the manufacturer's directions and pour it into an ultrasonic bath. Completely immerse the polishers in the cleaning solution.
- Sonicate the polishers in the ultrasonic bath for 1 min.
- Remove the polishers from the cleaning solution and thoroughly rinse them (30 seconds) under running water.
- Check for cleanliness. Repeat the steps mentioned above if there is any visible residual contamination.

##### – Disinfection

- To disinfect, use an ultrasonic unit to immerse the polishers entirely in a disinfectant that is classified as suitable for rubber and silicone polishers and synthetics (e.g. ID 212/Dürr Dental) by the disinfectant manufacturer.
- Clean the polishers in the ultrasonic bath for 2 min.
- Additional exposure time (5 min) and concentration (2%) recommended by the manufacturer must be observed.
- Perform a final rinse of the polishers using distilled water (30 sec).
- Only use deionized, low-germ (max. 10 cfu/ml) and low-endotoxin (max. 0.25 EU/ml) water for the final rinse.

#### Drying

Dry the polishers with a fresh, clean, lint-free cellulose tissue.

#### Visual inspection:

- After cleaning and disinfecting, visually inspect the polishing instruments for damage and cleanliness. Repeat the cleaning and disinfecting procedure if there is evidence of any macroscopically visible residual contamination.

#### Packaging

- For packaging the polishers, a sterile barrier system (e.g. foil paper packaging, steriCLIN) according to DIN EN ISO 11607 needs to be used, which is designated by the manufacturer for steam sterilization. A suitable sealing seam system (e.g. HAWO, type 880 DC-V) is to be used.
- The polishers are double bagged.
- The bags must be big enough, so that the sealing seam is not stretched.
- Note: After the heat sealing procedure, the sealing seam must be visually checked for any defects. If a defect is found, the packaging must be re-opened and the instrument bagged and sealed again.

#### Sterilization

Polishers must be sterilized immediately before use. Ivoclar Vivadent AG recommends one of the following sterilization procedures:

- For countries outside the United States:
  - To pack the items for sterilization, only use a sterile barrier system compliant with ISO 11607-1 (e.g. steriCLIN) made of paper/film which is designated by the manufacturer for steam sterilization. The sterile barrier system being used must be large enough. The filled sterile barrier system must not be stretched.

- For the United States:

Insert the products in a perforated basket with lid and wrap in two layers of 1-ply polypropylene wrap using sequential envelope folding technique prior to sterilization. Note: Users in the United States must make sure that the sterilizer and any sterilization accessories (e.g. sterilization wraps, pouches, basket, biological or chemical indicators) are approved by the FDA for the intended sterilization.

Steam-sterilize with a fractionated pre-vacuum process according to ISO 17665 in a steam sterilizer (e.g. WEBECO A65-1) under the following conditions:

|   | Method                                                    | Conditions                 | Drying time    |
|---|-----------------------------------------------------------|----------------------------|----------------|
| 1 | Steam sterilization (autoclave)<br>Fractionated vacuum    | 134 °C (273 °F) for 5 min  | Local practice |
| 2 | Steam sterilization (autoclave)*<br>Fractionated vacuum   | 132 °C (270 °F) for 3 min  | 10 min         |
| 3 | Steam sterilization (autoclave)**<br>Fractionated vacuum  | 134 °C (273 °F) for 3 min  | Local practice |
| 4 | Steam sterilization (autoclave)***<br>Fractionated vacuum | 134 °C (273 °F) for 18 min | Local practice |

\* recommended for the USA

\*\* recommended for the UK

\*\*\* recommended for Switzerland and France

Note: Do not use a chemiclave or hot-air sterilizer.

#### Storage

Sterilized products packed in a sterile barrier system (e.g. sterilization pouch) are intended for immediate use and must not be stored for longer than 48 hours.

## 4 Safety information

- Do not use polishers that have been damaged.
- Wear protective goggles and face mask.
- Use only handpieces and contra-angles that are technically and hygienically impeccable and regularly cleaned and serviced.
- The handpieces and contra-angles used for the instruments must ensure accurate rotation.
- The polisher must be firmly secured in the handpiece.
- Jamming or leveraging should be avoided during rotation to reduce the risk of breaking the instrument.
- Do not leave the polishing instrument rotating in the same spot for too long since this may cause overheating.
- Polishers tend to vibrate if the maximum admissible rotary speed is exceeded. This may lead to damage of the polisher, deformation of the shank and/or fracture of the instrument.
- Excessive contact pressure or insufficient water cooling may lead to damage of the restoration and the surrounding tissue and it may accelerate the wear of the polishing instrument.
- Make sure to disinfect and sterilize the polishers properly.
- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: [www.ivoclar.com](http://www.ivoclar.com) and your responsible competent authority.
- The current Instructions for Use are available on the website ([www.ivoclar.com](http://www.ivoclar.com)).
- Explanation of the symbols: [www.ivoclar.com/eIFU](http://www.ivoclar.com/eIFU)

#### Warnings

- Observe the Safety Data Sheet (SDS) (available at [www.ivoclar.com](http://www.ivoclar.com)).
- Strong acids and strong alkalines may oxidize the stainless steel shank.
- Avoid temperatures >150 °C / 302 °F.
- Do not immerse the polishers for more than 1 h in disinfectant.

#### Disposal information

Remaining stocks must be disposed of according to the corresponding national legal requirements.

#### Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following clinical residual risks are known:

- Aspiration of material (extraoral and intraoral use)
- Ingestion of material (intraoral use)

## **5 Shelf life and storage**

- Do not use the product after the indicated expiry date.
- Expiry date: See note on the packages.
- Before use, visually inspect the packaging and the product for damage. In case of any doubts, please contact Ivoclar Vivadent AG or your local dealer.

## **6 Additional information**

Keep material out of the reach of children!

The product has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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