SR Ivocap® High Impact

[en] Instructions for Use

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1 Intended use

Intended purpose

- Fabrication of bases for removable dentures
- Fabrication of hard occlusal splints

Patient target group

- Patients with permanent teeth
- Adult patients with dental implants
- Edentulous adult patients

Intended users / Special training

- Dental laboratory technicians (fabrication of restorations in the dental laboratory)
 - Denturists (fabrication of restorations in the dental laboratory)
- No special training required.

Use

For dental use only.

Description

Heat-curing denture base material for the injection technique

Technical specifications

Characteristics	SR Ivocap High Impact	SR Ivocap Clear
Flexural strength	≥ 65 MPa	≥ 65 MPa
Flexural modulus	≥ 2000 MPa	≥ 2000 MPa
Residual MMA content	≤ 2.2 %	≤ 2.2 %
Water absorption (7 days)	≤ 32 µg/mm ³	≤ 32 µg/mm ³
Solubility (7 days)	≤ 1.6 µg/mm ³	≤ 1.6 µg/mm ³
Critical stress intensity factor	≥ 1.9 MPa m ^{1/2}	--
Overall work of fracture Wf	≥ 900 J/m ²	--

Indications

- Partial edentulism in the anterior and posterior region
- Complete edentulism
- Dental dysfunctions and parafunctions

Areas of application:

- Complete denture prosthetics
- Partial denture prosthetics
- Implant prosthetics
- Occlusal splints

Contraindications

The use of this product is contraindicated if the patient is known to be allergic to any of its ingredients.

Limitations of use

- Direct contact of uncured material with intraoral tissues
- The shade guide cannot be sterilized and is thus not suitable for use on the patient.
- High heat generation (exceeding 110 °C) should be avoided during the grinding, polishing and cleaning process. Avoid contact with solvents.
- The final denture or hard occlusal splint is for use by one patient only.
- Use only in conjunction with Cap Vibrator (in combination with either the SR Ivocap equipment or IvoBase Injector)

Side effects

In individual cases, allergic reactions to methyl methacrylate materials have been reported.

Interactions

There are no known interactions to date.

Clinical benefit

- Reconstruction of chewing function
- Restoration of esthetics

Composition

SR Ivocap High Impact polymer:

High impact modified polymethyl methacrylate, polymethyl methacrylate, dibenzoyl peroxide

SR Ivocap High Impact monomer:

Methyl methacrylate, ethylene glycole dimethacrylate, styrol ethylene/butylene styrol block copolymer

SR Ivocap Clear polymer:


Polymethyl methacrylate, poly(methyl methacrylate ethylacrylate), dibenzoyl peroxide, tributyl o-acetyl citrate

SR Ivocap Clear monomer:

Methyl methacrylate, ethylene glycole dimethacrylate

2 Usage

List of parts

A) Pressure apparatus ID2	C) Flask
A1) Manometer	C1) Flask half M
A2) Extrusion piston	C2) Flask half W
A3) Vise shank	C3) Flask cover
A4) Safety loop	C4) Centering insert
B) Clamping frame	D) Capsule plunger
B1) Ratchet	E) Investment aid
B2) Pressure ring with lever	F) Funnel
B3) Pressing plate	G) Timer

2.1 Fabricating the denture using the SR Ivocap equipment

Preparing the model

1. Soak the model in water for 5-10 min.
2. Isolate the model.
3. Allow the model to dry for 5 min.
4. Repeat the steps 2 to 3.

Preparing the flask

1. Isolate the inner surfaces of both flask halves using petroleum jelly.
2. Place the investment aid (E) and the flask cover (C3).

Investing the model and the waxed-up denture

1. Fill flask half M (C1) with stone.
2. Position the model 1 cm away from the flask margin.
3. Press the model down to the height of the flask margin.
4. Skim off the stone at the stone edge in such a way that it is flush with the investment aid (E).
5. Remove excess stone.
6. Allow the stone to set.
7. Remove the investment aid (E).
8. Push the funnel (F) into the centering insert (C4) until it snaps into place.
9. Place the centering insert (C4) with the inserted funnel (F).

Placing the wax components

NOTICE! The injection channels must have a diameter of 3–5 mm.

A Maxillary complete dentures:

1. Place an injection channel, fanning it out from the funnel tip towards the palate.

B Mandibular complete dentures:

2. Place two injection channels, fanning them out from the funnel tip towards the dorsal extremities of the dentures.

C Partial dentures:

3. Place an injection channel at each saddle, fanning it out from the funnel tip towards the dorsal ends.

Casting the counter model

1. Cover the model and denture teeth with a heat-resistant silicone (Shore A hardness of >80).
2. Allow the silicone to dry.
3. Thinly isolate the stone-to-stone contact surfaces.
4. Position the flask half W (C2).
5. Cover the incisal edges and occlusal surfaces with stone, ensuring there are no bubbles.
6. Place a moistened round paper insert.
7. Fill stone to the flask brim.
8. Place the flask cover (C3).
9. Press the flask cover down by hand until it snaps into place.
10. Allow the stone to set.

Boiling out and cleaning the flask

1. Immerse the flask into the water bath (90° C, 5-8 min).
2. Open the flask.
3. Remove coarse wax.

NOTICE! The flask edges and thermal insulation component must be completely free of wax.

4. Boil out the wax.
5. Clean the flasks using clean, boiling water.
6. Roughen the basal tooth surfaces using tungsten carbide burs.

Isolating the stone surfaces

1. Allow the flask halves to cool down to <30°C.
2. Isolate the surfaces of the plaster model, ensuring there is no pooling.
3. Wait for 5 min.
4. Repeat the steps 2 and 3.

Assembling the flask

1. Push the funnel (F) into the centering insert (C4) until it snaps into place.
2. Place the centering insert with the funnel (F) into the flask half M (C1).
3. Position the flask half W (C2).
4. Place the flask cover (C3).
5. Push the flask into the clamping frame (B) until it snaps into place.
6. Secure the clamping frame (B) with the inserted flask in a hydraulic dental press (3 t / 80 bar).
7. Push the ratchet (B1) to the right.
8. Remove the clamping frame (B) with the flask from the press.

Mixing the material

NOTICE! You can store mixed capsules refrigerated in the dark for up to 5 days. Allow previously mixed capsule to regain room temperature before reuse. Extend the injection time to 10 min when using previously mixed capsules.

1. Remove the monomer container from the capsule.
2. Remove the capsule seal of the polymer.
3. Twist off the sealing cap of the monomer container.
4. Pour the monomer into the open polymer capsule.
5. Screw the sealing cap on the empty monomer container.
6. Reinsert the closed monomer container into the capsule.
7. Shortly mix the monomer container.
8. Secure the capsule in the Cap Vibrator using a rubber band.
9. Mix the capsule for 5-7 min.
10. Remove the monomer container from the capsule.
11. Press the denture base material upwards with light rocking movements using the capsule plunger (D).
12. Open the capsule seal.

NOTICE! There should be no air between the plunger and the material.

13. Press the capsule plunger upwards.

Injecting the material

NOTICE! If the red mark on the plunger rod is visible, the capsule is empty. Replace it with a new, freshly mixed capsule.

1. Remove the capsule seal.
2. Push the capsule into the flask until it snaps into place.
3. Fully extend the extrusion piston (A2) of the pressure apparatus when the compressed air locking valve is open.
4. Lift the safety loop (A4).
5. Press the extrusion piston (A2) into the capsule.
6. Pull down the safety loop (A4) until it snaps into place.
7. Connect the pressure apparatus (A) with the compressed air locking valve closed to the compressed air supply (6 bar).
8. Slowly open the locking valve.
9. Leave the flask under pressure for 5 min.

Polymerizing the flask

NOTICE! Do not interrupt the boiling process by introducing other flasks.

1. Place the pressure apparatus (A) in the polymerization bath.
- NOTICE! Maintain the water level and water temperature throughout the polymerization process.**
2. Fill the polymerization bath with boiling water up to the red mark on the clamping frame.

NOTICE! Make sure that no plastic floaters are trapped beneath the clamping frame.

3. Cover the water surface with plastic floaters.
4. Polymerize the pressure apparatus (A).

Polymerization time 35 min = Residual monomer content <2.2 %

Polymerization time 90 min = Residual monomer content <1 %

Cooling the flask

1. Remove the flask with the pressure apparatus (A) using thermal gloves.
2. Allow the flask in the pressure apparatus to cool down in a water bath <15 °C, 20 min.
3. Remove the pressure apparatus from the flask.
4. Allow the clamping frame with the flask to cool down in a water bath <15 °C, 10 min.

Divesting the denture

1. Secure the clamping frame in the dental press (3 t, 80 bar).
2. Move the ratchet (B1) to the left.
3. Release the pressure.
4. Remove both flask covers (C3).

NOTICE! The injection area of the flask must be placed towards the back.

5. Secure the flask with the divesting aid on top in the dental press.
6. Separate the upper flask half from the stone by levering it upwards.
7. Remove the flask from the press.
8. Turn the flask around.
9. Repeat the steps 5 to 7.
10. Remove the stone core.
11. Divest the denture.

Finishing the denture

NOTICE! Prevent any contact of the polymerized denture with solvents or monomer.

1. Check the occlusion and articulation.
2. Finish the denture.
3. Polish the denture to a high gloss.

2.2 Fabricating the denture using SR IvoBase Injector

1.  Observe the instructions of the manufacturer.

2.3 Extending/relining/repairing the denture

1. Use cold-curing polymers.  Observe the instructions of the manufacturer.

3 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclar.com and your responsible competent authority.
- The current Instructions for Use are available on the website (www.ivoclar.com).
- Explanation of the symbols: www.ivoclar.com/eIFU

Warnings

- The monomer contains methyl methacrylate (MMA).
- Methyl methacrylate is easily flammable and irritating (flash point +10 °C).
- Irritating to eyes, skin and respiratory system. Wear protective goggles during processing. A mask should be worn and a suction device used during grinding.
- Skin contact may cause sensitization.
- Avoid contact of the skin with monomer or uncured material. Commercially available medical gloves do not offer protection against the sensitizing effect of methacrylates.
- Do not inhale vapours.
- Keep away from sources of ignition. Do not smoke.
- Do not empty into drains.
- Take precautionary measures against electrostatic discharge.
- The processing of the material can result in sharp edges, so that there is a risk of injury.
- Observe the Safety Data Sheet (SDS; can be found in the download section of the Ivoclar Vivadent AG website www.ivoclar.com).

Disposal information

Remaining stocks must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following known residual clinical risks exist:

- Ingestion of fragments
- Occurrence of denture stomatitis

4 Shelf life and storage

- Storage temperature 2-28 °C
- Do not use the product after the indicated expiry date.
- Expiry date: See note on packaging
- Protect the product from direct sunlight.
- Check to make sure that the packaging and the product are intact and undamaged before use. If in doubt, contact Ivoclar Vivadent AG or your local sales partner.

5 Additional information

Keep material out of the reach of children!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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