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1 Intended use

Intended purpose

Veneering of metal-based and zirconium oxide-based frameworks in the fabrication of indirect anterior and posterior restorations. Inlays/onlays, veneers and anterior crowns. Modification/characterization of resin teeth and denture base materials.

Patient target group

- Patients with permanent teeth
- Adult patients with dental implants
- Edentulous adult patients

Intended users / special training

- Dental laboratory technicians (fabrication of restorations in the dental laboratory)
- Dentists (clinical procedure)
- No special training required.

Use

For dental use only.

Description

SR Nexco® is a purely light-curing lab composite with micro-opal fillers for framework-based and framework-free dental restorations.

Product name	Product description
SR Nexco Liner	Protective layer between the residual tooth structure and restorative material
SR Nexco Opaquer	Covering metal frameworks, zirconium oxide frameworks and establishing the basic shade
SR Nexco Pontic Fill	Supplementing the cervical areas of pontics
SR Nexco Paste Margin	Creating composite margins on metal frameworks
SR Nexco Paste Dentin	Veneering metal frameworks and zirconium oxide frameworks; creating inlays, onlays, veneers, anterior single-tooth crowns; modifying denture teeth
SR Nexco Paste Incisal	Veneering occlusal and incisal areas on metal frameworks and zirconium oxide frameworks; inlays, onlays, veneers, anterior single-tooth crowns; modifying denture teeth
SR Nexco Paste Effect SR Nexco Paste Effect Mamelon SR Nexco Paste Effect Transpa SR Nexco Paste Effect Opal	Characterizing veneered restorations on metal frameworks and zirconium oxide frameworks; inlays, onlays, veneers, anterior single-tooth crowns; modifying denture teeth
SR Nexco Stains	Characterizing and staining restorations
SR Nexco Paste Gingiva	Veneering the gingival areas of metal frameworks and zirconium oxide frameworks; customizing complete dentures and partial dentures
SR Nexco Paste Intensive Gingiva SR Nexco Paste Basic Gingiva	Characterizing gingival areas
SR Modelling Liquid	Wetting the dental technician's instruments during modelling and as a modelling aid
SR Link	Metal/zirconium oxide/composite bonding agent that provides a covalent bond between the metal/zirconium oxide framework and SR Nexco Paste

Technical data

SR Nexco is a Type 2, Class 2 polymer-based crown and veneering material (ISO 10477:2020) with the following technical properties:

Properties	Specification			Example values ³⁾
	SR Nexco Paste	SR Nexco Opaquer	SR Nexco Liner, Stains	
Sensitivity to ambient light [s]	≥ 65			SR Nexco Paste > 420 ⁴⁾
Flexural strength [MPa]	≥ 80	N/A		SR Nexco Paste 100 ± 6
Shear bond strength [MPa] ¹⁾	≥ 5	≥ 5	N/A	SR Nexco Paste/ SR Nexco Opaquer 21.8 ± 2.6 ⁵⁾
Hardness [HV 0.2] ²⁾	≥ 18	N/A		–
Water absorption (7 days) [µg / mm ³]	≤ 40	N/A		SR Nexco Paste 14.72 ± 0.5
Solubility (7 days) [µg / mm ³]	≤ 75			SR Nexco Paste 0 ± 0.44

1) on metal frameworks or zirconium oxide frameworks, with SR Link as bonding agent

2) according to JIS T 6517

3) from Verification Report, according to ISO 10477:2020

4) from Development Report, measured at 8000 lux, according to ISO 10477: 2004

5) on Pisces Plus, with SR Link

Polymerization table**Polymerization parameters SR Nexco**

Polymerization method	Material*	PrograPrint Cure program	Curing time	Peak wavelength	Light intensity
Intermediate curing / precuring (per layer)	Opaquer Liner Dentin PonticFill Incisal Effect Margin Gingiva Stains	SR Nexco	40 s	405 nm / 460 nm	86 mW/cm ² / 68 mW/cm ²
Polymerization	Opaquer Gingiva Margin	SR Nexco	5 min.		
Final polymerization	Final restoration	SR Nexco	5 min.		

* (layer thicknesses, see "Curing depths", Chapter 2)

The product is purely light-curing and can be light cured with commercially available curing units. The curing units must have at least a light spectrum with two peaks at 405 nm (86 mW / cm²) and at 460 nm (68 mW / cm²). At this light intensity, a curing time of at least 40 seconds is required for precuring and a curing time of 5 minutes for polymerization and final polymerization. During the curing process, temperatures exceeding 50 °C should be avoided on the object.

Indications

Missing tooth structure in anterior and posterior teeth, partial edentulism in the anterior and posterior region, complete edentulism

Areas of application:

Fixed denture prosthetics

Framework-based

- Veneering of metal-supported restorations
- Veneering of zirconium oxide (ZrO₂)-supported restorations
- Veneering of combination dentures (e.g. veneering of telescope crowns)
- Veneering of fixed/removable implant superstructures
- Veneering of gingiva portions in fixed/removable implant superstructures
- Veneering of CAD/CAM-manufactured metal frameworks
- Veneering of CAD/CAM-manufactured zirconium oxide (ZrO₂) frameworks
- Masking of model cast frameworks with SR Nexco Opaquer pink

Framework-free

- Inlays/onlays/veneers
- Anterior crowns

Modification/Characterization

- Surface characterization of Ivoclar resin teeth with SR Nexco Stains in conjunction with an MMA-based light-curing conditioner for bonding light-curing materials to PMMA, followed by layering with SR Nexco Paste layering materials
- Shape and shade modifications of Ivoclar resin teeth with SR Nexco Paste layering materials in conjunction with an MMA-based light-curing conditioner for bonding light-curing materials to PMMA
- Modification and characterization of Telio® CAD with SR Nexco Stains, Dentin, Incisal and Effect materials in conjunction with an MMA-based light-curing conditioner for bonding light-curing materials to PMMA

Contraindications

The use of the product is contraindicated if the patient is known to be allergic to any of its ingredients.

Limitations of use

- Posterior crowns without framework support
- Untreated bruxism (a splint is indicated after incorporation)
- If the stipulated working technique cannot be applied
- Conventional cementation of fixed metal-free restorations
- Framework-free long-term temporaries worn for longer than 12 months
- Patients with inadequate oral hygiene and substantial medication intake (e.g. medications that reduce saliva flow)
- Repair of chipped denture teeth
- The final restoration must not be reused.

Side effects

There are no known side effects to date.

Interactions

There are no known interactions to date.

Clinical benefit

- Reconstruction of chewing function
- Restorations of esthetics

Composition

SR Nexco Pasten (Margin, Dentin, Incisal, Gingiva, Intensive Gingiva, Effect)

Copolymer, highly dispersed silicon dioxide, aromatic-aliphatic UDMA, D3MA, DCP, UDMA

Total content of inorganic fillers: 60–63 wt.%.
Particle size of the inorganic fillers: between 0.05 µm and 0.80 µm.

SR Nexco Opaquer

Aromatic-aliphatic UDMA, zirconium oxide, copolymer, D3MA, highly dispersed silicon dioxide, UDMA, iron oxide

Total content of inorganic fillers: 28–47 wt.%.
Particle size of the inorganic fillers: between 0.03 µm and 0.80 µm.

SR Nexco Liner

Barium glass, aromatic-aliphatic UDMA, D3MA, DCP, highly dispersed silicon dioxide, UDMA, hexanediol-UDMA copolymer

Total content of inorganic fillers: 56–58 wt.%.
Particle size of the inorganic fillers: between 0.03 µm and 0.80 µm.

SR Nexco Stains

Aromatic-aliphatic UDMA, highly dispersed silicon dioxide, copolymer, D3MA, DCP, UDMA, hexanediol-UDMA copolymer, rosin ester, titanium oxide

Total content of inorganic fillers: 39–44 wt.%.
Particle size of the inorganic fillers: between 0.05 µm and 0.80 µm.

SR Modelling Liquid

D3MA, EGDMA

SR Link

Butyl alcohol, ketone, UDMA, methacrylate phosphoric acid ester

2 Application

2.1 Working times / Curing depths

Working times

SR Nexco materials are sensitive to light. The working time depends on the layer thickness, shading and the existing lighting conditions. Light shades react faster than dark shades. The times listed below represent mean values at a light intensity of 3000 lux, which corresponds to the light encountered in a well-lit working space. Bear the maximum time limits in mind, when extruding the veneering material from the syringe.

SR Nexco		Time
Low viscosity	SR Nexco Liner	2-25 min
	SR Nexco Opaquer	
	SR Nexco Stains	
High viscosity	SR Nexco Margin	4-25 min
	SR Nexco Dentin	
	SR Nexco Incisal	
	SR Nexco Effect	
	SR Nexco Gingiva	

Curing depths

Due to their light sensitivity, the curing depth of the SR Nexco materials (precuring with a handheld curing light or in a curing unit) depends on the shade and above all the layer thickness. Light and translucent shades cure more readily because light can penetrate them more easily compared to darker and more opaque shades. These values must be taken into consideration during the layering of the different materials.

SR Nexco	Curing depths*
SR Nexco Opaquer	0.05 mm
SR Nexco Stains	0.2-0.8 mm
SR Nexco Paste Incisal, Dentin, Effect	2.0 mm
SR Nexco Paste Margin, Gingiva, Intensive Gingiva	1.0 mm

* Curing depths

Precuring: 40 seconds in the PrograPrint Cure / Polymerization and final polymerization: 5 minutes in the PrograPrint Cure

2.2 Framework-free restorations (inlays/onlays/anterior crowns)

Sealing the dies and adjoining parts of the model

Apply SR Model Separator in two coats. Apply the first coat slightly more generously onto the die and allow it to react for 3 minutes. Apply the second coat thinly and allow the die to dry for 3 minutes, flipped upside down. Additionally, apply SR Model Separator to adjoining model surfaces, including the opposing bite, allow to react for a short time, and then disperse excess material using oil-free compressed air.

SR Nexco Liner combination table

	BL		A					B				C				D		
Tooth shade	BL3	BL4	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Liner	BL	BL	1	2	2	3	4	1	2	3	3	1	5	5	4	5	5	5

Applying the liner

Dispense some liner paste from the syringe and spread it out slightly on a mixing pad using a disposable brush. First, apply a thin coat of liner onto the surfaces of the die. Make sure to fully cover all areas, as the liner provides an essential bonding surface to the luting composite. The liner must be applied in a thickness of at least 150 µm. Precure each segment with a handheld curing light or in a curing unit.

In the case of non-vital, discoloured cavities and preparations, coat the entire die with liner. This way, the dark shade can be blocked from the depth and, at the same time, sufficient brightness can be achieved.

Refer to the Polymerization table for the details on polymerization.

Removing the inhibition layer

Thoroughly remove the resulting inhibition layer using a disposable sponge (do not use any solvent) so that no residue remains on the liner surface. Make sure that the liner shows a slightly shiny surface.

Layering

Adapt the first layer firmly (press into place) to ensure an effective bond between the lab composite and liner surface. Generally, precure each layer as you go along, using a handheld curing light or a curing unit.

Inlay/onlay layering

Note: Apply SR Nexco Stains only in very thin layer thicknesses of 0.2 to 0.8 mm – the curing depth of dark shades is very low. Enhance the shade effect in the interdental and cavity area with Occlusal Dentin. First, build up the marginal ridges and cusps using Dentin. Apply characterizations using SR Nexco Stains, precure, and then cover with Incisal and Transpa materials. Design a true-to-nature occlusal morphology.

Anterior crown layering

Add Dentin material to the marginal areas. Enhance the shade effect in the palatal area using Stains and Occlusal Dentin materials. Outline the mesial and distal ridges using Opal Effect materials. Outline the ridges from the palatal aspect with Opal Effect 2 and cover with Incisal material. Apply Mamelon materials or Stains and precure. Complete the labial tooth contours using Incisal and Transpa materials. As a modelling aid, SR Modelling Liquid can be used on a brush to wet the instruments. Refer to the Polymerization table for the details on polymerization.

Final polymerization

After the layering procedure has been completed, ensure that all areas have been precured. Next, apply a covering layer of SR Gel to the entire veneering surface, ensuring that all areas are fully covered without the layer being too thick. Refer to the Polymerization table for the details on polymerization.

Finishing/polishing

Completely remove SR Gel using running water or a steamer and carefully remove the restoration from the die. Remove the inhibition layer and finish the surface using cross-cut burs.

Carefully smooth out the surface with rubber polishers and silicone polishing wheels. For pre-polishing and high-gloss polishing, use a goat hair brush, cotton or leather buff and Universal Polishing Paste. Use low speed and light contact pressure when polishing.

Preparing for cementation

Note: Adhesive cementation is mandatory for framework-free SR Nexco restorations.

To achieve a bond with the luting composite, the cavity side of the restoration has to be carefully blasted with Al_2O_3 (80–100 μm) at 1 bar/14.5 psi pressure in the laboratory. After the try-in in the dental practice and subsequent cleaning, the cavity side should be roughened again using a 50–100 μm diamond immediately before adhesive bonding. After that, conditioning is carried out using a bonding agent to enable chemical bonding.

2.3 Fixed metal-supported/zirconium oxide-supported restorations

Framework design

The framework must reflect the final anatomical shape in reduced form to ensure uniform support and consistent composite layer thicknesses. Create the framework according to the specified parameters provided by the alloy or zirconia manufacturer.

Alloy frameworks

Applying retention beads is advisable as they provide additional mechanical retention on top of the chemical bond with SR Link. SR Retention Adhesive is applied thinly to prevent the SR Micro Retention beads from being completely immersed in adhesive and to provide sufficient surface area for mechanical retention. If possible, create a tapered chamfer or wraparound design in the cervical area of the framework. Polish those parts of the restoration which will not be veneered before you start with the veneering process.

Conditioning the framework

Note:

- When using SR Link, do not clean the framework with steam or with an air gun after blasting!
- Do not touch the surfaces after they have been cleaned.

Compatibility of alloys / SR Link

Alloys with a gold, palladium and platinum content of less than 90%; alloys with a copper and silver content of less than 50%, base metal alloys, titanium and titanium alloys.

Conditioning alloy frameworks with SR Link

Carefully blast the framework with aluminium oxide Al_2O_3 (80–100 μm) at 2–3 bar (29–43.5 psi) pressure and tap off with an instrument (see instructions for use of the respective manufacturer). Apply SR Link bonding agent using a clean disposable brush and allow to react for 3 minutes.

Conditioning zirconium oxide frameworks (ZrO_2) with SR Link

Carefully blast the previously sintered framework using aluminium oxide Al_2O_3 (80–100 μm) at a maximum pressure of 1 bar (14.5 psi) and tap off with an instrument (see instructions for use of the respective manufacturer). Apply SR Link bonding agent using a clean disposable brush and allow to react for 3 minutes.

Applying the opaquer

1st opaquer layer

Apply the first opaquer layer (wash) thinly using a brush. Level and/or fill any roughness and retentions and then precure using the curing unit.

2nd opaquer layer

Apply the second opaquer layer in a well covering layer and precure. Then polymerize in the curing unit.

Procedure for the bridge pontic with Pontic Fill

Apply the second opaquer layer in such a way that the framework is completely covered and precure segment by segment with a handheld curing light or precure in the curing unit. Next, build up the hollowed out space on the bridge pontic to the level of the abutment teeth and precure for 40 seconds. Then, apply a layer of opaquer directly onto the inhibition layer of the Pontic Fill, precure for 40 seconds and then polymerize in the curing unit.

Refer to the Polymerization table for the details on polymerization.

Removing the inhibition layer

Thoroughly remove the resulting opaquer inhibition layer using a disposable sponge (do not use any solvent) so that no residue remains on the opaquer surface. Make sure that the opaquer shows a slightly shiny surface.

Sealing the model

Before applying the dentin and incisal layers, isolate all areas of the model which may come into contact with SR Nexco using SR Model Separator.

Cervical, dentin and incisal layering

Note: Before use, gently knead the SR Nexco Pastes on a mixing pad or mixing plate using a modelling instrument to achieve a smoother consistency. During layering, tightly adapt the SR Nexco Paste materials to prevent air trappings.

Adapt the first layer firmly (press into place) to ensure an effective bond between the lab composite and the opaquer surface. Generally, precure each layer as you go along, using a handheld curing light or a curing unit.

The layering procedure of the individual SR Nexco Paste materials is carried out either in accordance with the layering diagram (shade guide layering) or individually.

Margin material may be applied in a half-moon shape to cervical areas, pontics and crown margins that are thinning towards the metal to stabilize the shade effect. Pontic Fill is particularly suitable for pontic areas.

After that, build up the dentin layers step by step and precure each segment. The shade effect in the interdental area may be enhanced by applying chromatic materials, such as Occlusal Dentin orange. Design the dentin core in such a way that the mamelon shape remains outlined in the dentin. Make sure that enough space is left for the subsequent application of the Incisal and Transpa materials. The mamelons can be individually designed using either the Mamelon materials or SR Nexco Stains. Next, complete the restoration step by step using Incisal and Transpa materials.

As a modelling aid, SR Modelling Liquid can be used on a brush to wet the instruments.

Refer to the Polymerization table for the details on polymerization.

Final polymerization

After the layering procedure has been completed, ensure that all areas have been precured. Next, apply a covering layer of SR Gel to the entire veneering surface, ensuring that all areas are fully covered without the layer being too thick.

Refer to the Polymerization table for the details on polymerization.

Finishing/Polishing

Completely remove SR Gel using running water and/or a steamer. Remove the inhibition layer and finish the surface using cross-cut burs. Carefully smooth out the surfaces with rubber polishers and silicone polishing wheels. For prepolishing and high-gloss polishing, use a goat hair brush, cotton or leather buff and Universal Polishing Paste. Use low speed and light contact pressure when polishing.

2.4 Framework-supported combination dentures

Masking of model cast retentions with Gingiva Opaquer

Conditioning with SR Link

After finishing the model cast, carefully blast the retentions with aluminium oxide Al_2O_3 (80–100 μm) at 3 bar (43.5 psi) pressure (see instructions for use of the respective manufacturer). Remove blasting medium residue from the framework by tapping off.

Apply SR Link using a clean disposable brush and allow to react for 3 minutes.

Applying Gingiva Opaquer

1st layer of Gingiva Opaquer (wash)

Apply the first opaquer layer (wash) in a thin coat using a disposable brush. Make sure to thoroughly smooth out/fill any roughness as the wash layer represents the most important bond between the metal and composite. Subsequently, precure each segment of the wash using a handheld curing light or a curing unit.

2nd layer of Gingiva Opaquer

Apply the second opaquer layer in such a way that all retentions are adequately and completely covered and precure each segment using a handheld curing light or a curing unit. When positioning the model cast on the object holder, ensure adequate light exposure is provided (no shadow casting). Then, polymerize in the curing unit.

Refer to the Polymerization table for the details on polymerization.

Preparing for completion

After polymerization, check the curing depth with a probe. Then, remove the inhibition layer using the monomer of the respective denture base material and a disposable sponge.

It is recommended to use a cold-curing denture base resin material to secure the denture teeth to the model cast. Heat-curing may negatively affect the bond between the metal framework and the SR Nexco veneering material (opaquer).

2.5 Modifying and characterizing denture teeth, denture base materials or Telio® CAD

Denture teeth, denture base materials or Telio CAD restorations can be modified and characterized using an MMA-based light-curing conditioner to bond light-curing materials to PMMA (cold- or heat-curing polymers and resin teeth) and SR Nexco.

Conditioning the surface

After reducing the denture teeth, denture base material or Telio CAD restorations, blast the surface to be characterized using aluminium oxide Al_2O_3 (80–100 µm) at 2 bar (29 psi) pressure. Remove residue with oil-free air. Do not clean with steam! Apply the MMA-based light-curing conditioner in a thin layer, allow it to react for 2–3 minutes and then polymerize in a curing unit according to the manufacturer's instructions. Do not destroy the inhibition layer! Now, apply the SR Nexco veneering material.

Characterizing and customizing denture teeth and Telio CAD

Applying characterizations and customizations and performing shade and shape adjustments using Dentin and Effect materials. Additionally, intense characterizations can be applied using the SR Nexco Stains materials. Cover with Incisal material.

Characterizing and customizing denture base materials

Applying characterizations and customizations and performing shade and shape adjustments using Basic Gingiva B634, Gingiva and Gingiva Intensive. Additionally, intense characterizations can be applied using the SR Nexco Stains materials, which are covered with a final layer of Gingiva Paste.

Final polymerization

After the layering procedure has been completed, all layers must have been precured. Next, apply SR Gel on the entire veneering surface ensuring that all areas are fully covered without the layer being too thick. Please refer to the Polymerization table for details on polymerization.

Finishing/polishing

Completely remove SR Gel using running water and/or a steamer. Remove the inhibition layer and finish the surface using cross-cut burs. Carefully smooth out the surfaces using rubber polishers and silicone polishing wheels. For prep polishing and high-gloss polishing, use a goat hair brush, cotton or leather buff and Universal Polishing Paste. Use low speed and light contact pressure when polishing.

2.6 Restorations with gingiva portions

The SR Nexco Paste Gingiva shades are coordinated with the Ivoclar Gingiva Concept and the Gingiva Solution shade guide. In this way, natural-looking gingiva areas, particularly in conjunction with implant superstructures, can be designed using the same principles across all veneering systems.

Framework design

The framework design should be meticulously planned and fabricated by means of a wax-up and silicone keys fabricated from it. This ensures that the SR Nexco veneering material features an even layer thickness. Make sure that the soft-tissue contacts consist entirely of SR Nexco, so that SR Nexco Paste Gingiva can be used to supplement the restoration in case of later tissue recession.

Conditioning the alloy framework

Conditioning the alloy framework using SR Link

Carefully blast the framework with aluminium oxide Al_2O_3 (80–100 µm) at 2–3 bar (29–43.5 psi) pressure and tap off with an instrument (see instructions for use of the respective manufacturer). Apply SR Link bonding agent using a clean disposable brush and allow to react for 3 minutes.

Applying the opaquer and layering

Cover the dental portions of the framework with two layers of the tooth-coloured opaquer paste and precure. Polymerization is carried out in a curing unit. After removing the inhibition layer with a disposable sponge, complete the dental portions using the SR Nexco Paste materials.

Applying the opaquer for veneering the gingival portions

Applying the Gingiva Opaquer

Apply the first opaquer layer (wash) thinly using a brush and precure using a handheld curing light or a curing unit.

Completely cover the gingival area with the second opaquer layer and polymerize.

Please refer to the Polymerization table for details on polymerization.

Conditioning ceramic surfaces (e.g. IPS Style Ceram, IPS InLine, IPS e.max)

Blast the ceramic portions of the teeth to be veneered and the opacified ceramic gingiva portions of the restoration or implant superstructure using Al_2O_3 (80–100 µm) at a pressure of 2 bar (29 psi). As an alternative, IPS Ceramic Etching Gel can be applied to the ceramic portions using a disposable brush. Allow the etching gel to react for 60 seconds. After that, thoroughly rinse off the etching gel using a steamer and dry with oil-free air.

A bonding agent must be applied for conditioning (e.g. Monobond Plus, allowing a reaction time of 60 seconds, or Heliobond, light-cured with a handheld curing light or in a curing unit for 20 seconds).

Gingiva layering

First, apply Basic Gingiva B634 as the basic material to the entire area from the papilla contours to the model. To achieve an adequate depth effect, apply layers of Gingiva materials in a more intensive shade. In the process, papillae and the spaces between the alveoli can be layered in a true-to-nature fashion. To achieve a natural-looking appearance, use brighter, more translucent materials towards the top surface. Precure each individual layer segment by segment using a handheld curing light or in a curing unit.

Final polymerization

After the layering procedure has been completed, all layers must have been precured. Next, apply SR Gel on the entire veneering surface ensuring that all areas are fully covered without the layer being too thick. Please refer to the Polymerization table for details on polymerization.

Finishing/polishing

Completely remove SR Gel using running water and/or a steamer. Remove the inhibition layer and finish the surface using tungsten carbide burs.

Carefully smooth out the surfaces using rubber polishers and silicone polishing wheels. For prepolishing and high-gloss polishing, use a goat hair brush, cotton or leather buff and Universal Polishing Paste. Use low speed and light contact pressure when polishing.

2.7 Cementation

Depending of the indication, SR Nexco restorations must be seated using an adhesive, self-adhesive or conventional cementation method. Framework-free restorations must be seated adhesively. Framework-supported restorations must be seated in line with the instructions for use of the framework material being used.

3 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclar.com, and your responsible competent authority.
- The current Instructions for Use are available on the website at: www.ivoclar.com
- Explanation of symbols: www.ivoclar.com/eIFU
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>.
- Basic UDI-DI: 76152082ACOMP004LJ

Warnings

- Observe the Safety Data Sheet (SDS) (available at www.ivoclar.com).
- Contact of uncured material (pastes) with the skin, mucous membrane or eyes must be avoided. Contact with uncured material may have a slight irritating effect and may lead to a sensitization against methacrylates.
- Customary medical gloves do not provide protection against the sensitizing effect of methacrylates.
- Do not inhale grinding dust.

The safety notes on the individual primary packages and labels must be observed.

Disposal information

Remaining stocks or removed restorations must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

Some of these risks are listed below:

- Chipping / fracture / decementation of the restorative material may lead to accidental swallowing or inhalation and dental re-treatment.
- Cement excess may lead to irritation of the soft tissue / gingiva. Progressive inflammation may lead to periodontitis and bone resorption.

4 Shelf life and storage

- Storage temperature 2–28 °C/36–82 °F: SR Nexco Paste, SR Nexco Opaquer, SR Nexco Liner, SR Nexco Pontic Fill, SR Nexco Stains, SR Link, SR Modelling Liquid.
- Close syringes immediately after use. Exposure to light may lead to premature polymerization.
- Dispose of the cannula immediately after use and attach a corresponding cap to seal the syringe.
- Keep out of sunlight.
- Expiry date: See note on the packaging.
- Do not use the product after the indicated date of expiration.
- Before use, visually inspect the packaging and the product for damage. In case of any doubt, consult Ivoclar Vivadent AG or your local trade partner.

5 Additional information

Keep material out of the reach of children!

The product has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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