SR Triplex® Cold

[en] Instructions for Use

Self-curing denture base material

[de] Gebrauchsinformation

Autopolymerisierendes Prothesenbasismaterial

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1 Intended use

Intended purpose

- Fabrication of bases for removable dentures
- Fabrication of hard occlusal splints
- Fabrication of orthodontic appliances

Patient target group

- Patients with permanent teeth
- Adult patients with dental implants
- Edentulous adult patients

Intended users

- Dental laboratory technicians (fabrication of restorations in the dental laboratory)
- Denturists (fabrication of restorations in the dental laboratory)

Special training

None

Use

For dental use only.

Description

Self-curing denture base material for the pouring technique.

Technical specifications

Property	Value
Flexural strength	≥ 60 MPa
Flexural modulus	≥ 1500 MPa
Residual MMA content	≤ 4.5 %
Water absorption (7 days)	≤ 32 $\mu\text{g}/\text{mm}^3$
Solubility (7 days)	≤ 8.0 $\mu\text{g}/\text{mm}^3$

Indications

- Partial edentulism in the anterior and posterior region
- Complete edentulism
- Dental dysfunctions and parafunctions
- Malalignment of teeth

Types of restorations:

- Complete dentures
- Partial dentures
- Implant-supported dentures
- Relines
- Repairs
- Occlusal splints
- Orthodontic appliances

Contraindications

The use of this product is contraindicated if the patient is known to be allergic to any of its ingredients.

Limitations of use

- Direct contact of uncured material with intraoral tissues
- Use of the non-sterilizable shade guide on the patient
- Exceeding the maximum processing temperature of 110 °C / 230 °F
- Contact of polymerized PMMA with solvents
- Reuse of the restoration by other patients

Side effects

Allergic reactions to methacrylate in individual cases

Interactions

None known to date

Clinical benefit

- Reconstruction of chewing function
- Restoration of esthetics

Composition

Polymer: Polymethyl methacrylate, dibenzoyl peroxide
Monomer: Methyl methacrylate, ethylene glycol dimethacrylate

2 Usage

2.1 Fabricating the denture

Fabricating the plaster or silicone key

1. Mix the plaster or silicone material.  Observe the manufacturer's instructions.
2. Form the plaster or silicone key.
3. Allow the plaster or silicone key to set.
4. Remove the denture teeth from the waxed-up denture.

NOTICE! Tooth necks and the basal tooth surfaces must be completely free of wax.

5. Clean the denture teeth.
6. Roughen the basal tooth surfaces.
7. Position the denture teeth in the plaster or silicone key.
8. Carve a pouring opening into the plaster or silicone key.

Preparing the model

1. Soak the model in water for 5-10 min.
2. Isolate the model.
3. Allow the model to dry for 5 min.
4. Repeat the steps 2 to 3.
5. Secure the plaster or silicone key on the plaster model.

Mixing and pouring the material

NOTICE! Do not mix the material with other materials.

NOTICE! Observe the following processing parameters:

Mixing ratio	Dough time	Working time at 23 °C / 73 °F	Polymerisation in water bath
13 g polymer : 10 ml monomer	15 s	Max. 9 min (flow phase: 2.5-3 min transition phase: 3 min modelling phase: 3 min)	15 min 40 °C / 104 °F 2-6 bar

Table 1

1. Thoroughly mix the monomer and polymer to a homogeneous consistency for 10-15 s in a mixing cup using a spatula.
2. Allow the material to swell according to Table 1.
3. Pour the prepared material into the plaster or silicone key.
4. Polymerize the denture according to Table 1 with the pouring opening facing upwards.

Finishing the denture

NOTICE! The processing of the material can result in sharp edges so that there is a risk of injury.

1. Remove the plaster or silicone key.
2. Check the occlusion in the articulator.
3. Finish the denture.
4. Polish the denture to a high gloss.

2.2 Extending/relining/repairing the denture

NOTICE! You can also use other cold-curing polymers.  **NOTICE! Observe the instructions of the manufacturer.**

1. Roughen the contact area of the denture.
2. Wet the roughened surface with monomer.
3. Mix the material according to Table 1.
4. Process the mixed material according to Table 1.
5. Polymerize the material according to Table 1.
6. Finish the denture (See chapter "Finishing the denture")

3 Safety information

- The product has been developed solely for use in dentistry. Process according to the Instructions for Use.
- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, www.ivoclar.com and your responsible competent authority.
- Before use, visually inspect the packaging and the product for damage. In case of any doubts, please contact Ivoclar Vivadent AG or your local dealer.

Warnings

CAUTION

Product contains methacrylates.

Irritation of eyes, respiratory tract and skin
Sensitization of skin

- a) Avoid skin contact with monomer and unpolymerized material (commercially available medical gloves do not provide adequate protection).
- b) Use an extraction unit and wear protective goggles and a face mask.

CAUTION

Product contains methyl methacrylate (MMA). MMA is highly flammable.

Fire hazard (flash point + 10 °C / 50 °F)

- a) Keep away from sources of ignition.
- b) Take precautionary measures against electrostatic discharge.

Supporting documents

Document	To be found at:
Current version of Instructions for Use	www.ivoclar.com/eIFU
Structure of Instructions for Use and Warnings	www.ivoclar.com/eIFU
Explanation of symbols	www.ivoclar.com/eIFU
Safety Data Sheet (SDS)	www.ivoclar.com

Disposal information

Dispose contaminated products, remaining stocks or removed restorations according to the corresponding national legal requirements.
The product or product residues must not reach the sewage system.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following known residual clinical risks exist:

- Ingestion of fragments
- Occurrence of denture stomatitis

4 Shelf life and storage

- Storage temperature 2–28 °C (36–82 °F)
- Do not use the product after the indicated expiry date.
- Protect the product from direct sunlight.

5 Additional information

Keep material out of the reach of children!

The manufacturer accepts no liability for damages resulting from a failure to observe the Instructions for Use. Furthermore, the user is responsible for testing the products for their suitability for any purpose not explicitly stated in the Instructions.

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