IvoBase® High Impact

[en] Instructions for Use

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1 Intended use

Intended purpose

- Fabrication of bases for removable dentures

Patient target group

- Patients with permanent teeth
- Adult patients with dental implants
- Edentulous adult patients

Intended users / Special training

- Dental laboratory technicians (fabrication of restorations in the dental laboratory)
 - Denturists (fabrication of restorations in the dental laboratory)
- No special training required.

Use

For dental use only.

Description

Self-curing, impact-resistant denture base material for the injection technique.

Technical specifications

Characteristics	Value	
Flexural strength	≥ 65 MPa	
Flexural modulus	≥ 2000 MPa	
Critical stress intensity factor	≥ 1.9 MPa m ^{1/2}	
Overall work of fracture Wf	≥ 900 J/m ²	
Residual MMA content	≤ 1.5 %	≤ 1.0 using RMR*
Water absorption (7 days)	≤ 32 µg/mm ³	
Solubility (7 days)	≤ 1.6 µg/mm ³	

* RMR: Using the optional function "Residual Monomer Reduction"

Indications

- Partial edentulism in the anterior and posterior region
- Complete edentulism

Areas of application:

- Complete denture prosthetics
- Partial denture prosthetics
- Implant prosthetics
- Relines
- Repairs

Contraindications

The use of this product is contraindicated if the patient is known to be allergic to any of its ingredients.

Limitations of use

- Direct contact of uncured material with intraoral tissues
- The shade guide cannot be sterilized and is thus not suitable for use on the patient.
- High heat generation (exceeding 110 °C) should be avoided during the grinding, polishing and cleaning process. Avoid contact with solvents.
- The final denture is for use by one patient only.
- Use only in combination with IvoBase Injector (exception: repairs and relines).

Side effects

In individual cases, allergic reactions to methyl methacrylate materials have been reported.

Interactions

There are no known interactions to date.

Clinical benefit

- Reconstruction of chewing function
- Restoration of esthetics

Composition

IvoBase High Impact polymer:



High-impact-modified polymethyl methacrylate, poly(methyl methacrylate-ethylacrylate), polymethyl methacrylate, tributyl 0-acetyl citrate, dibenzoyl peroxide

IvoBase High Impact monomer:

Methyl methacrylate, 1,4-butanediol dimethacrylate

2 Application

List of parts

A) IvoBase capsule	B) IvoBase wax components	C) Spatula
A1 Capsule seal	B1 Filter wax component	D) IvoBase funnel
A2 Capsule	B2 Injection wax component	E) Aeration filter
A3 Plunger	B3 Aeration wax component	F) Separating Fluid

2.1 Fabricating the denture

Preparing the model

1. Soak the model in water for 5-10 min.
2. Isolate the model.
3. Allow the model to dry for 5 min.
4. Repeat the steps 2 to 3.

Preparing the flask

1. Isolate the inner surfaces of both flask halves using petroleum jelly.
2. Place the flask base, access former half and the filter wax component (B1).

Investing the model and the waxed-up denture

1. Fill one of the two identical flask halves with stone.
2. Position the model 1 cm away from the flask margin.
3. Press the model down to the height of the flask margin.
4. Skim off the stone at the stone edge in such a way that it is flush with the access former half.
5. Remove excess stone.
6. Allow the stone to set.
7. Remove the access former half.
8. Place the access former full.

Placing the wax components (B)

A Maxillary complete dentures:

1. Place the injection wax component (B2), fanning it out at the palatal aspect and the lateral aspect of the maxillary tuberosity.
2. Press all channels of the aeration wax component (B3) onto the anterior aspect of the model and filter wax component.

B Mandibular complete dentures:

3. Place the injection wax component (B2), fanning it out at the lingual aspect and the second molar.
4. Press all channels of the aeration wax component (B3) onto the anterior aspect of the model and filter wax component.

C Partial dentures:

5. Place the injection wax component (B2), fanning out all injection channels at the dorsal end of each denture saddle.
6. Position all channels of the aeration wax component (B3) onto the anterior stops of the denture saddles and the filter wax component (B1).

Casting the counter model

1. Cover the model and denture teeth with a heat-resistant silicone (Shore A hardness of >80).
2. Allow the silicone to dry.
3. Thinly isolate the stone-to-stone contact surfaces.
4. Position the upper flask half.
5. Close the locking clasps.
6. Fill stone into the opening of the flask.
7. Skim off excess stone so that no stone protrudes from the flask opening.
8. Allow the stone to set.

Boiling out and cleaning the flask

1. Immerse the flask into the water bath (90° C, 5-8 min).
2. Remove the locking clasps.
3. Roughly remove the wax.
4. Remove the access former full.

NOTICE! Tooth necks and the basal tooth surfaces must be completely free of wax.

5. Boil out the wax.
6. Clean the flasks using clean, boiling water.
7. Roughen the basal tooth surfaces using tungsten carbide burs.

Isolating the stone surfaces

1. Allow the flask halves to cool down to <30°C.
2. Check the flask temperature using an IR thermometer.
3. Isolate the surfaces of the plaster model without pooling.
4. Wait for 5 min.
5. Repeat the steps 3 and 4.

Assembling the flask

1. Place the aeration filter (E).
2. Push the Ivobase funnel (D) into the centring insert until it snaps into place.
3. Place the centring insert with the funnel (D) in the lower flask half.
4. Place the upper flask half.
5. Close the locking clasps.

Mixing and injecting the material

NOTICE! One capsule is sufficient for the injection of a maxillary or mandibular complete denture. With very large denture bodies, use two capsules. Fill the content of the two separately mixed capsules into one capsule. Post-injection is not possible.

1. Start the Ivobase Injector.  Observe the instructions of the manufacturer.
2. Select the program.
3. Remove the monomer container (A4) from the capsule (A).
4. Remove the capsule seal (A1) of the polymer.
5. Twist off the sealing cap of the monomer container (A4).
6. Pour the monomer into the open polymer capsule (A2).
7. Mix the monomer and polymer material using the spatula (C) to a homogeneous mixture. (20-30 s)
8. Press the flask onto the capsule. Do no longer tilt the flask.
9. Slide the flask with the capsule via the flask holder into the polymerization chamber until it stops.
10. For minimal residual monomer content: Activate the RMR key.
11. Start the program.

Cooling the flask

1. Remove the flask using thermal gloves.
2. Allow the flask to cool down in a water bath <15 °C, 15-25 min.

Divesting the denture

1. Open the locking clasps.
2. Secure the flask with the divesting aid on top in the dental press.
3. Separate the upper flask half from the stone by levering it upwards.
4. Remove the flask from the press.
5. Turn the flask around.
6. Repeat the steps 2 to 5.
7. Remove the stone core.
8. Separate the capsule from the injection channel.
9. Divest the denture.

Finishing the denture

NOTICE! Prevent any contact of the polymerized denture with solvents or monomer.

1. Check the occlusion and articulation.
2. Finish the denture.
3. Polish the denture to a high gloss.

2.2 Extending/relining/repairing the denture

NOTICE! Observe the following processing parameters:

Material	Mixing ratio	Dough phase	Processing time at 23 °C	Polymerization in the pressure pot
IvoBase Hybrid	9 g polymer : 5 g monomer	30 s	max. 10 min (pouring phase 1-2 min; modelling phase 3-5 min)	55 °C 2.5 bar 20 min
IvoBase High Impact	8 g polymer : 5 g monomer			
ProBase Cold	 NOTICE! Observe the instructions of the manufacturer.			

1. Mix the material according to the processing parameters.
2. Allow the material to react according to the processing parameters.
3. Apply the prepared material. Observe the processing time.
4. Polymerize the denture according to the processing parameters.
5. Finish the denture. (See: Finishing the denture)

3 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclar.com and your responsible competent authority.
- The current Instructions for Use are available on the website (www.ivoclar.com).
- Explanation of the symbols: www.ivoclar.com/elFU

Warnings

- The monomer contains methyl methacrylate (MMA).
- Methyl methacrylate is easily flammable and irritating (flash point +10 °C).
- Irritating to eyes, skin and respiratory system.
- Wear protective goggles during processing. A mask should be worn and a suction device used during grinding.
- Skin contact may cause sensitization.
- Avoid contact of the skin with monomer or uncured material. Commercially available medical gloves do not offer protection against the sensitizing effect of methacrylates.
- Do not inhale vapours.
- Keep away from sources of ignition. Do not smoke.
- Do not empty into drains.
- Take precautionary measures against electrostatic discharge.
- The processing of the material can result in sharp edges, so that there is a risk of injury.
- Observe the Safety Data Sheet (SDS; can be found in the download section of the Ivoclar Vivadent AG website www.ivoclar.com).

Disposal information

Remaining stocks must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following known residual clinical risks exist:

- Ingestion of fragments
- Occurrence of denture stomatitis

4 Shelf life and storage

- Storage temperature 2-28 °C
- Do not use the product after the indicated expiry date.
- Expiry date: See note on packaging
- Protect the product from direct sunlight.
- Check to make sure that the packaging and the product are intact and undamaged before use. If in doubt, contact Ivoclar Vivadent AG or your local sales partner.

5 Additional information

Keep material out of the reach of children!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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- Jei abejojate, kreipkitės į Ivoclar Vivadent AG arba vietinį prekybos partnerį.