Palodent® 360

Circumferential Matrix System

INSTRUCTIONS FOR USE ENGLISH 2

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Distributed by
Dentsply Canada
525 Rowntree Dairy Rd
Woodbridge, Ontario
L4L 7K6 Canada

Palodent® 360

Circumferential Matrix System

INSTRUCTIONS FOR USE – ENGLISH

Caution: This is a medical device. For dental use only.

USA: Rx only.

1. PRODUCT DESCRIPTION

The Palodent® 360 Circumferential Matrix System is designed to be used as a matrix to contain direct-placement restorative materials in posterior cavities involving proximal surfaces. The Palodent® 360 matrix band with integrated tightener/retainer allows matrix placement and retention without the need for bulky retainers or accessory placement and tightener instruments, allowing for ease of placement, better access, clear view of the operating field and greater patient comfort. Access for wedge placement (where indicated) is not hindered by an external retainer. The single-use, disposable Palodent® 360 system bands allow proper subgingival and proximal contact adaptation while protecting against cross-contamination.

1.1 Indication

Palodent® 360 Circumferential Matrix System is indicated where placement of a matrix band around a prepared tooth will facilitate contouring of a direct placement restorative material (e.g., Class II preparations).

1.2 Contraindications

None known.

1.3 Delivery forms

The Palodent® 360 system matrices are available in one thin gauge (thickness) in 3 heights: 4.5 mm, 5.5 mm, 6.5 mm to fit a wide range of tooth sizes and shapes.

1.4 Composition

- Matrix Bands: Polished stainless steel.
- Integrated tightener/retainer: Polyacetal Thermoplastic.

2. GENERAL SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other chapters of these Instructions for Use.

Safety alert symbol.

This is the safety alert symbol. It is used to alert you to potential personal injury hazards.

Obey all safety messages that follow this symbol to avoid possible injury.

2.1 Warnings

- Palodent® 360 System matrix band consists of small, sharp parts. As with all intraoral instruments, care must be exercised to ensure stable handling and placement. Placement of rubber dam or other suitable isolation is advised prior to operative procedures. Wear suitable protective eyewear and gloves.
 - Aspiration:** If aspiration of band or any part occurs seek immediate medical attention.
 - Ingestion:** If accidental swallowing occurs, encourage patient to seek medical assistance to locate swallowed part. If nausea or illness develop, seek medical attention immediately.
 - Laceration:** Palodent® 360 system matrix bands are sharp, and may cut unprotected skin and periodontal tissues. Avoid excessive apical pressure. Handling bands with gloved hands is recommended. In the event a clean band causes a skin laceration, employ usual first aid measures. Discard and properly dispose of band.
- Used matrix bands are contaminated. Follow proper exposure control plans. In the event of an accidental exposure, follow recommended post-exposure procedures. Following use of Palodent® 360 system matrix bands, proper handling and disposal techniques are required. Used matrix bands are contaminated and sharp. When handling or disposing, use precautions as when handling or disposing of other contaminated sharps.
- Product is provided non-sterile. Components should be sterilized prior to clinical use (see Precautions, Step-by-Step Instructions, Hygiene and Disposal).
- Care should be taken to have good control on any instrument to protect patients from injury.

2.2 Precautions

- This product is intended to be used only as specifically outlined in the Instructions for Use. Any use of this product inconsistent with the Instructions for Use is at the discretion and sole responsibility of the practitioner.
- Devices marked "single use" on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
- Do not use matrices which have not been properly pre-use processed. Before use, all components must be properly autoclaved to reduce the risk of cross contamination (see Hygiene and Disposal below).
- Avoid placing excessive pressure in an apical direction when placing or on placed bands. Apical force may result in injury to periodontal tissues (see Warnings).
- Do not use accessory tightening tools with Palodent® 360 system bands. Injury may result (see Step-by-Step Instructions).

2.3 Adverse Reactions

- Product misuse may cause failure resulting in respiratory, digestive, eye, skin, or mucosa injury (see Contraindications, Warnings and Precautions).
- Failure to properly employ accepted handling and disposal techniques may present bio-hazard (see Warnings and Precautions).

2.4 Storage

Store Palodent® 360 system bands at 2 °C-36 °C/35 °F-96 °F, away from moisture or excessive humidity. Allow components to reach room temperature prior to use. Do not freeze.

3. STEP-BY-STEP INSTRUCTIONS

3.1 Pre-operative Sterilization

Matrices are intended for single use, and are not supplied sterile. These components should be sterilized one time by steam autoclaving prior to clinical use. Paper/plastic steam sterilization pouches may be used, but are not required. Follow pre-vacuum steam autoclave manufacturer's Instructions for Use for loading and operation cycle (see Hygiene and Disposal Section).

3.2 Cavity Preparation

Complete cavity preparation as usual utilizing rubber dam or other suitable isolation technique. **Technique Tip:** Pre-wedging or pre-operative Palodent® Plus/Palodent® V3 retaining ring placement to achieve slight separation may facilitate matrix placement and acceptable restoration proximal contact. Once the preparation is complete, remove the ring, the wedge or the metal shield of the WedgeGuard (if applicable).

3.3 Matrix Selection/Preparation-Adjusting Matrix Circumference

There are three heights and one gauge of matrix band. To properly fit the tooth, matrix should be slightly larger than the occlusal circumference.

Using glove-protected hands, grasp the retainer body of a clean, new Palodent® 360 system matrix band. Twist the thumb wheel counterclockwise to unroll the coil slightly and enlarge the loop, or clockwise to take up the coil and reduce the loop, if necessary.

3.4 Positioning Matrix

Carefully slide matrix around tooth, through the interproximals using a downward oscillating motion. Matrix may be positioned with the thumb wheel tightener at the center of buccal or lingual or most secure remaining tooth wall, or in either interproximal embrasure.

Technique Tips:

- To seat matrix in tight interproximal spaces – Pre wedging or use of a Palodent® Plus or Palodent® V3 retaining ring or an Elliot Separator on posteriors may gain temporarily required space.
- To clear debris and facilitate seating – Use floss or a soft wedge between tight interproximals.
- For hard to reach molars – Use locking college pliers to grip matrix along the side and facilitate matrix seating. Release pliers when matrix is securely seated on tooth.

3.5 Tightening the Matrix

- Stabilize the matrix with gentle finger pressure apically until the matrix is fully tightened. Check to be sure gingival tissue is not trapped between band and tooth surface. If so, slightly lift and reposition band prior to tightening.
- Grasp tightening thumb wheel with thumb and forefinger. Rotate barrel clockwise.

NOTE: Palodent® 360 matrix system is designed for finger tightening ONLY. When band is snug, and resistance is felt, STOP tightening. The band is properly subgingivally tightened. **NOTE:** The integrated tightener/retainer keeps matrix locked tightly until removal is desired.

Technique Tip: Where space or grip is limited, the Palodent® 360 matrix system bands may be tightened by grasping the thumbwheel and moving the retainer body forward in a ratchet motion, releasing the thumbwheel and allowing the body to return back to center, re-grasping thumbwheel and repeating until band is tight. Reverse the direction of the retainer body movement to loosen.

Technique Tip: Palodent® 360 matrix system bands are compatible with Palodent® Plus or Palodent® V3 retaining rings to further stabilize and assist with proper embrasure form. When restoring an MOD preparation, two rings may be placed facing opposite directions.
- Place wedges as desired (buccal and/or lingual). If desired, tighter contacts may be achieved by creating additional separation using Palodent® Plus or Palodent® V3 retaining rings in combination with preferred wedges.
- Burnish contact areas.

3.6 Placement of Restorative

Place desired restorative material as per restorative material manufacturer's instructions for use.

3.7 Matrix Removal

- Following appropriate restorative curing/setting, remove retaining rings and wedges, open and remove the matrix by reversing (counterclockwise) rotation on the retainer thumb wheel until matrix is loose. This spool will usually remain attached to the band upon removal, but may become dislodged. If spool has separated from band, remove it from the mouth and properly discard.
- Support the restoration marginal ridge(s) with a condenser. Slide loosened matrix buccally or lingually, while removing matrix occlusally. Properly discard matrix. Finish contouring and verify appropriate proximal contact. **Technique Tip:** It is advisable to "post-cure" light cured restorative materials following matrix removal from all surfaces, including proximo-facial and lingual line angles.

Small parts – To reduce risk of ingestion/inhalation.

- Ensure matrix is grasped securely with fingers or locking college pliers during withdrawal from the mouth.
- Recover all pieces from the mouth immediately.

4. HYGIENE AND DISPOSAL

Cross-contamination – To reduce risk of infection.

- Do not reuse single use products. Dispose of in accordance with local regulations.
- When not in use, all supplies and instruments should be in covered storage, such as drawers and cabinets, and away from potential contamination.
- During treatment clinicians with patient contact should only handle items to be used.

NOTICE

Wrong cleaning or disinfection method.

Damage to components.

Pre-Use Cleaning – Matrices

Pre-use cleaning and disinfection may be omitted for single-use matrices stored in their original packaging and protected from contamination.

Proceed to packaging and steam autoclaving according to instructions below.

Process		✓ Acceptable
Cleaning	Mechanical (Ultrasonic)	
	Manual	
	Automated (Washer/Disinfector)	
Disinfection	Automated (Washer/Disinfector)	
	Manual	
	Immersion	
Sterilization	Steam Autoclave	✓
	Immersion	
Single Use (Matrices)	Dispose of after Use	✓
Reusable/Not Reprocessable	Dispose of if contaminated	

4.1 Sterilization of matrices pre-use should be accomplished (one time) as outlined below:

Instructions for sterilizing Palodent® 360 Matrices	
Warnings	- High-level disinfection has not been validated as a terminal process for the Palodent® 360 Matrix. Disinfection is not recommended as a terminal process. Proceed to Sterilization following any optional cleaning/intermediate disinfection process. - Steam autoclaving sterilization is appropriate and recommended for the Palodent® 360 Matrix. - Follow proper infection prevention activities, such as proper hand washing and donning new puncture and chemical gloves at pertinent steps. - Do not allow the Palodent® 360 matrices to exceed 137 °C. - Do not use phenol-based glutaraldehyde. - Sterilize metal and plastic components in separate pouches to avoid damage to components. - If the integrity of the pouch has been compromised, the device should be processed again before use. - Matrix should be sterilized before use. Cleaning and Disinfection of new, uncontaminated matrices is not required prior to sterilization. Cleaning and disinfection have not been validated.
Limitations on Reprocessing	Palodent® 360 Matrix is single use only, and may be processed one time before use. Do not reprocess for re-use. Cold liquid immersion disinfection/sterilization, automated washer/disinfection, manual or mechanical cleaning or disinfection have not been tested or validated for efficacy and are not recommended for use.
Initial Treatment at the Point of Use	No particular requirements if device has been stored: - in its original packaging. - in covered drawer or cabinet. - outside the patient treatment operator. - where contamination is avoided. Contaminated matrices must be disposed of according to local and national regulations. Do not use.
Preparation before cleaning	N/A
Cleaning and Disinfection: Mechanical, Manual, Automated	No validated process.
Maintenance, Inspection and Testing	- Visually inspect the Palodent® 360 Matrix prior to use. If the device is contaminated, discolored, cracked, worn, distorted, etc., discard and do not use. - Ensure moving parts operate smoothly.
Packaging	- Place each Palodent® 360 Matrix in a separate appropriately sized steam-sterilization pouch. - Uses of FDA-cleared, ISO 11607 compliant Paper/plastic steam sterilization pouches (e.g., AssurePlus® Sterilization Pouches) are recommended.
Sterilization	The Palodent® 360 matrices may be sterilized by Dynamic-Air-Removal/Pre-Vacuum or Gravity Displacement steam autoclave. The following Pre-Vacuum cycles may be used for Palodent® 360 matrices: 132 °C (270 °F) for 4 minutes with a minimum drying time of 20 minutes 134 °C (273 °F) for 3 minutes with a minimum drying time of 18 minutes 135 °C (275 °F) for 3 minutes with a minimum drying time of 16 minutes The following Gravity cycles may be used for Palodent® 360 matrices: 132 °C (270 °F) for 15 minutes with a minimum drying time of 30 minutes 135 °C (275 °F) for 10 minutes with a minimum drying time of 30 minutes Follow manufacturer's instructions for loading and operation cycle. An alternate method of sterilization is to place non-bagged instruments into the steam sterilizer and run at one of the above listed cycles.
Storage	- Do not remove from the pouch until ready for use to prevent contamination. - To prevent contamination, store processed device outside of the patient treatment operator and/or in covered storage, such as drawers and cabinets, protected from contaminated hands, spatter, spray or aerosol of body fluids. - Prior to use, inspect the sterilization pouch and the matrix. If the integrity of the pouch has been compromised, the matrix should be disposed of in accordance with local regulations. - Matrices steam autoclaved unwrapped should be used immediately.

Additional Information	The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for use. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and trained personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process. Likewise, any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse consequences. Use of other reprocessing procedures/methods is at the discretion and sole responsibility of the practitioner.
Manufacturer Contact	Within the United States, call Dentsply Sirona at 1-302-422-4511. For areas outside the United States, contact your local Dentsply Sirona representative.

4.2 Disposal

- Dispose of in accordance with national and local regulations.
- During incineration, matrix retainer body may produce irritating and highly toxic gasses.
- Used matrices are contaminated and sharp. When handling or disposing, use precautions as when handling or disposing of other contaminated sharps.

5. LOT NUMBER AND CORRESPONDENCE

- The following numbers should be quoted in all correspondence:
 - Reorder number
 - Lot number on the package
 - Copy of proof of purchase
- Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.
- Device Identification (Basic UDI-DI): ++D002MATRIXSYSTEMSPX

© 2025 Dentsply Sirona Inc. All Rights Reserved. **559600WEB** (R 1/9/25)

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