AutoMatrix®

Retainerless Matrix System

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AutoMatrix® Retainerless Matrix System

For use with AutoMate® III Tightening Device and AutoMatrix® Snippers+

INSTRUCTIONS FOR USE – ENGLISH

CAUTION: This is a medical device. For dental use only by dental professionals.
USA: Rx only.

1. PRODUCT DESCRIPTION

The AutoMatrix® Retainerless Matrix System is designed to be used as a matrix to contain direct-placement restorative materials in posterior cavities involving proximal surfaces. The AutoMatrix® system allows matrix placement and retention without the need for bulky retainers, allowing for ease of placement, better access, clear view of the operating field and greater patient comfort. Access for wedge placement (where indicated) is not hindered by an external retainer. The single-use, disposable AutoMatrix® system bands allow proper subgingival and proximal contact adaptation while protecting against cross-contamination.

1.1 Indications

AutoMatrix® Retainerless Matrix System is indicated where placement of a matrix band around a prepared tooth will facilitate contouring of a direct placement restorative material (e.g., Class II preparations).

1.2 Contraindications

None known.

1.3 Delivery forms¹

The AutoMatrix® system matrices are available in four height/width and gauge (thickness) configurations (see table below). All AutoMatrix® system matrix sizes fit teeth of any circumference.

Type	Medium Thin ²	Narrow Regular	Medium Regular	Wide Regular
Height & Gauge	6.35 mm × 0.038 mm (1/4" × .0015)	4.76 mm × 0.051 mm (3/16" × .002)	6.35 mm × 0.051 mm (1/4" × .002)	3.36 mm × 0.051 mm (5/16" × .002)
Typical Clinical Uses	Bicuspid/molars	Short Bicuspid and molars	Bicuspid/molars	Long molars

¹ Some delivery forms may not be available in all countries.

² Thin-gauge bands identified by notch on auto-lock loop.

1.4 Composition

- Matrix Bands: Polished stainless steel
- AutoMate® III tightening device: Stainless steel
- AutoMatrix® Snippers+: Stainless steel, with plastic catch-tray

2. GENERAL SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other chapters of these Instructions for Use.

Safety alert symbol.

This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury.

2.1 Warnings

- AutoMatrix® System matrix bands contain Nickel (metallic) which may be irritating to skin and oral mucosa and may cause allergic contact dermatitis or mucositis in susceptible persons.
 - Avoid skin contact** to prevent irritation and possible allergic response. In case of contact, red-dish rashes may be seen on the skin. In case of skin or mucosa sensitization or rash, discontinue use and seek medical attention.
- AutoMatrix® System matrix band consists of small, sharp parts. As with all intraoral instruments, care must be exercised to ensure stable handling and placement. During matrix removal procedures a segment of the AUTO-LOCK LOOP (retention clip) may separate from the matrix band. Prevent this segment from being swallowed or aspirated. Remove any separated piece immediately from the mouth. Placement of rubber dam or other suitable isolation is advised prior to operative procedures. Wear suitable protective eyewear and gloves.

Aspiration: If aspiration of matrix band occurs seek immediate medical attention.

Ingestion: If accidental swallowing occurs, encourage patient to seek medical assistance to locate swallowed part. If nausea or illness develop, seek medical attention immediately.

Laceration: AutoMatrix® system matrix bands are sharp, and may cut unprotected skin and periodontal tissues. Avoid excessive apical pressure. Handling bands with gloved hands and/or locking college pliers or similar instrument is recommended. In the event a clean band causes a skin laceration, employ usual first aid measures. Discard and properly dispose of band.
- Used matrices and accessories are contaminated. Follow proper exposure control plans. In the event of an accidental exposure, follow recommended post-exposure procedures. Following use of AutoMatrix® system matrix bands, proper handling and disposal techniques are required. Used matrices are contaminated and sharp. When handling or disposing, use precautions as when handling or disposing of other contaminated sharps.
- Product is provided non-sterile. Components should be sterilized prior to clinical use (see Precautions, Step-by-Step Instructions, Hygiene).
- Care should be taken to have good control on any instrument to protect patients from injury.
- Do not use excessive force when inserting AutoMate® III tightening device into tightening coil. Do not allow rotating coil to contact oral tissues during tightening. When using removal snippers, employ a technique that prevents entrapping soft tissues (skin, gingiva, mucosa) within the snipper jaws which may result in injury to the tissue. If injury occurs, thoroughly wash the affected area with water and seek appropriate medical attention.
- AutoMatrix® system should not be used with patients who have a known hypersensitivity to any of the components.

2.2 Precautions

- This product is intended to be used only as specifically outlined in the Instructions for Use. Any use of this product inconsistent with the Instructions for Use is at the discretion and sole responsibility of the practitioner.
- Wear suitable protective eyewear, mask, clothing and gloves. Protective eyewear is recommended for patients.
- Excessive pressure on matrix in an apical direction may result in injury to periodontal tissues (see Warnings).
- Devices marked "single use" on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent crosscontamination.
- Use only clean and properly disinfected and sterilized, or ready for use accessories. Before use and each re-use, accessory tools must be properly cleaned, disinfected and autoclaved to reduce the risk of cross contamination (see Hygiene below).
- Matrices are provided ready for use but should be sterilized prior to use (see Hygiene below). Use only new, clean unused matrices, stored in covered storage, such as drawers and cabinets, away from potential contamination, and removed from their original packaging using either clean gloved hands or a clean suitable instrument at time of use. To reduce risk of infection, do not use matrices if packaging is damaged or compromised.
- Do not use accessory tools other than those designed, supplied, and/or outlined for use with the AutoMatrix® System. Injury may result (see Warnings, Step-by-Step Instructions).
 - Always make sure that the plastic-basket is properly connected to the AutoMatrix® Snippers+ before application.
 - When clipping AUTO-LOCK LOOP with removal snippers, shield open end with fingers to prevent clipped part from projecting into the oral cavity. Immediately retrieve any loose part and remove from the mouth (see Warnings).

2.3 Adverse reactions/Undesirable side effects

- Product may irritate the skin and oral mucosa. Skin contact: Irritation or possible allergic response. Reddish rashes may be seen on the skin. Mucous Membranes: Inflammation, sloughing or possible allergic response (see Warnings).
- Product misuse may cause failure resulting in respiratory, digestive, eye, skin or mucosa injury (see Warnings and Precautions).
- Failure to properly employ accepted handling and disposal techniques may present biohazard (see Warnings and Precautions).

2.4 Storage

Store AutoMatrix® system bands and accessory AutoMate® III tightening device and AutoMatrix® Snippers+ at room temperature, away from moisture or excessive humidity. Tightening device and snippers should be cleaned, disinfected/sterilized and dried prior to storage (see Hygiene).

3. STEP-BY-STEP INSTRUCTIONS

3.1 Pre-operative sterilization

Matrices are intended for single use and are not supplied sterile. These components should be steam autoclaved one time prior to clinical use. Tightening device and snippers are not supplied sterile and should be steam autoclaved prior to first and each subsequent use (see Hygiene section).

3.2 Place rubber dam or other suitable isolation technique.

3.3 Cavity preparation

Complete cavity preparation as usual utilizing rubber dam or other suitable isolation technique. Technique Tip: Pre-wedging, or placement of a Palodent® Plus ring (available separately, see complete Instructions for Use) to achieve slight separation may facilitate matrix placement and acceptable restoration proximal contact.

3.4 Matrix selection and preparation – Adjusting matrix circumference

Select the appropriate height and gauge matrix band. To properly fit the tooth, matrix should be slightly larger than the occlusal circumference. The height should extend slightly beyond the marginal ridge, but not obstruct the procedure.

- Using protected hands, or suitable instrument (e.g., college pliers), unroll the coil of a clean, new AutoMatrix® system matrix band slightly.
- While securely holding matrix band (with instrument or by side; do not hold by sharp edges) carefully slide the AUTO-LOCK against the coil. If necessary, repeat steps 3.4.1 and 3.4.2 above to further enlarge matrix to the appropriate size.

3.5 Positioning matrix

Carefully slide matrix around tooth, through the interproximals. Position matrix with the AUTO-LOCK LOOP at the center of buccal or lingual or most secure remaining tooth wall.

Technique Tips:

- To seat matrix in tight interproximal spaces – pre wedging, placement of a Palodent® Plus ring or use of an Elliot Separator on posteriors may gain temporarily required space.
- To clear debris and facilitate seating – use waxed floss or a soft wedge between tight interproximals.
- For hard to reach molars – use locking college pliers to grip matrix along the side and facilitate matrix seating. Release pliers when matrix is securely seated on tooth.

3.6 Tightening the matrix (using the AutoMate® III tightening device)

- Stabilize the matrix with gentle finger pressure apically until the matrix is fully tightened. Check to be sure gingival tissue is not trapped between band and tooth surface. If so, slightly lift and reposition band prior to tightening. Ensure matrix is apical to preparation margin with no gap.
- Insert AutoMate® III tightening device tip into the matrix coil. DO NOT USE FORCE.
- Grasp AutoMate® III tightening device with thumb and forefinger on barrel, palm grip in palm. Rotate knurled barrel clockwise. **NOTE:** AutoMate® III tightening device has an audible torque indicator, thereby indicating the rotational (tightening) force required. After 2 to 3 "clicks", STOP tightening. The band is properly subgingivally tightened. **NOTE:** The AUTO-LOCK LOOP keeps matrix locked tightly until removal is desired.
- Remove the tightening device. Rotate AutoMate® III tightening device continuously counter-clockwise (backwards) while simultaneously and continuously withdrawing the AutoMate® III tightening device from the coil.

3.7 Placement of restorative

- Close the gingival margin by placing a contoured wedge or anatomic wooden wedge (buccal and/or lingual).
- Ensure matrix is adapted to the adjacent tooth contact area. Lightly burnish the matrix against the adjacent tooth both buccally and lingually to further adapt it to form proper contours. **NOTE:** Over burnishing may distort matrix and compromise contours.
- Place desired restorative material as per restorative material manufacturer's directions for use.

3.8 Removing the matrix

Following appropriate restorative curing/setting, open and remove the matrix using shielded AutoMatrix® Snippers+.

3.8.1 Preparation of AutoMatrix® Snippers+:

- Align catch-tray with the instrument axis of the snipper.
- Clip catch-tray straight onto the snippers by sliding the u-shaped cavity into the nut below the bolt head.
- Make sure that the pin on the lower surface of the catch-tray clicks in the hole on the upper side of the blade edge.

Displacement of catch-tray. Small parts – To reduce risk of ingestion/inhalation.

- Do not open Snippers+ more than 60° to prevent displacement of catch-tray, which may cause leakage of the catch-tray and release of cut pieces.
- Before application, make sure that the catch-tray is properly connected to the snippers. Turn instrument and check that there is no gap between catch-tray and snippers. Make sure that the catch-tray is empty.

3.8.2 Clip any part of the projecting end of the AUTO-LOCK LOOP.

3.8.3 After clipping immediately retrieve any loose part and remove from the mouth.

3.8.4 Keep snippers closed containing the cut piece during withdrawal from the mouth.

3.8.5 Empty catch-tray and properly discard cut piece.

Small parts – To reduce risk of ingestion/inhalation.

- Do not open Snippers+ widely during use. Widest opening is designed to dislodge the catch-tray.
- After clipping the AUTO-LOCK LOOP into the catch-tray, keep Snippers+ closed during withdrawal from the mouth.
- If tray is dislodged during use, recover any pieces from the mouth immediately.
- Ensure matrix is grasped securely with fingers or locking college pliers during withdrawal from the mouth.
- Recover all pieces from the mouth immediately.

3.8.6 Insert an explorer, scaler or suitable instrument through the lockrelease hole located on the band near the coil. Gently pull lock-release hole occlusally until the matrix slips open.

NOTE: Locate the gingival segment of the AUTO-LOCK LOOP. This part will usually remain attached to the band upon removal but may become dislodged. If segment of AUTO-LOCK LOOP has separated from band, remove it from the mouth and properly discard.

3.8.7 Support the restoration marginal ridge(s) with a condenser. Slide loosened matrix buccally or lingually, while removing matrix occlusally. Properly discard matrix.

4. HYGIENE, PROCESSING AND DISPOSAL

Cross-contamination – To reduce risk of infection.

- Do not reuse single use products. Dispose of in accordance with local regulations.
- Reprocess reusable products as described below.
- When not in use, all supplies and instruments should be in covered storage, such as drawers and cabinets, and away from potential contamination.
- During treatment clinicians with patient contact should only handle items to be used.
- Dispensing materials with clean/disinfected gloves or hands or suitable transfer forceps in a separate room, bringing only what will be used into the operatory is strongly recommended.

NOTICE

Wrong cleaning or disinfection method.

Damage to components.

- Always disassemble Snippers+ catch-tray before final reprocessing.
- Do not clean or disinfect with aggressive agents (e.g., solutions based on orange oil or acetone).

For manual cleaning and disinfection:

- Discard used gloves according to local regulations.
- Disinfect hands with an appropriate bactericidal, virucidal, and fungicidal hand disinfectant solution according to local regulations. Use according to disinfectant solution manufacturer's Instructions for Use.
- Use a new clean pair of examination gloves.

Process		✓ Acceptable
Cleaning	Mechanical (Ultrasonic) Excluding AutoMate® III Tightening Device	✓
	Manual	✓
	Automated (Washer/Disinfector)	✓
Disinfection	Automated (Washer/Disinfector)	✓
	Manual Excluding AutoMate® III Tightening Device	✓
	Immersion	
Sterilization	Steam Autoclave	✓
	Immersion	
Single Use (Matrices)	Dispose of after Use	✓
Reusable/Not reprocessible	Dispose of if contaminated	

Instructions for cleaning, disinfecting and sterilizing AutoMate® III Tightening Device, AutoMatrix® Snippers+ and catch-tray, AutoMatrix® matrices		
Instructions for processing AutoMate® III Tightening Device ("X" applicable)		
Instructions for processing AutoMatrix® Snippers+ and catch-tray, AutoMatrix® matrices ("X" applicable)		
Warnings	- High-level disinfection has not been validated as a terminal process for the AutoMate® III Tightening Device, AutoMatrix® Snippers+ and catch-tray or AutoMatrix® matrices. Disinfection is not recommended as a terminal process. Proceed to Sterilization following cleaning and any optional/intermediate disinfection process. - Steam autoclaving sterilization is appropriate and recommended for the AutoMate® III Tightening Device, AutoMatrix® Snippers+ and catch-tray and AutoMatrix® matrices. - Follow proper infection prevention activities, such as proper hand washing and donning new puncture and chemical resistant gloves at pertinent steps. - Do not use hard wire brushes to clean the Instruments as scratching and mechanical degradation can occur. - Use only a disinfecting solution which is approved for its efficacy, EPA registered (and/or Health Canada approved) and use in accordance with the IFU of the disinfecting solution manufacturer. - Do not allow the device to exceed 137 °C (279 °F). - Do not use phenol-based glutaraldehyde. - Always use a pH neutral instrument cleaning solution if authorized. - Sterilize metal and plastic components in separate pouches to avoid damage to components. - If the integrity of the pouch has been compromised, the device should be processed again before use. - In the USA, all devices should be cleaned and sterilized before each use. For the countries that require a three-step process, in which disinfection is required prior to sterilization, all devices must be cleaned, disinfected, and sterilized before use.	X X
Warnings	- Do not mechanically clean and/or disinfect AutoMate® III Tightening Device in ultrasonic cleaner. - Do not immerse AutoMate® III Tightening Device in liquid. Cleaning/disinfection using static immersion and/or an ultrasonic unit is not validated and should not be used.	X
Limitations on reprocessing	- The AutoMatrix® Retainerless Matrix System matrix band is intended for single use only. Do not attempt to reprocess and/or re-use. - Repeated reprocessing has minimum effect on the accessory instruments. End of life is normally determined by wear and damage due to use. - Cold liquid immersion disinfection/sterilization, chemical vapor sterilization, and dry heat sterilization methods have not been tested or validated for efficacy and are not recommended for use.	X X
Initial treatment at the point of use	- Do not allow residue or any form of contamination to dry on the device. Remove excess soil/gross debris by wiping, brushing, and/or rinsing under water. - Transport to processing area and follow validated processing instructions below. - It is recommended that the device be processed as soon as is reasonably practical following use. - Start reprocessing within 1 hour after use.	X X
Preparation before cleaning	- Remove the catch-tray from the AutoMatrix® Snippers+: 1. For loosening of catch-tray, open AutoMatrix® Snippers+ as widely as possible. 2. Pull off catch-tray from the AutoMatrix® Snippers+. - Do not attempt to disassemble the AutoMate® III Tightening Device, AutoMatrix® Snippers+ or AutoMatrix® Retainerless Matrix System matrix band.	X X
Cleaning: Mechanical (Ultrasonic)	- Do not mechanically clean AutoMate® III Tightening Device in ultrasonic cleaner. No mechanical cleaning process has been validated. See "Manual cleaning" or "Automated cleaning and disinfection" below. For AutoMatrix® Snippers+ and catch-tray, AutoMatrix® matrices: - The chemical manufacturer's instructions for the cleaning and disinfection solution should be followed, observing concentration rates, contact times, fill level, degassing processes, etc. - The recommended solution is Resurge® Instrument Cleaning Solution mixed one packet 14.8 ml (1/2 fl. oz.) with 3.8 l (1 gallon) of warm water, at 35 °C-40 °C (95 °F-104 °F) or equivalent enzymatic cleaning solution. - Immerse device in an ultrasonic bath containing the prepared solution for at least 20 minutes unless otherwise specified by the time specified by the solution manufacturer. - Rinse under tap water of at least drinking quality for a minimum of 30 seconds. - Allow to air dry. - Visually inspect for visible soil, and repeat the cleaning and disinfection procedure if necessary.	X
Cleaning: Manual	- Clean instruments and components with a soft brush and tap water until free of visible contamination. Pay special attention to device seams and insertions. - The recommended solution for manual cleaning is Resurge® Instrument Cleaning Solution mixed one packet 14.8 ml (1/2 fl. oz.) with 3.8 l (1 gallon) of warm water, at 35 °C-40 °C (95 °F-104 °F) or equivalent enzymatic cleaning solution. - The recommended solution for manual cleaning of the AutoMate® III Tightening Device in EU countries is Dürr Dental ID215 Instrument Cleaning Solution mixed 20 ml with 1 l or equivalent enzymatic cleaning solution. - Scrub with a lint-free disposable wipe pre-impregnated with instrument cleaning solution ensuring that the fluid penetrated all seams and crevices. - Wrap in a second lint-free disposable wipe pre-impregnated with instrument cleaning solution. Keep the surface of the instrument wet for 5 minutes with Resurge®, 15 minutes with 2% Dürr Dental ID215. - Remove cleaning solution residue by rinsing under warm tap water for 2 minutes. - Dry instruments and components with a lint-free, single-use cloth or allow to air dry in a clean area. - Visually inspect for visible soil, and repeat the cleaning procedure if necessary.	X X

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© 2026 Dentsply Sirona Inc. Tous droits réservés. **563030WEB** (R 6/19/26)
Dentsply DeTrey GmbH