KerrTM

HaweTM

Centring Rings | X-Ray Holders

X-Ray Centring Ring

**X-Ray Film and Phosphor Plate Holders
and X-Ray Sensor Holders**

HaweTM

X-Ray Centring Devices

Centring Devices

IFUs

Instructions for Hawe™ Centring Rings, X-Ray Holders and Hawe™ X-Ray Centring Devices

SYMBOLS DESCRIPTION Full explanation of symbols used on Kerr packaging located at: <http://www.kerrdental.com/symbols-glossary>

	Catalogue number
	Manufacturer
	Use-by date
	European conformity (CE) mark for Class I medical devices
	Consult instructions for use
	Caution
	Authorized representative in the European Community
	Batch code
 for US law only	Caution: Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner. 1. for US law only
 135°C 275°F	2. Sterilizable in a steam sterilizer (autoclave) at temperature specified
	3. Ultrasonic bath
	4. Washer-disinfector for thermal disinfection
	5. Medical Device
	Importer

en - ENGLISH

Hawe™ Centring Rings | X-Ray Holders

X-Ray Centring Ring | X-Ray Film and Phosphor Plate Holders and X-Ray Sensor Holders.

General Information

X-Ray Film and Phosphor Plate Holders and X-Ray Sensor Holders are available with index or ring and intended to hold films, phosphor plates or digital sensors respectively. Together with Centring Device, they facilitate the adjustment of the X-ray tube positioning in order to get the radiography right the first time, reducing the exposure of the patient to radiation.

The X-Ray Holders are available for paralleling techniques, as well as for endodontic imaging where endodontic files and instruments may otherwise impede proper positioning of the receptor behind the tooth.

The X-Rays Holders are made of plastic materials.

Indications for Use

To assist the direct or indirect imaging of the bitewing, periapical or endodontic radiographs for human diagnosis purposes.

Contraindications

None known.

Warnings

None known.

Precautions

- These products are intended to be used by qualified dental professionals on the general patient population.
- Wear personal protective equipment when using and handling these products.
- This device is sold non-sterile and needs to be processed prior to first use according to section "Cleaning, Disinfection, and Sterilization".
- Discard the device if corrosion, deformation, breakages, discoloring or missing parts are noticed.
- To avoid cross contamination, this product must be reprocessed after each patient following this instruction for Use.

Adverse Events

If a serious incident occurs with this medical device, report it to the manufacturer and to the competent authority for that country in which the user and/or patient is established.

Step-by-Step Instructions

The X-Ray Film and Phosphor Plate Holders and X-Ray Sensor Holders can be used without or in conjunction with the X-Ray Centring Devices System as per info cards:


- Paro-Bite and Kwik-Bite (p. 86-89)
- Super-Bite Anterior and Super-Bite Posterior (p. 90-93)
- Endo-Bite Anterior and Endo-Bite Posterior (p. 94-97)
- Kwik-Bite Senso (p. 98-101)
- Super-Bite Senso Anterior and Super-Bite Senso Posterior (p. 102-105)
- Endo-Bite Senso Anterior and Endo-Bite Senso Posterior (p. 106-109)

en

Cleaning, Disinfection, and Sterilization

Instructions for Cleaning / Disinfecting / Sterilizing X-Ray Holders and Hawe Centring Rings

WARNINGS	- Start the reprocessing procedure immediately after use because dried contaminants may negatively influence the cleaning effort. - During sterilization, avoid direct contact with the metallic part of the autoclave as the temperature there can exceed the pre-set temperature.
Limitations on processing	Do not use chemical autoclave or cold sterilization solutions.
Initial treatment at the point of use	- Remove the Film / Phosphor Plate / Sensor. - Remove centring ring from the holder stem. - Remove the centring device.
Preparation before cleaning	Twist holder plate into the loose intermediate position (Super-Bite, Endo-Bite) or open position (Kwik-Bite Senso). 
Cleaning and Disinfection: Automated³	Use only a properly maintained, inspected, calibrated and approved washer-disinfector in compliance with ISO 15883-1. Perform a cleaning/disinfection program cycle: 4 min. pre-wash with cold water; 10 min. wash with alkaline detergent at 55°C ±2°C (131±3°F); 1 min. enzymatic neutralization in cold tap water, if recommended by the manufacturer of the cleaning agent used; 3 min. rinse in cold water; 5 min. thermal disinfection at 93°C (199°F) with demineralized water; 35 min. drying at 100°C (212°F).
Cleaning: Manual	Remove any visible debris or organic matter from the instrument by scrubbing in warm tap water (> 40°C-104°F) until no visible debris is observed. Never use metal brushes or wire wool.

Cleaning: Manual (Continued)	• With the naked eye, inspect instruments for debris. Repeat previous step if there is still debris. • Thoroughly scrub the instrument surfaces, for at least 1 minute, with a soft brush dampened with an enzymatic solution¹ following the manufacturer's instructions with regard to the concentration and temperature. Do not use detergents containing strong alkalis, strong acids, hydrogen peroxide or organic solvents. • Soak the instrument in an enzymatic solution following the manufacturer's instructions¹ with regard to the concentration, temperature and duration. Do not use detergents containing strong alkalis, strong acids, hydrogen peroxide or organic solvents. • Alternatively to the previous step, place the instrument in an ultrasonic cleaner with a cleaning or enzymatic agent following manufacturer's instructions¹ with regard to the concentration, temperature and duration.																				
Disinfection: Manual	• Rinse all parts for at least 30 sec with warm tap water (> 40°C -104°F). • Soak the instrument in a disinfectant solution following the manufacturer's instructions² with regard to the concentration, temperature and duration. Do not use disinfectant containing strong alkalis, strong acids, hydrogen peroxide or organic solvents. • Alternatively to the previous step, place the instrument in an ultrasonic cleaner with disinfecting agent following manufacturer's instructions² with regard to the concentration, temperature and duration.																				
Drying	Dry with clean debris-free soft cloth.																				
Maintenance, Inspection and Testing	Inspect and discard the instrument if there is evidence of plastic deformations or cracks, bent parts, brittle material, corrosion or missing parts. The examination can be made with the naked eye or with a magnifier under normal lighting conditions.																				
Packaging	Paper / plastic steam sterilization pouches may be used.																				
Sterilization	Place the instrument in a sterilization pouch suitable for steam sterilization following the instructions of the pouch manufacturer. Pay attention to the tip of the instruments as it may damage the packaging during handling, distribution or storage. Proper packaging protects the sterilized instruments from bacterial recontamination during handling and storage. Use a steam autoclave with distilled water. Do not exceed 135°C (275°F). Autoclave Parameters <table><tr><th>Region</th><th>Autoclave Type</th><th>Cycle Temperature</th><th>Exposure time</th><th>Minimum Dry Time</th></tr><tr><td>EU</td><td>Pre-vacuum</td><td>134°C (273°F)</td><td>3 minutes</td><td>14 minutes</td></tr><tr><td>US</td><td>Gravity (Wrapped)</td><td>121°C (250°F)</td><td>30 minutes</td><td>15-30 minutes</td></tr><tr><td>US</td><td>Pre-Vacuum (Wrapped)</td><td>135°C (275°F)</td><td>3 minutes</td><td>16 minutes</td></tr></table>	Region	Autoclave Type	Cycle Temperature	Exposure time	Minimum Dry Time	EU	Pre-vacuum	134°C (273°F)	3 minutes	14 minutes	US	Gravity (Wrapped)	121°C (250°F)	30 minutes	15-30 minutes	US	Pre-Vacuum (Wrapped)	135°C (275°F)	3 minutes	16 minutes
Region	Autoclave Type	Cycle Temperature	Exposure time	Minimum Dry Time																	
EU	Pre-vacuum	134°C (273°F)	3 minutes	14 minutes																	
US	Gravity (Wrapped)	121°C (250°F)	30 minutes	15-30 minutes																	
US	Pre-Vacuum (Wrapped)	135°C (275°F)	3 minutes	16 minutes																	

it

fr

es

de

da

pl

nl

no

sv

et

fi

pt

hu

lt

cs

ro

sr

sk

sl

ru

uk

bg

el

tr

ar

is

Storage	When the product is taken out of the autoclave and before use, check the package for damage. In case of perforated package, re-pack the product and repeat the sterilization procedure.
Manufacturer contact	See Manufacturer contact details on the IFU.

¹Process validated with Empower™ enzymatic solution:

- Concentration: 30ml/3L
- Temperature: 35°C ± 2°C (95°F±3°F)

²Process validated with undiluted CaviCide™ disinfectant solution:

- Concentration: undiluted
- Contact time: 3 minutes

³Process validated with Miele Vario TD program:

- cleaning agent: Neodisher MediClean

Storage and Disposal

Refer to outer package for expiration date. Do not use after expiration date. For proper disposal always follow local and national regulations.

Disclaimer

The instructions provided above have been approved by Kerr Corporation as being CAPABLE of preparing a medical device, including reusable ones. It remains the responsibility of the dental professional / user to ensure that the reprocessing actually performed by equipment, materials and personnel is done per the instructions provided above in order to achieve the desired result. Any deviation by the dental professional / user from the instructions provided above will be done at your own risk; the Kerr Corporation will not be able to address any requests for refund or exchanges under guarantee of any products that have not been handled or reprocessed in accordance with the above instructions.

Hawe™ X-Ray Centring Devices

Centring Devices.

General Information

X-Ray Centring Devices are intended to facilitate the exact positioning of rectangular or other special X-Ray tubes when using the X-Ray Film and Phosphor Plate Holders or X-Ray Sensor Holders.

Indications for Use

To assist the direct or indirect imaging of the bitewing, periapical or endodontic radiographs for human diagnosis purposes.

Contraindications

None known.

Warnings

None known.

Precautions

- These products are intended to be used by qualified dental professionals on the general patient population.
- Wear personal protective equipment when using and handling these products.
- This device needs to be processed prior to first use according to section « Cleaning, Disinfection, and Sterilization ».
- Discard the device if deformation, breakages, discoloring or missing parts are noticed.

- To avoid cross contamination, this product must be reprocessed after each patient use following this instruction for Use.

Adverse Events

If a serious incident occurs with this medical device, report it to the manufacturer and to the competent authority for that country in which the user and/or patient is established.

Step-by-Step Instructions

1. Remove centring ring from the X-Ray Holder's stem.
2. Assemble the X-Ray Centring Device on the X-Ray Film and Phosphor Plate Holder's or X-Ray Sensor Holder's centring ring. The X-Ray Centring Devices stay in place on the ring by a slight press-fit.
3. Mount the Centring ring on the X-Ray Holder's stem. (p. 110)

Cleaning and Disinfection

Instructions for Cleaning/Disinfecting the Hawe Centring Devices	
WARNINGS	Start the reprocessing procedure immediately after use because dried contaminants may negatively influence the cleaning effort.
Limitations on processing	Product is not sterilizable.
Initial treatment at the point of use	Remove the X-Ray Centring Devices from the X-Ray Film and Phosphor Plate Holders or X-Ray Sensor Holders' Centring ring.
Cleaning: Manual	• Remove any visible debris or organic matter from the instrument by scrubbing in warm tap water (> 40°C-104°F) until no visible debris is observed. Never use metal brushes or wire wool. • With the naked eye, inspect instruments for debris. Repeat previous step if there is still debris. • Thoroughly scrub the instrument surfaces, for at least 1 minute, with a soft brush dampened with an enzymatic solution¹ following the manufacturer's instructions with regard to the concentration and temperature. Do not use detergents containing strong alkalis, strong acids, hydrogen peroxide or organic solvents. • Soak the instrument in an enzymatic solution following the manufacturer's instructions¹ with regard to the concentration, temperature and duration. Do not use detergents containing strong alkalis, strong acids, hydrogen peroxide or organic solvents. • Alternatively to the previous step, place the instrument in an ultrasonic cleaner with a cleaning or enzymatic agent following manufacturer's instructions¹ with regard to the concentration, temperature and duration.

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Disinfection: Manual	• Rinse all parts for at least 30 sec with warm tap water (> 40 °C/104°F). • Soak the instrument in a disinfectant solution following the manufacturer's instructions² with regard to the concentration, temperature and duration. Do not use disinfectant containing strong alkalis, strong acids, hydrogen peroxide or organic solvents. • Alternatively to the previous step, place the instrument in an ultrasonic cleaner with disinfecting agent following manufacturer's instructions² with regard to the concentration, temperature and duration.
Drying	Dry with clean debris-free soft cloth.
Maintenance, Inspection and Testing	Inspect and discard the instrument if there is evidence of plastic deformation or cracks, bent parts, brittle material or missing parts. The examination can be made with the naked eye or with a magnifier under normal lighting conditions.
Storage	Store according to section "Storage and Disposal".
Manufacturer contact	See Manufacturer contact details on the IFU.

¹Process validated with Empower™ enzymatic solution:

- Concentration: 30ml/3L
- Temperature: 35°C ± 2°C (95°F±3°F)

²Process validated with undiluted CaviCide™ disinfectant solution:

- Concentration: undiluted
- Contact time: 3 minutes

³Process validated with Miele Vario TD program:

- cleaning agent: Neodisher MediClean

Storage and Disposal

Refer to outer package for expiration date. Do not use after expiration date. For proper disposal always follow local and national regulations.

Disclaimer

The instructions provided above have been approved by Kerr Corporation as being CAPABLE of preparing a medical device, including reusable ones. It remains the responsibility of the dental professional / user to ensure that the reprocessing actually performed by equipment, materials and personnel is done per the instructions provided above in order to achieve the desired result. Any deviation by the dental professional / user from the instructions provided above will be done at your own risk; the Kerr Corporation will not be able to address any requests for refund or exchanges under guarantee of any products that have not been handled or reprocessed in accordance with the above instructions.

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©Kerr Corporation, 2021. All rights reserved.
Kerr Australia Pty. Ltd.
Kerr Italia S.r.l.