




SOOOO COMFORTABLE
IS **IDENTIUM**®

PRECISION IMPRESSION



KETTENBACHDENTAL

Simply intelligent

IDENTIUM[®] – THE UNIQUE MATERIAL COMBINATION!

A new generation of precision impression material

Combines the advantages of two proven impression materials (A-silicone and polyether) - Vinylsiloxanether[®].




ONE-STEP IMPRESSION

- Monophase impressions
- Double-mix impressions




INDICATIONS

- Crowns and bridges
- Implants
- Inlays and onlays
- Combined restorations (fixation impression)
- Pick-up impressions
- Functional impressions



Double-mix impression




Implant impression



IDENTIUM® – COMFORTABLE IN ALL DIMENSIONS

Removal from the mouth - as light as a feather

Identium® vs. Impregum™: **2.6 times** lower withdrawal forces



By removing from the patient mouth without exerting force, you avoid deformations and irritations.


For your precision.



Advantages your patient will appreciate too

Fast final setting time

with sufficient processing time and short intraoral setting time



The sufficient processing time allows you to work calmly.
For your precision.

Outstanding flowability



With excellent flowability, you also accurately reproduce even the narrowest gaps.
For your precision.



- **Odorless and tasteless material** reduces retching and vomiting reflex
- **Quick setting and easy removal from the mouth** for a relaxed and pleasant feeling
- **Thixotropic flowability** avoids stains on skin and clothing

COMFORTABLE AND PRECISE IMPRESSIONS



EASY REMOVAL FROM THE MOUTH

SHORT INTRAORAL SETTING TIME

OUTSTANDING FLOWABILITY


Available in viscosities
Heavy, Medium and Light



Selectable setting times
regular (4:30 min.)
fast (3:30 min.)



Compatible with all common
mixing devices



IDENTIUM® – THE RANGE

Identium®

Identium® Heavy

Identium® Heavy Fast

Intro pack

380 mL cartridge, 50 mL cartridge Identium® Light, 10 mL adhesive,
10 dynamic mixers, 6 yellow mixing tips, 6 yellow intraoral tips

Refill pack

2 × 380 mL cartridges

Identium® Medium

Identium® Medium Fast

Intro pack

380 mL cartridge, 10 mL adhesive
10 dynamic mixers, 1 application syringe

Refill pack

2 × 380 mL cartridges

Identium® Light

Identium® Light Fast

Normal pack

2 × 50 mL cartridges,
8 yellow mixing tips



Kettenbach GmbH & Co. KG
Im Heerfeld 7
35713 Eschenburg, Germany
Phone +49 (0)2774 7050
info@kettenbach.com
www.kettenbach-dental.com


KETTENBACHDENTAL
Simply intelligent

022743_0919

